## **WITNESS FIV**

### **GENERALITES**

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## **WITNESS FIV**

### **GENERAL INFORMATION**

Feline Immunodeficiency Virus (FIV) is a lentivirus of worldwide distribution, first isolated in 1986 from several cats exhibiting signs of immunodeficiency. Prevalence rates vary from less than 1 % to more than 20 %, depending on the cat population (sex, age, behaviour), health status, and geographical area.

FIV is most commonly transmitted by biting. The infection is firstly expressed by a transient primary illness lasting several weeks, followed by a long period of apparent normal health which precedes a terminal immunodeficiency phase characterized by ailments such as generalized lymphadenopathy, stomatitis, gastro-enteritis, anaemia and neurological disorders. Affected cats will finally die from the development of a variety of secondary opportunistic infections or tumours.

Infection by FIV is generally accompanied by a rapid rise of antibody levels, particularly those directed against the transmembrane protein of the virus, gp40. The presence of anti-gp40 antibodies indicates that the cat has been exposed to the virus, even in the case of complete absence of symptoms.

### **TEST INDICATION**

WITNESS FIV test is recommended for use when history and/or clinical signs may suggest FIV infection.

### **TEST PRINCIPLE**

The WITNESS FIV test is a simple test, based on *Rapid Immunomigration (RIM)* technology, using a synthetic gp40 peptide from the FIV transmembrane region.

Sensitized colloidal gold particles bound to anti-FIV antibodies present within the sample (whole blood, plasma or serum) are allowed to migrate along a strip. The complex is then captured on a sensitized reaction line where its accumulation causes the formation of a clearly visible pink band. A control band, located at the end of the reading window ensures that the test was performed correctly.

### **SPECIMEN INFORMATION•**

- The test can be performed on unclotted whole blood anticoagulated with EDTA or heparin, serum or plasma.

- Samples should always be collected with a sterile needle and syringe.
- Haemolysis does not significantly interfere with the test, but strongly haemolyzed samples may partly obscure a weak positive line.

#### Storage

Anticoagulated whole blood samples should preferably be tested immediately after collection but not longer than 4 hours after collection if stored at room temperature.

If testing is further delayed, samples should be kept refrigerated (+2°C and 8°C) for up to 4 days.

For prolonged storage, samples (serum and plasma only) should be kept frozen (-20°C).

#### KIT CONTENTS

- 10 pouches, each containing 1 test device and desiccant.
- 1 Buffer dropper bottle (of 2.8 ml).
- Instructions for use.
- 10 pipettes.

#### GENERAL PRECAUTIONS

1. Do not use components after expiration date.
2. Store the test kit at 2°C - 25°C. Do not freeze.
3. Use the test immediately after opening the sealed pouch (within 10 minutes).
4. Avoid touching or damaging membrane at windows (1), (2), (3).
5. The WITNESS device should be placed on a flat, horizontal surface while performing the test.
6. Use a separate pipette for each sample.
7. Hold pipette and buffer bottle vertically when dispensing sample and buffer respectively.
8. Handle all samples as capable of transmitting FIV.
9. For veterinary use only.

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#### TEST PROCEDURE AND RESULTS

<p><b>1. Sample application</b></p> <ul style="list-style-type: none"> <li>• Tear open a pouch provided and place the test device on a flat horizontal surface.</li> <li>• Holding the provided pipette vertically, transfer one drop of sample to the sample well (1).</li> </ul>	<p><b>2. Buffer dispensing</b></p> <ul style="list-style-type: none"> <li>• Remove the cap from the buffer bottle, hold it vertically and add three drops of buffer to the sample well (1).</li> <li>• Leave the test device flat during migration of sample reagent complex through the reading window.</li> </ul>
<p><b>3. Reading test</b></p> <p>After 10 minutes, observe the presence or absence of pink bands in reading windows (2) and (3).</p> <p><i>Note :</i></p> <ul style="list-style-type: none"> <li>• <i>It is possible to read the test before 10 minutes if two pink bands (respectively in (2) and (3)) are clearly visible.</i></li> <li>• <i>The presence of only one band in reading window (3), prior to the end of the development time (10 mins), does not mean that the test is complete, as weak samples may appear slower than the control band.</i></li> </ul>	<p><b>RESULTS</b></p> <p><b>1. Validation</b></p> <p>Test is validated if a pink band is present in the reading window (3).</p> <p><b>2. Interpretation</b></p> <ul style="list-style-type: none"> <li>• No band in reading window (2), with one band in window (3) : sample is negative for FIV antibodies.</li> <li>• One band in reading window (2), with one band in window (3) : sample is positive for FIV antibodies.</li> </ul> <p><b>Note :</b></p> <ul style="list-style-type: none"> <li>• No band in control window (3) : invalid test.</li> <li>• A test result should always be interpreted in the context of all available clinical information and history for the cat being tested.</li> </ul>

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