

Clinispin[®] CT40 Centrifuge User Manual



WOODLEY
EQUIPMENT COMPANY LTD.

Table of Contents

| | |
|---|-----------|
| 1. Introduction | 4 |
| 2. Intended Use | 4 |
| 3. Symbols | 4 |
| 4. Features | 4 |
| 5. Standard Accessories | 4 |
| 6. Technical Specifications | 5 |
| 7. Standard Parts Listing | 7 |
| 8. Installation | 7 |
| 8.1 Location & Mounting | 7 |
| 8.2 Rotor Installation | 7 |
| 9. User Interface & Display | 9 |
| 10. Operating the Centrifuge | 10 |
| 10.1 Switch the Centrifuge On | 10 |
| 10.2 RPM/RCF Setting | 10 |
| 10.3 Timer Setting | 10 |
| 10.4 ACC/DCC..... | 11 |
| 10.5 PROG | 12 |
| 10.6 Quick Spin | 12 |
| 10.7 Rotor | 12 |
| 10.8 Start & Stop..... | 12 |
| 10.9 Imbalance Detection | 13 |
| 10.10 Emergency Stop/Lid Release..... | 13 |
| 11. Safety Precautions | 13 |
| 12. Maintenance and Cleaning | 14 |
| 13. Transportation & Storage | 15 |
| 13.1 Transport..... | 15 |
| 13.2 Storage | 15 |
| 14. Troubleshooting | 16 |
| 15. Warranty Statement | 17 |
| 16. Product Disposal | 18 |

1. Introduction

The Clinispin CT40 Centrifuge is equipped with a maintenance free drive, a large display and simple interface for efficient operation for daily lab usage. The programmable centrifuge can deliver up to 4500 RPM and can accommodate different types of rotors. It features various programmable modes to save time and add convenience.




2. Intended Use

The Clinispin CT40 is a large capacity centrifuge designed to separate, sediment, spin down aqueous solutions and solvent suspensions of differing densities in compatible sample containers.

NOTE:

Before using the centrifuge, please read this user manual carefully. This user manual is intended to assist with the operation and care of the unit. It is not a document which aids in repair. If repairs are needed, please contact Woodley Equipment Company.

3. Symbols

| | |
|---|--|
|  Warning |  Waste electrical |
|  Fuse | |

4. Features

- Delivers up to 4500 RPM for all compatible rotors
- Brushless DC maintenance free motor drive
- Imbalance detection with auto cutoff feature
- Lid lock safety feature – lid does not open during operation
- Program mode for customised operation
- Speed setting by RPM/RCF mode
- Countdown timer range from 1 to 99 minutes
- Last run memory feature
- Convenient and easy user interface
- Emergency lid release during power outage
- Automatic internal diagnosis and error display






5. Standard Accessories






- Power Cord
- T-Allen Key
- Grease Tube
- User Manual

6. Technical Specifications


| | | |
|--|--|---|
| Motor Type | Brushless DC Motor | |
| Max. Capacity | 360ml (24 x 15ml) | |
| Max. Speed | 4500 RPM | |
| Speed Accuracy & Speed Step | ± 100 RPM | |
| Run Time | 1 min to 99:59 mins & infinite mode | |
| Min. Acceleration Time | ≤35 seconds | ACC 9 / DCC 9 mode, 32 x 6ml swing out rotor (without load) |
| Min. Deceleration Time | ≤35 seconds | |
| Noise Level | 75dB ±10dB (A) | |
| Ambient Temperature | +5°C to +40°C | |
| Relative Humidity | ≤80% | |
| Dimensions | 553 (L) x 427 (W) x 247 (H) mm | |
| Weight | 16.2kg (without rotor) | |
| Power Supply | 230 VAC, 50 Hz | |
| Power Consumption | 391 W | |
| Safety Fuse Rating | 250 VAC, 3A | |
| Altitude | Use up to an altitude of 2000m above MSL (Mean Sea Level) | |
| Pollution Degree | 2 | |
| Environment | For indoor use only | |

| Rotor Type | Fixed Angle | Swing Angle |
|-------------------------|--------------------|--------------------|
| Bucket Type | 4 x 50ml | 4 x 50ml |
| | 24 x 15ml | 24 x 15ml |
| | 32 x 6ml | 32 x 6ml |
| | 16 x 15ml | 16 x 15ml |
| Max. RCF | 2766 xg | 3334 xg |
| Max. Volume | 360ml (24 x 15ml) | |
| Imbalance Cutoff | Automatic | |
| Display Type | 4 Digit, 7 Segment | |

| Swing Out Rotor (Plastic) | | | |
|---|---|--|---|
|  | Rotor Type | | RCF Value (xg) |
| | 4 x 50ml | | 3315 |
| | 16 x 15ml | | 3334 |
| | 24 x 15ml | | 3220 |
| | 32 x 6ml | | 3153 |
| Adaptor | | | |
| 50ml Adaptor (Aluminium) | 15ml Adaptor (Aluminium) | 15ml Glass Tube Adaptor (Aluminium) | 6ml Adaptor (Aluminium) |
|  |  |  |  |

| Fixed Angle Rotor (Plastic) | | | |
|---|---|--|---|
|  | Rotor Type | | RCF Value (xg) |
| | 4 x 50ml | | 2562 |
| | 16 x 15ml | | 2744 |
| | 24 x 15ml | | 2744 |
| | 32 x 6ml | | 2766 |
| Adaptor | | | |
| 50ml Adaptor (Plastic) | 15ml Adaptor (Plastic) | 15ml Glass Tube Adaptor (Plastic) | 6ml Adaptor (Plastic) |
|  |  |  |  |

7. Standard Parts Listing



8. Installation

The Clinispin CT40 is supplied in a packaging box. Upon receiving the centrifuge, gently remove the packaging and take out the centrifuge by holding it from the bottom. Ensure that all the packaging accessories are removed from the rotor chamber. Please keep all the packaging safely for at least two years for warranty purposes.

8.1 Location & Mounting

Place the centrifuge on a flat, solid and level surface. Ensure that all the four feet of the centrifuge stand on the surface firmly. Avoid installing on a slippery surface or surface prone to vibration.

- Ideal ambient temperature is $+25^{\circ}\text{C} \pm 5^{\circ}\text{C}$. Avoid placing the centrifuge in direct sunlight.
- Ensure clearance of at least 30cm from all sides for ease of use.
- Keep away from heat and water to avoid sample temperature issues or centrifuge failure.
- Do not place the centrifuge anywhere that makes it difficult to operate.

8.2 Rotor Installation



Rotor Removal and Replacement Process

To remove or replace the rotor, follow the instructions below.

Removing the Rotor

1. Using the T-Allen Key, loosen the rotor nut by turning it counterclockwise. Do not try to pull the rotor, the rotor will come up automatically.
2. Once the rotor nut has been loosened completely, pull up the rotor vertically.

Replacing the Rotor



Do not interchange the buckets of the fixed angle and swing rotors.


1. To replace or install the rotor, take the rotor and load vertically onto the motor shaft.
2. Place the rotor nut in the centre hole of the rotor onto the motor shaft.
3. Put the T-Allen Key into the rotor nut and turn clockwise to tighten and counterclockwise to loosen the rotor.

NOTE:

- Check the rotor is firmly tightened before running the next program.
- Ensure you do not interchange the buckets of the fixed angle and swing rotors. The warranty will be void if they are interchanged.

Balancing the Rotor


Swing Angle Rotor



Balanced Rotor

Imbalanced Rotor

Fixed Angle Rotor




Balanced Rotor

Imbalanced Rotor

1. Always balance the rotor before centrifugation. Above are examples of balancing the rotors.
2. The samples in the tubes should be of equal volume.
3. If the tubes are not loaded correctly – vibration or imbalance can occur which can cause serious damage to the centrifuge.
4. If the tubes are not loaded symmetrically, the imbalance detector will cut off the

running centrifuge for device and user safety. This will stop the centrifuge and Err 55 will be seen indicating tubes are not loaded symmetrically. To resume operation, load tubes symmetrically and restart the centrifuge.

9. User Interface & Display



| No. | Button | Function |
|-----|---------------------------|--|
| 1 | Start | Press once and the motor will start as per the selected RPM and time if the lid is closed. The START button only works if the lid is closed. |
| 2 | Stop/Open | Press once and the motor will stop. After the motor has stopped, the lid will open automatically. The OPEN button will only work when the unit has come to a complete stop. |
| 3 | RPM/RCF | Press once and the speed display will flash 5 times to set RPM. Select the RPM value to be set for centrifugation. Press and hold the button to switch between RPM and RCF. Select the RCF value to be set for centrifugation. RPM/RCF mode is indicated via LED on the panel. |
| 4 | Time | Press once and the minutes can be set. The timer can be set from 1 min to 99 mins 59 seconds, infinite mode and roll over mode. The display will flash 5 times and then set the selected time. |
| 5 | Increment/Decrement | Press Plus (+) or Minus (-) to increase or decrease the speed, time, program value or program number. |
| 6 | Acceleration/Deceleration | Press once to set acceleration. Press once again to set deceleration. |
| 7 | Rotor | Press once to select the rotor number. Press multiple times to change the rotor number. |
| 8 | PROG | Press and hold to enter program mode. The display will flash 5 times to set the program number. To change the program number, use the Plus (+) and Minus (-) buttons. After 5 flashes, specific program |

| | | |
|---|------------|--|
| | | values can be set by using the Plus (+) and Minus (-) buttons, for confirmation each value will flash 5 times. |
| 9 | Quick Spin | Press and hold the button to run the last set RPM program. |

LEDs for lid open/close indication, run indication and RPM/RCF functions.

Display Element: 4 digit, 7 segment.

10. Operating the Centrifuge

10.1 Switch the Centrifuge On

After connecting the power cord, switch the main power supply ON. Then, switch the power switch located on the rear side of the centrifuge to ON. Check the rotor fitment before use. The centrifuge will not operate when the lid is open.

NOTE:

Maintain a gap of 3 seconds between switching OFF and switching ON again. DO NOT switch OFF and ON again immediately.

10.2 RPM/RCF Setting

- After closing the centrifuge lid, press the RPM/RCF button to select speed setting in RPM/RCF mode. Press the (+) button to increase the RPM/RCF value or press the (-) button to decrease the RPM/RCF value.
- Minimum and Maximum RPM of the centrifuge is 500 RPM to 4500 RPM respectively. No roll over mode.
- The flashing of the RPM value will stop and the RPM value will be stored automatically if no buttons are pressed after 5 flashes of the adjustment.
- RPM/RCF can also be changed while the centrifuge is under operation. Press the RPM/RCF button and use the (+) and (-) buttons to change the value. Changing the RPM/RCF value between the ongoing centrifugation will run the centrifuge at the updated speed for the remaining time as indicated by the timer.
- To set the RCF value, press and hold the RPM/RCF button to change the RCF setting. Press the (+) and (-) buttons to increase or decrease the value. The flashing of the RCF value will stop and the RCF value will be stored automatically if no buttons are pressed after 5 flashes of the adjustment.
- The value will flash 5 times to indicate acceptance.

10.3 Timer Setting

- Operating time can be selected from 1 min to 99:59 mins, or continuous.
- Press the TIME button to change the time setting.
- Press the (+) and (-) buttons to increase or decrease the value. After the display flashes 5 times, the value will be set.
- The centrifuge time set for run between 1 min to 99:59 minutes or operated in infinite time mode or roll over mode.
- Infinite timer will be indicated by][on the display.
- The timer on the centrifuge's countdown timer in the display will be in Minutes/Seconds.

- The time will be shown after 5 seconds.
- The time will be indicated on the TIME LED light.
- The input will be accepted once the value has flashed 5 times. The value will flash 5 times to indicate acceptance.

10.4 ACC/DCC

This button will set the parameters of acceleration or deceleration of the rotor ramp. Press the ACC/DCC button once to set the acceleration from 1-9 by using the +/- buttons for increment/decrement respectively. Similarly, press the ACC/DCC button again to set deceleration from 1-9 by using the +/- buttons for increment/decrement respectively. The values of acceleration and deceleration timings are listed below.

| Acceleration & Deceleration Time | | | | |
|----------------------------------|----------|-----------|----------------|----------|
| Rotor Type: Swing Out | | | Max. RPM: 4500 | |
| Tube Capacity | 4 x 50ml | 16 x 15ml | 24 x 15ml | 32 x 6ml |
| Mode ACC/DCC | Time (S) | Time (S) | Time (S) | Time (S) |
| DL | 478 | 478 | 481 | 478 |
| D0 | 293 | 293 | 297 | 293 |
| A1/D1 | 181/178 | 184/177 | 181/183 | 188/180 |
| A2/D2 | 172/164 | 174/164 | 172/170 | 173/168 |
| A3/D3 | 152/143 | 155/145 | 151/158 | 152/146 |
| A4/D4 | 131/125 | 133/126 | 131/139 | 131/128 |
| A5/D5 | 110/103 | 112/107 | 111/109 | 113/108 |
| A6/D6 | 90/86 | 93/87 | 91/90 | 91/89 |
| A7/D7 | 71/65 | 73/68 | 72/71 | 74/69 |
| A8/D8 | 51/47 | 52/46 | 52/49 | 52/46 |
| A9/D9 | 34/30 | 38/30 | 33/32 | 35/30 |

| Acceleration & Deceleration Time | | | | |
|----------------------------------|----------|-----------|----------------|----------|
| Rotor Type: Fixed Angle | | | Max. RPM: 4500 | |
| Tube Capacity | 4 x 50ml | 16 x 15ml | 24 x 15ml | 32 x 6ml |
| Mode ACC/DCC | Time (S) | Time (S) | Time (S) | Time (S) |
| DL | 483 | 484 | 484 | 484 |
| D0 | 301 | 302 | 301 | 301 |
| A1/D1 | 182/187 | 181/188 | 182/188 | 181/188 |
| A2/D2 | 173/174 | 172/176 | 172/174 | 172/176 |
| A3/D3 | 152/153 | 150/153 | 151/152 | 151/153 |
| A4/D4 | 131/133 | 129/134 | 130/133 | 130/134 |
| A5/D5 | 111/112 | 110/113 | 111/113 | 110/113 |
| A6/D6 | 91/93 | 89/94 | 90/93 | 90/94 |
| A7/D7 | 71/73 | 68/74 | 71/72 | 70/73 |
| A8/D8 | 51/48 | 50/49 | 51/49 | 51/49 |
| A9/D9 | 33/34 | 32/34 | 33/34 | 32/34 |

NOTE:

The readings were taken under 'No Load' conditions so the results may vary when the equipment is loaded.

10.5 PROG

- Press the PROG button to select a program from P01 – P99.
- Each stored program can be used for direct operation.
- Different parameters of RPM/RCF, TIME and ACC/DCC can be set in each program.
- All the parameters save automatically after the value has flashed 5 times.
- RPM/RCF, TIME and ACC/DCC parameter values entered respectively.
- After entering each value, it will flash 5 times to indicate acceptance.

10.6 Quick Spin

Press and hold the QUICK SPIN button for immediate centrifugation of the samples up to previous run RPM. Releasing the button will stop the rotor gradually.

10.7 Rotor

Press the ROTOR button to select the rotor, as listed below. The rotor numbers can be changed by using the +/- buttons for increment or decrement of the rotor number.

| Swing Out Rotor | | | | |
|-----------------|-------------|----------|----------|-------------|
| Rotor No. | Max. Volume | Max. RPM | Max. RCF | Radius (cm) |
| 31 | 32 x 6ml | 4500 | 3152 | 13.9 |
| 32 | 24 x 15ml | 4500 | 3220 | 14.6 |
| 33 | 16 x 15ml | 4500 | 3334 | 14.7 |
| 34 | 4 x 50ml | 4500 | 3315 | 14.6 |

| Fixed Angle Rotor | | | | |
|-------------------|-------------|----------|----------|-------------|
| Rotor No. | Max. Volume | Max. RPM | Max. RCF | Radius (cm) |
| 51 | 32 x 6ml | 4500 | 2766 | 12.2 |
| 52 | 24 x 15ml | 4500 | 2744 | 12.3 |
| 53 | 16 x 15ml | 4500 | 2744 | 12.3 |
| 54 | 4 x 50ml | 4500 | 2562 | 11.5 |

| Swing Out Bucket | | Fixed Angle Bucket | |
|------------------|---------------|--------------------|---------------|
| No. | Max. Volume | No. | Max. Volume |
| CT40-Swing-BH | Bucket Holder | CT40-Fix-BH | Bucket Holder |
| 31 | 32 x 6ml | 51 | 32 x 6ml |
| 32 | 24 x 15ml | 52 | 24 x 15ml |
| 33 | 16 x 15ml | 53 | 16 x 15ml |
| 34 | 4 x 50ml | 54 | 4 x 50ml |

10.8 Start & Stop

Press the START button to start the operation and press the STOP button to stop the ongoing operation. When the centrifuge is running, the RUN LED light will be lit. Pressing the STOP/OPEN button will stop the operation. Once the centrifuge has stopped, the lid will automatically open.


Once the set time has been reached, the centrifuge will stop automatically and the lid will open automatically.

10.9 Imbalance Detection

The centrifuge is equipped with an imbalance detection safety feature. When the rotor is not loaded symmetrically, the imbalance detector activates and will stop the centrifugation. The error "Err 55" will be shown on the display. First, correct the imbalance load. Then, switch the centrifuge OFF and ON again. The values will be same as those set before imbalance. The imbalance detection feature cannot be deactivated, as it is factory fitted.

10.10 Emergency Stop/Lid Release

Disconnect the centrifuge from the main power supply. Wait until the rotor has come to a complete stop (this may take some time). Once the rotor has stopped, pull the lid release thread located on the bottom left (to your right) of the centrifuge. This will open the centrifuge lid.



NOTE:

This method of opening the lid should only be used in case of emergency or power failure.


11. Safety Precautions



Read all safety and usage information provided in this manual carefully before using the centrifuge.

- Do not use the centrifuge in any manner not specified in this manual.
- Only use rotors and spare parts supplied by Woodley Equipment Company for best results and product safety.
- The rotors must be loaded symmetrically. Each tube should be counter balanced by another tube of the same weight.
- Do not use a centrifuge or rotor that has not been correctly installed or shows any sign of damage.
- The rotor must always be fitted securely. If the centrifuge makes unusual noises during operation, the rotor fitment needs to be checked. Switch OFF the centrifuge immediately by pressing STOP. Check the rotor fitment and fasten it if required.
- When first unboxing the centrifuge, please plug the centrifuge in and close the lid.
- Never move the centrifuge whilst in operation.
- Prior to centrifugation, the tubes should be visually inspected for any damage. Damaged tubes must not be used. This is because broken tubes can result in contamination, sample loss and can create imbalance in the centrifuge which can result in further damage to the centrifuge and its accessories.

- Do not fill tubes while they are in the rotor. Liquids spillage may damage the centrifuge. If liquids are spilled in the rotor or rotor chamber, the centrifuge must be cleaned carefully and properly before being used again.
- The centrifuge must be used for the specified applications only. It must not be operated in a hazardous or flammable environment, and must not be used to centrifuge explosive or highly reactive material. Also, do not place any potentially hazardous material within the clearance area.
- Equipment if used in any manner not specified in this manual or by the manufacturer can result in voiding the product warranty.
- Repairs must only be performed by Woodley Equipment authorised personnel.
- Do not lean on the equipment. It may damage the equipment or harm the operator.
- In the event of contamination caused by aggressive agents, the rotor must be cleaned immediately using a natural cleaning liquid. This is particularly important for the bores of the tubes. If there are any signs of damage, contact Woodley Equipment Company.
- Before using cleaning or decontamination methods other than those stated in this manual, contact Woodley Equipment Company to ensure that the intended method will not damage the centrifuge.
- For safety, Woodley Equipment Company have provided protective earthing with the power supply. Please ensure the power supply is earthed.
- Be sure to close the tube lids tightly prior to centrifugation. Open tube lids can be torn off during centrifugation and can damage the rotor or centrifuge.
- Rotor and adaptors are high grade components which are subject to extreme mechanical strain. Even slight scratches and tears can lead to serious internal material damage. Ensure to check the rotor for any signs of damage before use. Rotors showing visible signs of corrosion or mechanical damage should not be used.
- The centrifuge's capacity must not be exceeded.
- **For safety and air-tight sealing to minimise the air noise, higher force needs to be applied on the lid to close the latch properly.**
- The Safety Fuse provided is of 3 Amp configuration which can be replaced by the operator. The same will protect the machine circuit during an electrical fault or overload.



12. Maintenance and Cleaning

- The rotor and the outside of the centrifuge should be cleaned regularly with a damp (with water) cloth.
- Ensure that the centrifuge is not plugged in while cleaning.
- Wear protective clothing and safety glasses while operating and cleaning the centrifuge.
- The brushless motor in the centrifuge requires no routine maintenance. Any required service should be performed by Woodley Equipment authorised personnel only. Repairs performed by unauthorised personnel may void the warranty.
- Always keep the centrifuge housing, rotor chamber and rotor clean. All parts should be wiped down periodically with a soft cloth.

NOTE: Liquid should not come into contact with the motor.

- After cleaning, ensure that all parts are dry before using.
- Regular cleaning of the rotor is important.

- If the rotor chamber needs cleaning, use a cloth or sponge dampened with a neutral detergent solution.
- Do not place the rotor into the cleaning solution.
- If corrosive, toxic or pathogenic bacteria are accidentally spilled in the rotor or rotor chamber, the centrifuge must be decontaminated thoroughly.
- A Grease Tube has been provided as a standard accessory. A small amount of grease has to be applied on the Motor Shaft (Threaded Parts, Grooves) every 2 weeks to prevent rotor and shaft jam. This will ease the installation process in case of regular rotor change.

13. Transportation & Storage

- Only use the original packaging during transportation.
- For longer distances, utilise transportation assistance such as robust trucks.
- Avoid knocking, harsh shaking or jolting the centrifuge.
- Always retain the packaging material and transportation protections for longer storage or transportation.

13.1 Transport



- Before transporting the centrifuge, the transport securing device must be installed.
- When the centrifuge and its accessories are transported, the following ambient conditions must be complied with:
 - Ambient Temperature: +5°C to +40°C
 - Relative Humidity: ≤80%, non-condensing

13.2 Storage

- The centrifuge and its accessories should only be stored in dry environments, away from direct sunlight.
- When the centrifuge and its accessories are stored, the following ambient conditions must be complied with:
 - Ambient Temperature: +5°C to +40°C
 - Relative Humidity: ≤80%, non-condensing

14. Troubleshooting

The centrifuge has a self-diagnostic function. If a problem occurs, an error/warning code will be displayed on the display screen and the operator can determine the error using the following table.

| Error | Problem | Solution |
|---|--|---|
| No display | No main power connection | Check the power and properly plug in the mains cable at both ends. |
| | Power failure | Check the mains fuse of the lab. |
| | Improper connection | Connect the adaptor properly. |
|   | Lid not closed properly | Close the lid correctly. |
| | Error with lid closing and opening mechanism | Contact Woodley Equipment Company. |
| Err 55 | Rotor not loaded symmetrically with equal weights | Load rotor symmetrically and restart the centrifuge. |
| Centrifuge lid cannot be opened | Rotor is still spinning | Wait for the rotor to come to a complete stop. |
| | Power failure | Use the emergency lid release after the rotor has come to a complete stop. |
| Centrifuge shakes during acceleration and unusual running noise | Rotor not loaded symmetrically | Stop the centrifuge, load the rotor symmetrically and restart operation. |
| | Either a broken tube, damage to the rotor or motor is the cause of the unusual run noise | Replace the broken tube. For damaged rotor/motor, contact Woodley Equipment Company. |
| | Rotor damaged | Remove and change the rotor. |
| Display error | Loose connection of display | Contact Woodley Equipment Company. |
| Err 1 | Latch damaged, latch jammed | Contact Woodley Equipment Company. |
| | Lid not closed/locked | Close/lock the lid. |
| Err 52 | Rotor stuck or incorrect operating voltage | Turn the centrifuge OFF. Check the rotor fitment or apply correct 230 VAC \pm 10 VAC operating voltage. |
| Power tripping | Cable not fit properly | Remove cable and connect properly. |

IMPORTANT NOTE:

- If the system gets hot or overheats due to over current, switch the centrifuge OFF and ON (restart), and check it again.
- Maintain a 3 second gap between turning the centrifuge OFF and ON. Immediately turning the centrifuge back on can lead to a reset, erasing last run memory.
- If the motor overheats, there will be fluctuation in speed value. Allow centrifuge to cool down for at least 30 minutes. Do not perform any operations for 30 minutes.

15. Warranty Statement

This product is warranted to be free from defects in material and workmanship for a period of one (1) year from the date of purchase. Your product will be duly repaired upon prompt notification in compliance with the following conditions:

This warranty is valid only if the product is used for its intended purpose and within the guidelines specified in this manual. This warranty does not cover damage caused by accidents, neglect, misuse, improper service, natural forces or other causes not arising from defects in original material or workmanship. This warranty does not cover any incidental or consequential damages, commercial loss or any other damages from the use of this product.

The warranty is invalidated by any non-factory modification, which will immediately terminate all liabilities on Woodley Equipment Company for the products or damages caused by its use. The buyer and its customer shall be responsible for the product or use of products as well as any supervision required for safety. If requested, the products must be returned well packed and all shipping charges must be paid.

Some states do not allow limitation on the length of implied warranties or the exclusion or limitation of incidental or consequential damages. This warranty gives you specific legal rights. This warranty is given expressly in lieu of all other warranties, expressed or implied.

The purchaser agrees that there is no warranty of merchantability or of fitness for any intended purpose and that there are no other remedies or warranties, expressed or implied, which extend beyond the description on the face of the agreement. This warranty is only applicable to the original purchaser.

Products received without proper authorisation will not be entertained. All items returned for service should be returned in the original packaging or other suitable carton, padded to avoid damage. Woodley Equipment Company is not responsible for damage incurred by improper packaging.

| | |
|---|----------------|
| For your reference, make a note of the serial number and date of purchase here. | |
| Serial No.: | Purchase Date: |

16. Product Disposal

In case the product is to be disposed of, the relevant legal regulations are to be observed.


Information on the disposal of electrical and electronic devices in the European Community:













The disposal of electrical devices is regulated within the European Community by national regulations based on EU Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE). According to these regulations, any devices supplied after 13.06.05 in the business-to-business sphere, to which this product is assigned, may no longer be disposed of in municipal or domestic waste. They are marked with the following symbol to indicate this.



As disposal regulations within the EU may vary from country to country, please contact your supplier if necessary.










Swing Out Rotor

| | | |
|---|---|---------------|
|  | Max. G Force/RCF | 3334 xg |
| | Max. Rotation Speed | 4500 RPM |
| | Acceleration & Deceleration Time | ≤35s |
| | Max. Load (Adapter, Tube Holder, Tube & Contents) | 16 x 161.2 gm |
| | Rotor Weight | 866.4 gm |

| No. | Tube | Tube Type | Adaptor | Adaptor Type | Holder | Holder | Tubes/ Holder | Max. G Force/ RCF | |
|-----|---|-------------------------------|---|--------------|---|-----------------|------------------|------------------------|--------|
| | | Capacity | | Diameter | | Type | | Max. Rotation Speed | |
| | | Diameter | | Material | | Material | | Tubes/ Rotor | Radius |
| | | Length | | Inner Depth | | Dia | | | |
| | | Weight | | Weight | | Weight | | | |
| 1 |  | Conical Centrifuge Tube w/Cap |  | Spherical |  | 4 x 50 ml | 1/1 | 3315 | |
| | | 50 ml | | Ø 31.5 mm | | Bucket Holder | | 4500 RPM | |
| | | Ø 29 mm | | SS | | Aluminium Alloy | 4/1 | 14.6 cm | |
| | | 116.6 mm | | 98.3 mm | | Ø 33 mm | | | |
| | | 12.8 gm | | 41.7 gm | | 152.7 gm | | | |
| 2 |  | Conical Centrifuge Tube w/Cap |  | Spherical |  | 16 x 15 ml | 4/1 | 3334 | |
| | | 15 ml | | Ø 18 mm | | Bucket Holder | | 4500 RPM | |
| | | Ø 15 mm | | SS | | Aluminium Alloy | 4/1 | 14.7 cm | |
| | | 120.5 mm | | 98.3 mm | | Ø 19.5 mm | | | |
| | | 6.6 gm | | 22.1 gm | | 131.4 gm | | | |
| 3 |  | Glass Tube |  | Spherical |  | 24 x 15 ml | 6/1 | 3220 | |
| | | 15 ml | | Ø 18 mm | | Bucket Holder | | 4500 RPM | |
| | | Ø 15 mm | | SS | | Aluminium Alloy | 4/1 | 14.6 cm | |
| | | 120.5 mm | | 98.3 mm | | Ø 19.5 mm | | | |
| | | 6.6 gm | | 22.1 gm | | 122.3 gm | | | |
| 4 |  | Conical Centrifuge Tube w/Cap |  | Spherical |  | 32 x 6 ml | 8/1 | 3152 | |
| | | 7 ml | | Ø 14 mm | | Bucket Holder | | 4500 RPM | |
| | | Ø 12.70 mm | | SS | | Aluminium Alloy | 4/1 | 13.9 cm | |
| | | 107 mm | | 89.5 mm | | Ø 16.5 mm | | | |
| | | 8 gm | | 15..7 gm | | 94.2 gm | | | |

Note: Rotor maximum service life 30,000 cycles (from first use).

Fixed Angle Rotor

|  | | | | Max. G Force/RCF | | 2766 xg | | |
|---|---|-------------------------------|---------|---|---|---------------|--------------|---------------------|
| | | | | Max. Rotation Speed | | 4500 RPM | | |
| | | | | Acceleration & Deceleration Time | | ≤35s | | |
| | | | | Max. Load (Adapter, Tube Holder, Tube & Contents) | | 32 x 154.2 gm | | |
| | | | | Rotor Weight | | 1465.0 gm | | |
| No. | Tube | Tube Type | Adaptor | Adaptor Type | Holder | Holder | Tubes/Holder | Max. G Force/RCF |
| | | Capacity | | Diameter | | Material | | Type |
| | | Diameter | | Diameter | | Dia | | Max. Rotation Speed |
| | | Length | | Material | | Inner Depth | | Radius |
| | | Weight | | Inner Depth | | Weight | | |
| | | Weight | | Weight | | Weight | | |
| 1 |  | Conical Centrifuge Tube w/Cap | - | - |  | 4 x 50 ml | 1/1 | 2562 |
| | | 50 ml | | | | Bucket Holder | | |
| | | ∅ 29 mm | | | | Plastic | 4/1 | 4500 RPM |
| | | 116.6 mm | | | | ∅ 32.47 mm | | |
| | | 12.8 gm | | | | 84 mm | | |
| | | | | 136.7 gm | | 11.5 cm | | |
| 2 |  | Conical Centrifuge Tube w/Cap | - | - |  | 16 x 15 ml | 1/1 | 2744 |
| | | 15 ml | | | | Bucket Holder | | |
| | | ∅ 15 mm | | | | Plastic | 4/1 | 4500 RPM |
| | | 120.5 mm | | | | ∅ 17.97 mm | | |
| | | 6.6 gm | | | | 84 mm | | |
| | | | | 145.2 gm | | 12.3 cm | | |
| 3 |  | Glass Tube | - | - |  | 24 x 15 ml | 6/1 | 2744 |
| | | 15 ml | | | | Bucket Holder | | |
| | | ∅ 15 mm | | | | Plastic | 4/1 | 4500 RPM |
| | | 120.5 mm | | | | ∅ 17.97 mm | | |
| | | 6.6 gm | | | | 84 mm | | |
| | | | | 151.8 gm | | 12.3 cm | | |
| 4 |  | Conical Centrifuge Tube w/Cap | - | - |  | 32 x 6 ml | 8/1 | 2766 |
| | | 7 ml | | | | Bucket Holder | | |
| | | ∅ 12.70 mm | | | | Plastic | 4/1 | 4500 RPM |
| | | 107 mm | | | | ∅ 14.98 mm | | |
| | | 8 gm | | | | 84 mm | | |
| | | | | 145 gm | | 12.2 cm | | |

Note: Rotor maximum service life 30,000 cycles (from first use).



Woodley Equipment Company Ltd.

UK Office:

Old Station Park Buildings,
St. John Street,
Horwich,
Bolton,
BL6 7NY, UK

Tel: +44 (0) 1204 669033

Email: sales@woodleyequipment.com

Woodley Equipment Inc.

USA Office:

Brooklyn Navy Yard,
63 Flushing Avenue,
Building 212, Suite 209,
Brooklyn, NY 11205

Toll Free: 1-800-471-9200

Tel: 1-315-618-7600

Email: hello@woodleytrialsolutions.com

Web: www.woodleyequipment.com