



## InSight V-IA fPL Evaluation vs. Vcheck V200

✉ Old Station Park Buildings, St Johns Street, Horwich, BL6 7NY

☎ +44 (0) 1204 669033

🌐 [www.woodleyequipment.com](http://www.woodleyequipment.com)

The InSight V-IA is an easy to use veterinary Immunoassay Analyser providing accurate and reliable results in 3-15mins. The InSight V-IA uses immunofluorescence technology for accurate results. A competitive binding assay is based upon the competition of labelled and unlabelled analytes for a limited number of antibody binding sites. Unbound antibodies and immunocomplexes migrate along the nitrocellulose membrane towards the test line. The unbound antibodies are then captured by antigens immobilised on the test line. The fluorescent signal intensity reflects the amount of analytes captured and is measured by the InSight V-IA.

### **Feline Pancreatic Lipase (fPL)**

fPL is used as an indicator for diagnosis of pancreatitis in cats.

### **Comparison Items**

#### **InSight V-IA:**

Test Item – fPL

Quantity of Samples – 32 Tests

Lot No. – H064200205

Instrument – InSight V-IA Veterinary Immunoassay Analyser

#### **Bionote (Vcheck):**

Test Item – fPL

Quantity of Samples – 32 Tests

Lot No. – F104D014

Instrument – Vcheck V200


### **Test Results**

Sample	Vcheck V200 (ug/L)	InSight V-IA (ug/L)	Mark
1	1	1.5	
2	1	2.33	
3	1	3.2	
4	1	1.78	
5	1	1.55	Critical
6	3.52	4.23	
7	1	2.06	
8	6.21	5.91	
9	6.85	8.66	
10	1	2.1	
11	1	1.35	
12	1	1.74	
13	5.23	5.55	
14	1	1.45	
15	1.2	1.56	
16	1	1.88	
17	10.2	10.05	
18	18.2	15.8	

19	15.2	15.7	
20	17.5	15.74	
21	1.14	1.62	
22	2.5	1.93	
23	6.1	5.32	
24	3.3	2.58	
25	3.2	3.12	
26	9.6	9.24	
27	6.5	4.4	
28	3	3.27	
29	4.2	4.23	
30	7.1	5.3	
31	10.2	9.86	
32	7.5	6.05	

**Linear Correlation**

<b>R<sup>2</sup></b>	0.9555
<b>R</b>	0.9775



**Accuracy**

fPL		InSight V-IA		
		Positive	Negative	Total
Vcheck V200	Positive	14	0	14
	Negative	1	17	18
	Total	15	17	32

Positive Coincidence Rate	93.3%
Negative Coincidence Rate	100.0%
Total Coincidence Rate	96.9%

### **Repeatability**

Sample Concentration	1	2	3	4	5	6	7	8	9	10	Mean Value	Bias Max %
1.9ug/L	2.18	1.99	2.06	1.87	2.01	1.82	1.74	2.02	1.94	2.31	1.99	5%
6.8ug/L	8	7.89	6.75	6.9	6.78	9.81	6.77	6.53	8	6.57	7.40	9%
15ug/L	18.44	21.38	12.76	16.82	13.25	16.6	14.68	14.73	14.65	15.57	15.89	6%