


# Vcheck Canine NT-proBNP

New Biomarker for the Diagnosis of Heart Failure



## What is NT-proBNP?

The pro-hormone (proBNP) is produced by cardiac myocytes and increases due to increased myocardial wall stress. Upon release into the blood, it is cleaved into BNP and NT-proBNP. Due to its longer half-life, NT-proBNP is better suited as a diagnostic parameter.



## What NT-proBNP levels tell us

In dogs, NT-proBNP is correlated with heart size and systolic function, suggesting that the concentrations can be used to detect animals with early disease.

### Distinguishes cardiac from respiratory disease

- In dogs with dyspnea requiring emergency care
- Differentiates cardiac and respiratory causes of respiratory signs

### Staging of Myxomatous Mitral Valve Degeneration (MMVD)


- Differentiates dogs with MMVD with and without congestive heart failure
- Chronic monitoring in dogs with MMVD

### Detects Dilated Cardiomyopathy in Doberman Pinschers

- Highly sensitive and specific for detecting occult DCM
- Predicts survival in Dobermans at high risk

**NT-proBNP should be interpreted with other diagnostic techniques, such as echocardiography, thoracic radiography, history and clinical signs, improving the accuracy of diagnosis.**


## NT-proBNP predicts the survival time of Dobermans



With low NT-proBNP,  
6.1 times longer !

## Clinical Algorithm

### NT-proBNP testing in dogs



*NT-proBNP provides insights to distinguish cardiac from non-cardiac respiratory distress*

# Vcheck Canine NT-proBNP


## Key Features

- **Quantitative Analysis**  
Actual results (instead of an estimate like others)
- **Proven Accuracy and Reproducibility**  
Correlated against an ELISA method from laboratories
- **Rapid and accurate results**  
Simple procedure and quick results within 15 min.



## Performance


- Strong correlation ( $R^2=0.952$ ) with an ELISA method (from company 'T' laboratories)



# Vcheck Canine NT-proBNP

## Performance

- Vcheck NT-proBNP levels based on stages of ACVIM System



	Control	Stage B1, B2	Stage C	Stage D
ACVIM Stage	Healthy	Asymptomatic (+ Left heart volume increased)	CHF signals	Terminal stage

## Ordering Information

Product No.	Product Name	Product Type	Packing Unit
VCF132DC	Vcheck Canine NT-proBNP	Device	5 Tests/Kit

# Vcheck Canine NT-proBNP


## Specifications

- Species : Dog
- Sample : Serum 100  $\mu$ l
- Testing Time : 15 minutes
- Measurement : Quantitative
- Measurement Range : 500 – 10,000 pmol/L
- Storage Condition : 2 - 8  $^{\circ}$ C



## Test Procedure


- 1 Add 100  $\mu$ l of the sample to the assay diluent tube



- 2 Mix well by using a 100  $\mu$ l pipette



- 3 Add the mixed sample (100  $\mu$ l) into the test device



**Samples should be centrifuged and tested immediately after collection.  
Alternatively, refrigerate and use within 24 hours or freeze.**

\* \* Degradation of NT-proBNP may occur if stored at room temperature or refrigerated for more than 24 hours, causing false negative results.

## Reference Ranges

< 900 pmol/L	900 – 1,800 pmol/L	> 1,800 pmol/L
Normal	Suspected* Additional diagnostics are recommended	Abnormal* Additional diagnostics are recommended

\* 'Abnormal' or 'Suspected' NT-proBNP test results should always be interpreted in combination and other diagnostic findings, such as an echocardiogram.

\*\* Concentration over 735 pmol/L in Doberman Pinschers indicates an increased risk for occult dilated cardiomyopathy.