

Vcheck Canine NT-proBNP

BIONOTE Marketing team

Apr. 2020



Vcheck Canine NT-proBNP

01 NT-proBNP

02 Product Introduction

- Vcheck Canine NT-proBNP


01 NT-proBNP

- What is NT-proBNP?
- What NT-proBNP levels tell us
- Algorithm: NT-proBNP testing in dogs

01 NT-proBNP

What is NT-proBNP?

- B-type natriuretic peptide (proBNP)
 - produced in the muscle cells of the heart
 - increases with excessive stretching of the cells
⇒ correlated to the severity of the underlying heart disease
- proBNP is cleaved into BNP and NT-proBNP
(※ NT-proBNP: N-terminal pro-B type natriuretic peptide)
- **NT-proBNP** is stable and has a long half-life, making it a more desirable biomarker.
 - ⇒ used to assess the magnitude of cardiac muscle stretching
 - ⇒ proportionate to the severity of cardiac disease



01 NT-proBNP

NT-proBNP in Dogs

What NT-proBNP levels tell us

- **Distinguishes cardiac from respiratory disease**
 - ✓ In dogs with dyspnea requiring emergency care
 - ✓ Differentiates cardiac and respiratory causes of respiratory signs
- **Staging of Myxomatous Mitral Valve Degeneration (MMVD)**
 - ✓ Differentiates dogs with MMVD with and without congestive heart failure
 - ✓ Chronic monitoring in dogs with MMVD
- **Detects Dilated Cardiomyopathy in Doberman Pinschers**
 - ✓ Highly sensitive and specific for detecting occult DCM
 - ✓ Predicts survival in Dobermans at high risk



▲ A dog presenting difficulty breathing (dyspnea)


Vcheck Canine NT-proBNP

What NT-proBNP levels tell us

- 01 Distinguishes cardiac from respiratory disease
- 02 Staging of Myxomatous Mitral Valve Degeneration (MMVD)
- 03 Detects Dilated Cardiomyopathy in Doberman Pinschers

What NT-proBNP levels tell us

01 Distinguishes cardiac from respiratory disease



- NT-proBNP concentration was higher in CHF dogs (5,110 pmol/L) compared to those with noncardiac respiratory distress (1,287 pmol/L).
- A cut-off >2,447 pmol/L : Discriminates CHF from non-cardiac respiratory distress (Sensitivity 81.1%, Specificity 73.1%)
- NT-proBNP biomarker is useful for discriminating CHF from non-cardiac respiratory distress.

✓check Canine NT-proBNP


What NT-proBNP levels tell us

- 01 Distinguishes cardiac from respiratory disease
- 02 Staging of Myxomatous Mitral Valve Degeneration (MMVD)
- 03 Detects Dilated Cardiomyopathy in Doberman Pinschers

What NT-proBNP levels tell us

02 Staging of Myxomatous Mitral Valve Degeneration (MMVD)

- NT-proBNP in staging of MMVD in dogs



◀ In dogs with MMVD categorized according to the CHIEF classification


CHIEF	Interpretation
A	Control
B1, B2	MMVD only
C2, C3, D	MMVD + CHF

⇒ NT-proBNP values of all disease stages (B1-D) were significantly higher than NT-proBNP of the control group.

What NT-proBNP levels tell us

02 Staging of Myxomatous Mitral Valve Degeneration (MMVD)

- NT-proBNP in staging of MMVD in dogs



Cut-off 1207.5

- Sens: 83.3 %
- Spec: 84.8 %

◀ Sensitivity and specificity of NT-proBNP at diagnostic cut-off values of 1207 pmol/L

CHIEF	Interpretation
A	Control
B1, B2	MMVD only
C2, C3, D	MMVD + CHF

In canine MMVD, NT-proBNP is useful to discriminate between asymptomatic dogs and dogs with CHF.

✓check Canine NT-proBNP

What NT-proBNP levels tell us

- 01 Distinguishes cardiac from respiratory disease
- 02 Staging of Myxomatous Mitral Valve Degeneration (MMVD)
- 03 Detects Dilated Cardiomyopathy in Doberman Pinschers

What NT-proBNP levels tell us

03 Detects Dilated Cardiomyopathy in Doberman Pinschers

- Diagnostic utility of NT-proBNP assay

▼ Diagnostic utility of NT-proBNP assay and Holter for the detection of ODCM in 155 Dobermans

Criteria	Se	Sp	PPV	NPV	Accuracy
NT-proBNP > 457	69.9	80.5	76.1	75.0	75.5
NT-proBNP > 900	32.9	93.9	82.8	61.1	65.2
Holter, NT-proBNP, or both > 457	94.5	87.8	87.3	94.7	91.0
Holter, NT-proBNP, or both > 900	84.9	95.1	93.9	87.6	90.3



The combination of NT-proBNP assay and Holter detects ODCM with high accuracy

- Combination of NT-proBNP and Holter to detect ODCM yielded **sensitivity of 94.5%**, **specificity of 87.8%**, and **accuracy of 91.0%**. (the presence of > 50 VPCs during Holter)




▲ 24-hour Holter monitoring

What NT-proBNP levels tell us

03 Detects Dilated Cardiomyopathy in Doberman Pinschers

- Diagnostic utility of NT-proBNP assay




- The median survival time of Dobermans:
 - NT-proBNP > 900 pmol/L: 284 days (6.1 times shorter)
 - NT-proBNP < 900 pmol/L: 1743 days

01 NT-proBNP

NT-proBNP in Dogs

Algorithm: NT-proBNP testing in dogs



02 Product Introduction

Vcheck Canine NT-proBNP

- Specifications
- Key Features
- Test Procedure
- Reference Range
- Performance

02 Product Introduction

Vcheck Canine NT-proBNP

- Specifications



- ✓ Species : Dog
- ✓ Sample : Serum 100 μ l
- ✓ Testing Time : 15 minutes
- ✓ Measurement : Quantitative
- ✓ Measurement Range : 500 – 10,000 pmol/L
- ✓ Storage Condition : 2 - 8 $^{\circ}$ C

Product No.	Product Name	Product Type	Packing Unit
VCF132DC	Vcheck Canine NT-proBNP	Device	5 Tests/Kit

02 Product Introduction

Vcheck Canine NT-proBNP

- Key Features

- ✓ **Quantitative Analysis**

- Actual results (instead of an estimate like others)

- ✓ **Proven Accuracy and Reproducibility**

- Correlated against an ELISA method from laboratories

- ✓ **Rapid and accurate results**

- Simple procedure and quick results within 15 min.




02 Product Introduction

Vcheck Canine NT-proBNP


- Test Procedure

- ✓ Samples should be centrifuged and tested immediately after collection.
Alternatively, refrigerate and use within 24 hours or freeze.


* Degradation of NT-proBNP may occur if stored at room temperature or refrigerated for more than 24 hours, causing false negative results.



Add 100 µl of the sample
to the assay diluent tube



Mix well
by using a 100 µl pipette



Add the mixed sample
(100 µl)

* Stored samples should be placed at room
temperature 30 min. before use.

02 Product Introduction

Vcheck Canine NT-proBNP

- Reference Range

< 900 pmol/L	900 – 1,800 pmol/L	> 1,800 pmol/L
Normal	Suspected* Additional diagnostics are recommended	Abnormal* Additional diagnostics are recommended

* 'Abnormal' or 'Suspected' NT-proBNP test results should always be interpreted in combination and other diagnostic findings, such as an echocardiogram.

** Concentration over 735 pmol/L in Doberman Pinschers indicates an increased risk for occult dilated cardiomyopathy.


02 Product Introduction

Vcheck Canine NT-proBNP

- Performance

Comparative evaluation of Canine NT-proBNP

- ✓ Vcheck Canine NT-proBNP has a strong correlation ($R^2=0.952$) with an ELISA method (from company 'I' lab)




$y = 0.9658x - 11.299$
 $R^2 = 0.9521$

02 Product Introduction

Vcheck Canine NT-proBNP

- Performance

NT-proBNP levels based on stages of ACVIM System



MMVD Stage	
By ACVIM (American College of Veterinary Internal Medicine)	
Control	Predisposed to heart disease
Stage B1, B2	Asymptomatic (+ Left heart volume increased)
Stage C	CHF signals
Stage D	Terminal stage

Thank you

BIONOTE Marketing team

Apr. 2020

