

V200 cCRP 2.0 Quick User Guide

lease read these instructions carefully before carrying out the test

sample type: Serum or plasma (Lithium heparin)

test Procedure

All reagents and samples must be at room temperature (15-30°C) before testing

1. Turn on Bionote V200 analyser by holding the power switch for a few seconds. When the main screen appears, select 'Standard Test'
2. Enter operator ID and patient ID, press enter then press Ok
3. Remove test cartridge from pouch and insert into the analyser when prompted. Push the cartridge into the test port until it clicks. The analyser will check the cartridge.
4. 'Apply Sample' will appear on the screen.
5. Using the 5ul micropipette, push down the plunger of the pipette before putting the pipette tip into the sample. Once the pipette tip is in the sample, release the plunger to take up the sample.
6. Add the sample to the diluent tube (supplied in the box).
7. Discard the 5ul micropipette tip, **using the 100ul pipette** and a new pipette tip, mix the sample and diluent a several times by pushing down the plunger and releasing at least 6 times. Draw up the diluted sample by pushing down the plunger of the pipette before putting the pipette tip into the diluted sample. Once the pipette tip is in the sample, release the plunger to take up the sample.
8. Dispense the diluted sample from the pipette onto the sample area of the test cartridge and press 'Start' to begin testing
9. After 5 minutes, the result will be displayed on screen and printed out with the reference ranges.

Reference range	Result
10~20 mg/L	Normal
20~30 mg/L	Equivocal
>30 mg/L	Abnormal

If you require assistance please contact your Woodley Equipment Area Account Manager or Woodley Equipment Company Technical Support Dept. on 01204 669033 option 1.