

# Vcheck Feline NT-proBNP

BIONOTE Marketing team

Jan. 2020



# Vcheck Feline NT-proBNP

## 01 NT-proBNP

## 02 Product Introduction

- Vcheck Feline NT-proBNP


# 01 NT-proBNP

- What is NT-proBNP?
- What NT-proBNP levels tell us
- How can I use the NT-proBNP test?
- Algorithm: NT-proBNP testing in cats
- Factors which increase NT-proBNP concentration

# 01 NT-proBNP

## What is NT-proBNP?

- B-type natriuretic peptide (BNP)
  - produced in the muscle cells of the heart
  - increases with excessive stretching of the cells
    - ⇒ correlated to the severity of the underlying heart disease
- proBNP is cleaved into BNP and NT-proBNP  
(※ NT-proBNP: N-terminal pro-B type natriuretic peptide)
- NT-proBNP is stable and has a long half-life, making it a more desirable biomarker.
  - ⇒ used to assess the magnitude of cardiac muscle stretching
  - ⇒ proportionate to the severity of cardiac disease



# 01 NT-proBNP

## NT-proBNP in Cats

What NT-proBNP levels tell us

- **To screen for occult heart disease**
  - ✓ Before anesthesia
  - ✓ In apparently healthy cats with heart murmurs
  - ✓ At risk breeds - Maine Coon, Ragdoll, Birman, Persian, etc.
- **To determine Cardiac or Respiratory disease**
  - ✓ In cats with respiratory signs such as dyspnea, tachypnea, cough
  - ✓ To differentiate cardiac and respiratory causes of dyspnea
- **To determine the severity of heart disease**
  - ✓ For monitoring stabilization of CHF during hospitalization
  - ✓ For predicting survival in cats with CHF

\* CHF: Congestive Heart Failure



# 01 NT-proBNP

## NT-proBNP in Cats

What NT-proBNP levels tell us


- **To screen for occult heart disease**
  - ✓ Before anesthesia
  - ✓ In apparently healthy cats with heart murmurs
  - ✓ At risk breeds - Maine Coon, Ragdoll, Birman, Persian, etc.
- **To determine Cardiac or Respiratory disease**
  - ✓ In Cats with respiratory signs such as dyspnea, tachypnea, cough
  - ✓ To differentiate cardiac and respiratory causes of dyspnea
- **To determine severity of heart disease**
  - ✓ For monitoring treatment in cats with CHF \*CHF: Congestive Heart Failure
  - ✓ For predicting prognosis

## HCM Hypertrophic Cardiomyopathy

- ✓ Most prevalent feline cardiac disorder
  - ✓ Affected Age: 2, 3 year-old (3 months ~ 10 years)
  - ✓ Process: Thick left ventricular muscle
    - Decreased volume size in left ventricle
- ⇒ Congestive Heart Failure (CHF), Thromboembolism



Normal heart  
(cross section)



Hypertrophic  
cardiomyopathy



# 01 NT-proBNP

## NT-proBNP in Cats

How can I use the NT-proBNP test?

### 01 Comprehensive Evaluation (In Cats without clinical signs)

- ✓ In cats at increased risk of having occult cardiomyopathy  
e.g. heart murmur, gallop heart sound, arrhythmia
- ✓ In high-risk cat breeds,  
e.g. Maine Coon, Ragdoll, Birman, Persian, American Short Hair, Himalayan, Siamese, Sphinx, Burmese, etc.
- NT-proBNP may help the veterinarian assess the likelihood of underlying heart disease in asymptomatic cats and guide the next diagnostic and therapeutic steps. ⇒ Echocardiogram should be performed for the definite diagnosis.
- NT-proBNP assay was 91.2% specific and 85.8% sensitive for detection of underlying heart disease based on findings from echocardiographic examination. *J Vet Intern Med 2011; in press.*
- +) Although NT-proBNP levels is less than 100 pmol/L, high-risk cats should be re-evaluated annually.



# 01 NT-proBNP

## NT-proBNP in Cats

How can I use the NT-proBNP test?

### 02 In Cats with Respiratory signs

- ✓ The presence of respiratory signs (dyspnea, tachypnea, cough)
  - Cardiac disease: underlying cardiomyopathy and congestive heart failure (CHF)
  - Primary respiratory disease: bronchitis/asthma, pneumonia, neoplasia, pleural space disease
- NT-proBNP testing is a useful option to rapidly assess the likelihood that heart disease is present in a cat with respiratory signs.
- In cats with respiratory signs such as dyspnea, tachypnea and cough, if the NT-proBNP is >270 pmol/L, CHF is the most likely cause of the clinical signs. (Concurrent respiratory disease cannot be still ruled out.)
- Cutoff value of 265 pmol/L: Sensitivity of 90.2% and Specificity of 87.9%

▽ a cat having difficulty in breathing






# 01 NT-proBNP

## NT-proBNP in Cats

### Algorithm: NT-proBNP testing in cats



- In cats with respiratory signs such as dyspnea, tachypnea and cough, If the NT-proBNP is >270 pmol/L, Congestive heart failure is the most likely cause of the clinical signs. (However, concurrent respiratory disease cannot be ruled out.)

# 01 NT-proBNP

## NT-proBNP in Cats

- ✓ NT-proBNP is not a stand-alone test
- ✓ NT-proBNP should be interpreted when used in conjunction with findings from the history, physical examination, ECG, auscultation, thoracic radiography and echocardiography ⇒ **Helping achieve a definite diagnosis**

Having a 'biomarker'-guided therapy may be very important as a complement to other testing in cats to improve the accuracy and confidence in diagnosing heart disease or making treatment decisions **for cats with CHF**.

# 01 NT-proBNP

## NT-proBNP in Cats

When NT-proBNP was measured,

More accurate the diagnosis is!




Fig 1. Accuracy of general practitioners' diagnosis in cats with respiratory signs

Much higher the confidence score is!




Fig 2. Confidence of general practitioners' diagnosis in cats with respiratory signs

# 01 NT-proBNP

## NT-proBNP in Cats

### Factors which increase NT-proBNP concentrations

- Feline Hyperthyroidism
- Renal insufficiency or Pre-renal azotemia  
⇒NT-proBNP levels should be interpreted in the context of renal function tests.  
(creatinine > 2.8)
- Systemic hypertension
- Severe arrhythmias
- Pulmonary hypertension (rare)

Elevated  
NT-proBNP concentrations

+

Normal  
Echocardiogram

||

Consider these factors  
which increase NT-proBNP levels!

## 02 Product Introduction

### Vcheck Feline NT-proBNP

- Specifications
- Key Features
- Test Procedure
- Reference Range
- Performance

# 02 Product Introduction

## Vcheck Feline NT-proBNP

- Specifications
  - ✓ Species : Cat
  - ✓ Sample : Serum, Plasma (Heparin, EDTA) 100 µl
  - ✓ Testing Time : 10 minutes
  - ✓ Measurement : Quantitative
  - ✓ Measurement Range : 50 – 1,500 pmol/L
  - ✓ Storage Condition : 1 - 30 °C





# 02 Product Introduction

## Vcheck Feline NT-proBNP

- Key Features

- ✓ Quantitative measurement

Quantifies the degree of elevation in NT-proBNP for an accurate evaluation

- ✓ High correlation with company 'I' lab

Vcheck Feline NT-proBNP has a high correlation ( $R^2 = 0.96$ ) with company 'I' laboratory.

- ✓ A wide range of measurement

Measures up to 1500 pmol/L

- ✓ A user-friendly procedure & Fast results


Simple one-step procedure, improving user convenience, and quick results within 10 min.



# 02 Product Introduction

## Vcheck Feline NT-proBNP


- Test Procedure
  - ✓ Samples should be centrifuged and tested immediately after collection. Alternatively, refrigerate and use within 24 hours or freeze.
    - \* Degradation of NT-proBNP may occur if stored at room temperature or refrigerated for more than 24 hours, causing false negative results.



Add 100 µl of the sample to the assay diluent tube



Mix well by using a 100 µl pipette



Add the mixed sample (100 µl) into the test device

\*Stored samples should be placed at room temperature 30 min. before use.

# 02 Product Introduction

## Vcheck Feline NT-proBNP

- Reference Range

< 100 pmol/L	≥ 100 pmol/L
Normal	Abnormal Additional diagnostics are recommended


- A positive NT-proBNP test result should always be interpreted in combination and other diagnostic findings.
- **In cats with respiratory signs**, if the NT-proBNP is >270 pmol/L, CHF is the most likely cause of the clinical signs.

# 02 Product Introduction

## Vcheck Feline NT-proBNP

- Performance

Comparative evaluation of feline NT-proBNP (N=37)



Reproducibility

CV < 15%

Accuracy

Bias < 20%

# Reference

## Reference

1. Mark Oyama. Cardiac Blood Tests in Cats: Another Tool for Detection of Heart Disease. Today's Veterinary Practice. September/October 2011
2. Natalie Stilwell, MVC 2018: Advances in Feline Heart Disease Diagnosis
3. Connolly DJ, Soares Magalhaes RJ, Fuentes VL, et al. Assessment of the diagnostic accuracy of circulating natriuretic peptide concentrations to distinguish between cats with cardiac and non-cardiac causes of respiratory distress. J Vet Cardiol 2009;11(Suppl 1):S41–S50
4. K.V. Pierce, J.E. Rush, V.K. Yang, et al. Association between Survival Time and Changes in NT-proBNP in Cats Treated for Congestive Heart Failure. J Vet Intern Med. 2017 May-Jun; 31(3): 678-684.
5. Dr. Sonya G Gordon. Cardiac Education Group. October 2014
6. Fox PR, Rush JE, Reynolds CA, et al. Multicenter evaluation of plasma N-terminal pro-brain natriuretic peptide (NT-pro BNP) as a biochemical screening test for asymptomatic (occult) cardiomyopathy in the cat. J Vet Intern Med 2011; in press.
7. Connolly, DJ , et al. The effect of protease inhibition on the temporal stability of NT-proBNP in feline plasma at room temperature. J Vet Cardiol 2011;13:13–19.

Product No.	Product Name	Product Type	Packing Unit
VCF130DC	Vcheck Feline NT-proBNP	Device	5 Tests/Kit

# Thank you

BIONOTE Marketing team

Jan. 2020

