

# Product Information **Vcheck Analyzer**

BIONOTE Marketing team

July. 2021



# Vcheck Analyzer

## 01 Vcheck Analyzer Principle

- About Vcheck Analyzer
- Available Test Parameters on Vcheck
- Vcheck Fluorescent Complexes

## 02 Product Introduction

- Vcheck V200 vs Vcheck V2400
- Product Specification (V200, V2400)
- V200 Test Modes
  - Standard Test Mode
  - Read Only Mode

## 03 Product Maintenance

- Patient Information Configuration
- Test Result Management
- Calibration and Control
- Troubleshooting
- Warranty Policy

# 01 Vcheck Analyzer Principle

- About Vcheck Analyzer
- Available Test Parameters on Vcheck
- Vcheck Fluorescent Complexes

# 01 Vcheck Analyzer Principle

## About Vcheck Analyzer

**Vcheck Analyzer** is a fluorescent immunoassay veterinary POCT which enables quantitative results only within a few minutes.



**UP TO 20+ BIOMARKERS AVAILABLE**



**INTUITIVE USER INTERFACE**



**2D BARCODE AUTOMATION**



**USER FRIENDLY PROCEDURES**



**DATA CONNECTIVITY**



**UP TO 24 TESTS AT ONCE FOR V2400 ONLY**



Experience the **Outstanding Efficiency** with Vcheck Analyzer!




# 01 Vcheck Analyzer Principle

## Available Test Parameters on Vcheck


Various test parameters available to test on one analyzer

\* Biomarkers under development



### Inflammatory

CRP, SAA



### Pancreatitis

cPL, fPL



### Hormone

Cortisol, T4, TSH  
Progesterone



### Ab Titers

CPV, CDV, CAV Ab  
FPV, FHV, FCV Ab



### Coagulation

D-dimer



### Cardiac

NT-proBNP (Canine, Feline)  
TnI (Canine, Feline\*)



### Tumor


TK1\*, AFP\*

Quantitative Test | Antibody Titer Test | Infectious Test

# 01 Vcheck Analyzer Principle


## Vcheck Fluorescent Complexes

Vcheck Analyzer utilizes **FLUORESCENT**, which is more powerful than colloidal gold material



**Non-parallel Value at High Concentration**

Superior Accuracy




**Parallel Value at High Concentration**

With **Fluorescent** particle, the result value remains parallel to the actual measured concentration and doesn't shed at a high concentration.

# 01 Vcheck Analyzer Principle

## Vcheck Fluorescent Complexes

Vcheck Analyzer utilizes **EUROPIUM**, a powerful fluorescent material



With **Europium**, the excitation wavelength and emission wavelength do not overlap, ensuring you receive precise diagnostics for improved care.



## 02 Product Information

- Vcheck V200 vs Vcheck V2400
- Product Specification (V200, V2400)
- V200 Test Modes
  - Standard Test Mode
  - Read Only Mode

# 02 Product Information

## Vcheck V200 vs Vcheck V2400

Two different fluorescent immunoassay analyzers with different needs

Compact and convenient analyzer,

### V200



Animal clinics and hospitals

### USER FRIENDLY

- 1 Test at a Time
- Compact Size (200x240x205 mm)
- 7" Color Touch Screen
- 2D Barcode Technology
- E-Chart Auto Sync
- Built-in Printer

Vcheck  
Analyzer

The best way to reduce turnaround time,

### V2400



Large hospitals, laboratory

### 24 TESTS AT ONCE

- 24 Tests at Once / 70 Tests Per Hour
- Dimension (510x566x297 mm)
- 10" Color Touch Screen
- 2D Barcode Technology
- E-Chart Auto Sync
- Built-in Printer

# 02 Product Information

## Product Specification

- **Vcheck V200**
  - ✓ **Test Capacity** : 1 test at a time
  - ✓ **Power** : AC/DC adaptor
  - ✓ **Display** : 7" Color Touch Screen
  - ✓ **Pinter** : Built-in
  - ✓ **Connectivity** : HL7 v2.6(PCD-01) / POCT1-A
  - ✓ **Dimension** : 200 x 240 x 205 mm
  - ✓ **Weight** : 2.5 kg



Product No.	Product Name	Product Type	Packing Unit
VC7402EA	V200	Equipment	1 EA

# 02 Product Information

## Product Specification

- **Vcheck V2400**
  - ✓ **Test Capacity** : 24 tests at once / 70 tests per hour
  - ✓ **Power** : AC/DC adaptor
  - ✓ **Display** : 10" Color Touch Screen
  - ✓ **Pinter** : Built-in
  - ✓ **Connectivity** : HL7 v2.6(PCD-01) / POCT1-A
  - ✓ **Dimension** : 510 x 566 x 297 mm
  - ✓ **Weight** : 20.0 kg




Product No.	Product Name	Product Type	Packing Unit
VC7403EA	V2400	Equipment	1 EA

# 02 Product Information


## V200 Test Modes

- **Standard Test Mode (1)**
  - Below is an example test process for **Vcheck Standard Test Mode**
  - Refer to the specific **Instructions for Use(IFU)** for each Test Kit.

### Coding



Select [**Standard Test**] mode



Enter the **Operator ID**,  
**Patient ID**, and **Order #**

\*If the ID is not entered, the analyzer will regard the operator as a guest user.




Insert the test device  
from the desired test kit

# 02 Product Information

## V200 Test Modes

- **Standard Test Mode (2)**
  - Below is an example test process for **Vcheck Standard Test Mode**
  - Refer to the specific **Instructions for Use(IFU)** for each Test Kit.

Apply



Apply the prepared sample to the sample hole of the test device

\*If the [START] button is not pressed within 20 seconds after applying the sample, the analyzer will display an error message.

Analysis



The analyzer will automatically scan the test device and analyze the test results

\*The scanning time can vary based on the type of the test device.

Result



Once the analysis is finished, **Result** screen will be displayed

# 02 Product Information

## V200 Test Modes

- **Read Only Mode (1)**
  - Below is an example test process for **Vcheck Read Only Mode**
  - Refer to the specific **Instructions for Use(IFU)** for each Test Kit.

### Coding




Select **[Read Only]** mode



Enter the **Operator ID**,  
**Patient ID**, and **Order #**

\*If the ID is not entered, the analyzer will regard the operator as a guest user.



Apply the prepared sample to test device,  
then immediately set the timer to the  
designated processing time as  
indicated in the specific IFU

\*Add the prepared sample to additional test devices  
for the same reagent at an interval of one minute.

# 02 Product Information


## V200 Test Modes

- **Read Only Mode (2)**
  - Below is an example test process for **Vcheck Read Only Mode**
  - Refer to the specific **Instructions for Use** for each Test Kit.

Apply

Analysis

Result



When the timer goes off, insert the first test device to have applied the sample



Once the analysis is finished, **Result** screen will be displayed within 10 seconds

\*The scanning time can vary based on the type of the test device.



Eject the test device, then continue to insert the additional test devices in the order of the sample application

\*Press [**Cancel**] at any time to stop the test progress


## 03 Product Maintenance

- Patient Information Configuration
- Test Result Management
- Calibration and Control
- Troubleshooting
- Warranty Policy

# 03 Product Maintenance


## Patient Information Configuration

Easily organize patient information and the test result settings



Print out the test result with **your handwriting** on the test device

The print-out includes **reference range** to help interpretation



Easily enter the operators and patients' information using a **USB keyboard**\*

\*Material not provided.

# 03 Product Maintenance

## Test Result Management

Organize the test results automatically through the integrated data connectivity system



### Printer

Internal printer with V200, V2400



### Memory

3,000 test with V200  
5,000 test with V2400



### Connectivity

V200 and V2400 analyzer support LIS/HIS connectivity  
User can transfer the data to their information system

The screenshot shows the Vcheck software interface with a table of analysis results. The table has columns for Analyzer, Observation Date, Panel, Result, Patient(ID), Gender, Order No, Operator, and Doctor. A teal arrow points from the table to the right-hand side of the slide.

Analyzer	Observation Date	Panel	Result	Patient(ID)	Gender	Order No	Operator	Doctor
V200 #1	2018-06-19 오전 10:29	cCortisol	19.50ug/dL				Boss	
V200 #1	2018-06-19 오전 9:55	cCortisol	1.73ug/dL				Boss	
V200 #1	2017-06-28 오후 2:42	Canine T4	7.26ug/dL				Boss	
V200 #1	2017-06-28 오후 2:41	Canine T4	0.75ug/dL				Boss	
V200 #1	2017-06-28 오후 2:41	CPV Ag	Negative(53.76%)				Boss	
V200 #1	2017-06-28 오후 2:40	CPV Ag	Positive(335.5%)				Boss	
V200 #1	2017-06-28 오후 2:40	CPV Ab	Titer value 4				Boss	
V200 #1	2017-06-28 오후 2:40	CPV Ab	Titer value 3				Boss	
V200 #1	2017-06-28 오후 2:39	CPV Ab	Titer value 2				Boss	
V200 #1	2017-06-28 오후 2:39	CPV Ab	Titer value 1				Boss	
V200 #1	2017-06-28 오후 2:38	CPV Ab	Titer value 0				Boss	
V200 #1	2017-06-28 오후 2:38	cCRP	169.60mg/L				Boss	
V200 #1	2017-06-28 오후 2:37	cCRP	16.70mg/L				Boss	



yyyy/mm/dd 12:00

Operator ID:

Patient ID:  
PUPPY

Canine CRP = 12.9 mg/L  
Procedural Control = Valid

Interpretation : Normal

< 20 mg/L : Normal  
20-30 mg/L : Equivocal  
> 30 mg/L : Abnormal

Note

# 03 Product Maintenance

## Calibration and Control

Maintain Vcheck Analyzer at its best with additional supplements



### Calibration Set

- Recommended to use once per month
- 2 years use with one set
- Consists of 3 calibration devices :
  - Cal-1 for White Calibration
  - Cal-2 for UV Sensor
  - Cal-3 for RGB Sensor

> Calibration for Analyzer only



### 2D Barcode Auto Calibration

- Automatic recognition of the test items  
> Doesn't need to insert ID chip
- Temperature calibration  
> Stable result from 15 – 30 °C
- Lot variation calibration  
> Doesn't need to calibrate for every batch



### Quantitative Test Control

- Separate package
- Individual tablet pipette for each control
- Different levels of control (Low, High)
- Freeze-drying sample

# 03 Product Maintenance

## Troubleshooting


Most common troubleshooting error messages on Vcheck Analyzer



### Cause

Start button pushed too late

Already used kit




### Cause

Small volume of sample

Too much foam

Device error




### Cause

Device temperature is not between 15~30°C

Analyzer temperature is not between 15~30°C

Defective temperature sensor


Analyzer error



# 03 Product Maintenance

## Troubleshooting

Most common troubleshooting error messages on Vcheck Analyzer



### Cause

Damaged barcode




### Cause

C-line is not colored



### Cause

Calibration overdue



### Cause

Unsupported device has been used

**Thank you**  
**Any Questions?**

BIONOTE Marketing team

July. 2021

