

www.woodleyequipment.com

g-PetTM

GLUCOSE
Veterinary Monitoring System



Quick, Easy & Economical



Key Features

- Calibrated specifically for cats, dogs & horses
- Biosensor auto code technology
- Accurate Results
- Easy to use blood glucose monitor
- Vets can train owners how to home test
- Monitor cat, dog & horse diabetes
- Small sample (1.5µl) of whole blood is required
- Confirmation window to show sufficient blood volume
- Large LCD screen
- Accurate results displayed within 8 seconds
- Readable range 1.1 – 33.3 mmol/L
- 300 result automatic memory
- Adjustable alert range alarms
- mmol/L or mg/dL

*Species Specific -
Cat, Dog & Horse*

Follow us:  

Tel: +44 (0)1204 669033
Fax: +44(0)1204 669034
Email: sales@woodleyequipment.com
Web: www.woodleyequipment.com



WOODLEY
EQUIPMENT COMPANY LTD.

OLD STATION PARK BUILDINGS, ST. JOHNS STREET, HORWICH, BOLTON, LANCASHIRE, BL6 7NY, UK

Why use the g-Pet veterinary glucose monitor rather than a human glucose meter?

Human glucometers are not always accurate when measuring glucose levels in animals. This is because red blood cell size and the distribution of glucose within the blood differs from species to species. A proportion of glucose is contained within the red blood cells and the remainder in the surrounding plasma. Humans have a higher haematocrit compared to cats, dogs and horses, as well as larger red blood cells. Therefore, the distribution of glucose within plasma is proportionally much lower in humans compared to these animals, as a larger portion of glucose is contained within the red blood cells. Subsequently, these differences in glucose distribution have a bearing on the glucose result when testing animals with a meter designed for human use, as false decreases are reported in some instances.

Species Specific Calibration

The g-Pet veterinary glucose monitor is calibrated specifically for use in cats, dogs and horses. Each g-Pet meter and each pack of g-Pet strips are supplied with three species chips (cat, dog and horse) which are individually inserted into the g-Pet meter prior to testing. The chips contain all relevant calibration information to ensure an accurate result is reported for each species. Vets no longer need to persist with inaccurate human meters as there is now a low cost, veterinary specific alternative.

Testing Made Easy



Kit Includes

- Meter
- 10 Test Strips
- Control
- User Guide
- Log Book
- Carry Case
- Lancing device and 10 Lancets



Tel: +44 (0)1204 669033
Fax: +44(0)1204 669034
Email: sales@woodleyequipment.com
Web: www.woodleyequipment.com

WOODLEY
EQUIPMENT COMPANY LTD.

OLD STATION PARK BUILDINGS, ST. JOHNS STREET, HORWICH, BOLTON, LANCASHIRE, BL6 7NY, UK