Vision V5000

 **OPERATION MANUAL**

 **GUIDA UTENTE**

 **MANUAL DE INSTRUCCIONES**

 **MANUEL D'INSTRUCTIONS**

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1.0 DESCRIPTION	page 4
2.0 INTRODUCTION	page 6
3.0 UNPACKING AND ASSEMBLY	page 6
4.0 USING THE MICROSCOPE	page 7
5.0 MAINTENANCE	page 8
6.0 TECHNICAL SPECIFICATIONS	page 9
7.0 RECOVERY AND RECYCLING	page 9

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use.

We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users.

We decline any responsibility deriving from instrument uses that do not comply with this manual.

Safety guidelines

This manual contains important information and warnings regarding safety about installation, use and maintenance of the microscope. Please read this manual carefully before using the equipment. To ensure safe use, the user must read and follow all instructions in this manual. The products are designed for safe use in normal operating conditions. The equipment and accessories described in the manual are manufactured and tested according to industry standards for safety instrumentation laboratory. Misuse can cause personal injury or damage to the instrument. Keep this manual at hand close to the instrument, for an easy consultation.

Electrical safety

Before connecting the power cord to wall outlet, ensure that your mains voltage for your region corresponds to the voltage supply of the instrument, and that the illuminator's switch is in position OFF. The user must observe the safety regulations in force in his region. The instrument is equipped with CE safety marking, in any case the user has full responsibility concerning the safe use of that instrument.

Warning/Caution symbols used in this manual

The user should be aware of safety aspects when using the instrument. Warning or hazard symbols are shown below. These symbols are used in this manual.

 **DANGER** The instructions on this symbol to avoid possible severe personal injuries.

 **WARNING** Warning of use; the incorrect operation on the instrument can cause damages to the person or instrument.

 **WARNING** Possibility of electric shock.

 **HOT!** Attention: high temperature surfaces. Avoid direct contact.

 **NOTE** Technical notes or usage tips.

1.0 DESCRIPTION

REAR:

TENSION
ADJUSTMENT KNOB

ON/OFF SWITCH

POWER SUPPLY
INPUT JACK

2.0 INTRODUCTION

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use.

Woodley reminds you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users.

Woodley declines any responsibility deriving from instrument uses that do not comply with this manual.

3.0 UNPACKING AND ASSEMBLY

- 3.1 The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.
- 3.2 Place the observation head onto the top of the arm and tighten the lock-screw. Insert the eyepieces into the eye tubes.
- 3.3 Connect the provided 7,5Vdc power supply to the power supply input jack on the rear of the microscope.

Connect the mains plug into the socket at the base

WARNING

Make sure, before you turn the illumination on, that the voltage selector is set to the mains voltage for your region.

WARNING

The power cord should be used only on network sockets equipped with adequate grounding. Contact a technician to check the state of your electrical system. If there is no need to install additional accessories, the instrument is now ready for use. Once positioned and installed with the necessary components, the microscope is ready to be used. Your microscope is a laboratory instrument designed to last. Handle it always carefully and avoid abrupt vibrations or shocks. Always disconnect the power cable from the microscope when not in use for long time, while you clean it or when you perform any maintenance.

WARNING

AVOID DISASSEMBLING THE INSTRUMENT

Do not disassemble the instrument. This entails the cancellation of the warranty and may cause malfunction.

4.1 Adjust the observation head

Loosen the lock-screw, turn the observation head to a comfortable position for observation, and then lock the lock-screw.

4.2 Place the specimen on the stage

Lock the specimen slide on the mechanical stage using the slide clamp. Ensure that the specimen is centred over the stage opening by adjusting the coaxial knobs of the stage.

4.3 Lamp settings

The microscope is fitted with power white LED illuminator. Before turning on the light, read the section 5.3 about electrical safety precautions. Insert the plug of the 7,5Vdc power supply into the jack on the rear and turn on the ON/OFF switch. Turn the brightness adjustment knob to a brightness suitable for observation.

4.4 Adjust interpupillary distance

Hold the right and left parts of the observation head by both hands and adjust the interpupillary distance by turning the two parts until one circle of light can be seen.

4.5 Focus and diopter adjustment

Rotate the coarse focusing knob to bring the slide into focus using an objective with a low magnification. Adjust the fine focusing knob to get the image sharp and clear while observing with your right eye, then turn the dioptic adjustment ring on the left eyepiece until a sharp and clear image can be observed also with the other eye. Turn the tension-adjust-knob to get a suitable tension for the focus system. When the image is in focus, turn the nosepiece to choose the objective you need.

4.6 Condenser adjustment

Hold the notched ring at the base of the condenser and rotate it up or down until you reach a clear and uniform sample illumination.

4.7 Set the numerical aperture

Adjust the aperture of the iris diaphragm to set the numerical aperture of the illuminator, thus controlling image contrast and resolution.

4.8 Additional filters

The green glass filter can be inserted in the flip-out filter holder underneath the condenser.

4.9 Solar panel

If a solar panel is needed to work with microscope, please refer to the instructions enclosed with solar panel kit. Use the DC jack from back of microscope to plug the solar panel.

4.10 Fitting the mirror

If the microscope is supplied with a plano-concave mirror, this can be used with the microscope instead of the electrical light, please follow these instructions:

- 1) Make sure the equipment is not plugged into the electrical power and the switch is in the OFF position.
- 2) Unscrew the light house from the illuminator, (black part).
- 3) Screw the mirror in place of the light house.

HOLD THIS NOTCHED RING AND ROTATE IT

5.1 Microscopy environment

This microscope is recommended to be used in a clean, dry and shock free environment with a temperature of 0-40°C and a maximum relative humidity of 85 % (non condensing). Use a dehumidifier if needed.

5.2 To think about when and after using the microscope

- The microscope should always be kept vertically when moving it and be careful so that no moving parts, such as the eyepieces, fall out.
- Never mishandle or impose unnecessary force on the microscope.
- Never attempt to service the microscope yourself.
- After use, turn off the light immediately, cover the microscope with the included dust-cover, and keep it in a dry and clean place.

5.3 Electrical safety precautions

- Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.
- Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

5.4 Cleaning the optics

- If the optical parts need to be cleaned try first to: use compressed air.
- If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent.
- As a final option: use the Microscope Lens Cleaner (Product Code WD0384) available from Woodley Equipment.

Note: For regular maintenance and cleaning we recommend using our Microscope Service Kit (Product Code WD0383) which includes: Lens Cleaning Fluid (WD0384), Air Duster (WD0385), Cotton Applicators (WD0386), Microglide Lubricant (WD0387), Lint-free lens paper (WD0388) and illustrated service instructions. All these products are also available individually


- 5.5 If you need to send the microscope to Woodley for maintenance, please use the original packaging.

HEAD:	V5000: binocular, 360° rotating, 30° inclined
EYEPIECES:	WF 10X/ 18mm
NOSEPIECE:	Quadruple, reversed
OBJECTIVES:	Achromatic 4x, 10x, 40x, 100x (oil immersion)
WORKING STAGE:	Double layer with mechanical sliding stage, 125x115mm, moving range 50x30mm
CONDENSER:	1.25 N.A. Abbe type, with iris diaphragm
FOCUSING SYSTEM:	Coaxial coarse and fine focusing system, with tension adjustment mechanism
ILLUMINATION:	Extra-power White LED, non-rechargeable, with brightness control.
MAINS POWER:	Input Voltage of the power supply: 100-240Vac / 50-60Hz

7.0 RECOVERY AND RECYCLING

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste.

The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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HALTEN SIE DEN RING AM KONDENSORFUSS FEST UND DREHEN IHN AUF ODER AB
Woodley Equipment Company Ltd