

# Abbott i-STAT Portable Clinical Analyser

Cartridges Available



Granted Waived Status for lithium heparin and sodium heparin whole blood venous samples only.

## Blood Gas & Electrolytes



**i-STAT<sup>CG8+</sup>**  
Sodium (Na)  
Potassium (K)  
Ionised Calcium (Ca)  
Glucose (Glu)  
Hematocrit (Hct)  
pH  
PCO<sub>2</sub>  
PO<sub>2</sub>  
TCO<sub>2</sub>\*  
HCO<sub>3</sub>\*  
BE<sub>ecf</sub>  
sO<sub>2</sub>\*  
Haemoglobin\* (Hgb)



**i-STAT<sup>EG7+</sup>**  
Sodium (Na)  
Potassium (K)  
Ionised Calcium (iCa)  
Haematocrit (Hct)  
pH  
PCO<sub>2</sub>  
PO<sub>2</sub>  
TCO<sub>2</sub>\*  
HCO<sub>3</sub>\*  
BE<sub>ecf</sub>  
sO<sub>2</sub>\*  
Haemoglobin\* (Hgb)



**i-STAT<sup>6+</sup>**  
Sodium (Na)  
Potassium (K)  
Chloride (Cl)  
Urea Nitrogen (BUN)/Urea  
Glucose (Glu)  
Haematocrit (Hct)  
Haemoglobin\* (Hgb)



**i-STAT<sup>EG8+</sup>**  
Sodium (Na)  
Potassium (K)  
Chloride (Cl)  
pH  
PCO<sub>2</sub>  
Urea Nitrogen (BUN)/Urea  
Glucose (Glu)  
Haematocrit (Hct)  
TCO<sub>2</sub>\*  
HCO<sub>3</sub>\*  
BE<sub>ecf</sub>  
Anoin Gap\* (Agap)  
Haemoglobin (Hgb)



**i-STAT<sup>EG6+</sup>**  
Sodium (Na)  
Potassium (K)  
Haematocrit (Hct)  
pH  
PCO<sub>2</sub>  
PO<sub>2</sub>  
TCO<sub>2</sub>\*  
HCO<sub>3</sub>\*  
BE<sub>ecf</sub>  
sO<sub>2</sub>\*  
Haemoglobin\* (Hgb)



**i-STAT<sup>EC4+</sup>**  
Sodium (Na)  
Potassium (K)  
Glucose (Glu)  
Haematocrit (Hct)  
Haemoglobin\* (Hgb)



**i-STAT<sup>E3+</sup>**  
Sodium (Na)  
Potassium (K)  
Haematocrit (Hct)  
Haemoglobin\* (Hgb)

## Blood Gas



**i-STAT<sup>G3+</sup>**  
pH  
PCO<sub>2</sub>  
PO<sub>2</sub>  
TCO<sub>2</sub>\*  
HCO<sub>3</sub>\*  
BE<sub>ecf</sub>  
sO<sub>2</sub>\*



**i-STAT<sup>CG4+</sup>**  
pH  
PCO<sub>2</sub>  
PO<sub>2</sub>  
TCO<sub>2</sub>\*  
HCO<sub>3</sub>\*  
BE<sub>ecf</sub>  
sO<sub>2</sub>\*  
Lactate

## Chemistry



**i-STAT<sup>Crea</sup>**  
Creatine (Crea)



**i-STAT<sup>G</sup>**  
Glucose (Glu)



**CHEM8+**  
Sodium (Na)  
Potassium (K)  
Chloride (Cl)  
Ionised Calcium (iCa)  
TCO<sub>2</sub>  
Glucose (Glu)  
Urea Nitrogen (BUN)/Urea  
Creatine (Crea)  
Haematocrit (Hct)  
Haemoglobin (Hgb)

## Coagulation



**i-STAT<sup>Celite ACT</sup>**  
Celite ACT



**i-STAT<sup>PT/INR</sup>**  
Prothrombin Time



**i-STAT<sup>Kaolin ACT</sup>**  
Kaolin ACT

## Cardiac Markers



**cTnl**  
Troponin I



**CK-MB**  
CK-MB



**BNP**  
BNP



\* Calculated

Tel: +44 (0) 8456 777001 | Fax: +44 (0) 8456 777002  
Email: sales@woodleyequipment.com  
www.woodleyequipment.com

**WOODLEY**  
EQUIPMENT COMPANY LTD.