

Change To QBC Product Code 424271

Please be informed that there has been a change to QBC Product Code 424271. The fibre optic cable, which is part of the Paralens system, will no longer be supplied with the black adaptor as per the images below.

Image 1 shows the configuration of the product previously available.



Image 2 shows the new configuration of the product.



Your customers will need to retain the black adaptor from their existing fibre optic cable to ensure the new fibre optic cable is compatible with their light source.

If the adaptor has been lost, the new fibre optic cable will not be compatible and a new light source will need to be purchased along with a new cable.