

www.woodleyequipment.com

Ready Whenever and Wherever
You Need It.

Introducing the New and Improved
i-STAT 1 Portable Critical Care Analyser.

Updated
& upgraded.

Delivering
more tests,
more cartridges
& more
capabilities.



**Clinical Chemistry, Electrolyte &
Blood Gas Analysis.**
Patient Side or in the Field.



Tel: +44 (0)1204 669033

Fax: +44 (0)1204 669034

Email: sales@woodleyequipment.com

Web www.woodleyequipment.com

WOODLEY
EQUIPMENT COMPANY LTD.

OLD STATION PARK BUILDINGS . ST. JOHNS STREET . HORWICH . BOLTON . LANCASHIRE . BL6 7NY

Accurate Results Within Easy Reach!

ANALYSER TEST MENU		CG4+	CG8+	G	Glu*	Crea	E3+	6+	CHEM 8+*	EC8+	ACT Celite	cTnl*
Haematology	Haematocrit (Hct)		X				X	X	X	X		
	Haemoglobin (Hgb)		X				X	X	X	X		
Chemistry	Blood Urea Nitrogen (BUN)							X	X	X		
	Creatinine					X			X			
	Ionized Calcium (iCa)		X						X			
	Glucose (Glu)		X	X	X			X	X	X		
Electrolytes	Chloride (Cl)						X	X	X	X		
	Sodium (Na)		X				X	X	X	X		
	Potassium (K)		X				X	X	X	X		
Acid Base	pH	X	X							X		
	PCO ₂	X	X							X		
	HCO ₃	X	X							X		
	TCO ₂	X	X						X	X		
	Anion Gap								X	X		
	Base Excess	X	X							X		
Blood Gas	PCO ₂	X	X									
	sCO ₂	X	X									
Extras	Lactate	X										
	ACT (Celite)										X	
	Troponin I (cTnl)											X

* = only available on the i-STAT 1 analyser

** = Chloride on the EC3 cartridge is available on the original i-STAT analyser only.

Features:

- ◇ Handheld
- ◇ Results in 2 minute
- ◇ Auto calibration and QC
- ◇ Laser barcode reader
- ◇ 5000 test storage ability
- ◇ Glucose strip reader
- ◇ Larger screen
- ◇ Cardiac Troponin I test
- ◇ Quick and accurate
- ◇ New ergonomic design
- ◇ Wide test menu
- ◇ CHEM8+ cartridge.

Tel: +44 (0)1204 669033

Fax: +44 (0)1204 669034

Email: sales@woodleyequipment.com

Web www.woodleyequipment.com

WOODLEY
EQUIPMENT COMPANY LTD.