

www.woodleyequipment.com

Ready Whenever and Wherever You Need It.

Introducing the New
and Improved
**i-STAT 1 Portable
Critical Care Analyser.**

**Updated
and Upgraded.**

Delivering more tests,
more cartridges
& more capabilities.

**Clinical Chemistry,
Electrolyte & Blood Gas
Analysis.**

Patient Side or in the Field.



Tel: +44 (0)1204 669033

Fax: +44 (0)1204 669034

Email: sales@woodleyequipment.com

Web www.woodleyequipment.com

**WOODLEY
EQUIPMENT COMPANY LTD.**

OLD STATION PARK BUILDINGS • ST. JOHNS STREET • HORWICH • BOLTON • LANCASHIRE • BL6 7NY.

Accurate Results Within Easy Reach!

ANALYSER TEST MENU

		CG4+	CG8+	G	Glu*	Crea	E3+	6+	CHEM 8+ *	EC8+	ACT Celite	cTnl *
Haematology	Haematocrit (Hct)		X				X	X	X	X		
	Haemoglobin (Hgb)		X				X	X	X	X		
Chemistry	Blood Urea Nitrogen (BUN)							X	X	X		
	Creatinine					X			X			
	Ionized Calcium (iCa)		X						X			
	Glucose (Glu)		X	X	X			X	X	X		
Electrolytes	Chloride (Cl)						X	X	X	X		
	Sodium (Na)		X				X	X	X	X		
	Potassium (K)		X				X	X	X	X		
Acid Base	pH	X	X							X		
	PCO ₂	X	X							X		
	HCO ₃	X	X							X		
	TCO ₂	X	X						X	X		
	Anion Gap								X	X		
	Base Excess	X	X							X		
Blood Gas	pO ₂	X	X									
	sO ₂	X	X									
Extras	Lactate	X										
	ACT (Celite)										X	
	Troponin I (cTnl)											X

* = only available on the i-STAT 1 analyser

** = Chloride on the E3+ cartridge is available on the original i-STAT analyser only.

Features:

- Handheld
- Results in 2 minutes
- Auto calibration and QC
- Laser barcode reader
- 5000 test storage ability
- Glucose strip reader
- Larger screen
- Cardiac Troponin test
- Quick and accurate
- New ergonomic design
- Wide test menu
- CHEM8+ cartridge

Tel: +44 (0)1204 669033

Fax: +44 (0)1204 669034

Email: sales@woodleyequipment.com

Web www.woodleyequipment.com

WOODLEY
EQUIPMENT COMPANY LTD.

OLD STATION PARK BUILDINGS • ST. JOHNS STREET • HORWICH • BOLTON • LANCASHIRE • BL6 7NY.