WOODLEY EQUIPMENT COMPANY

PRODUCT GUIDE

Issue 3.0

















"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"





WELCOME TO OUR PRODUCT GUIDE

Woodley Equipment Company was **established in 1989** and prides itself on being a market leading supplier of specialist **laboratory equipment**, **diagnostic and critical care** products to the **veterinary industry worldwide**.

A dedicated on-the road Account Manager team provides **in-practice support**, training and on-going advice. This team is supported at our headquarters with a technical support and service department staffed with **Biomedical Scientists** and **Service Engineers** who deliver proven high quality after-sales care.

As well as a **UK nationwide sales force**, Woodley has an international distribution network offering a rapidly growing portfolio of **high quality products** with worldwide distribution rights. Woodley also has a proven track record and outstanding **customer loyalty** in the form of repeat business and client referrals.

Woodley seeks to supply veterinarians and distributors with the equipment and services to best support Practices, thus enabling them to respond to the needs of their clients, while operating a more **efficient and profitable business.**

Quantum Vet Diagnostics, a division of Woodley Equipment, is the **UK's number 1** provider of veterinary wet chemistry solutions.



"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"









Insight V3®PLUS

Veterinary Haematology Analyser

A fully automated multi-species • Quick and easy to use **Haematology Analyser**

3 PART DIFFERENTIAL WITH EOSINOPHIL FLAG.

Parameters Available

WBC	RBC	RDW_SD
LYM%	HGB	PLT
MID%	HCT	MPV
GRAN%	MCV	PDW
LYM#	MCH	PCT
MID#	MCHC	P_LCR
GRAN#	RDW_CV	P_LCC

- Small sample volume 13μl whole blood
- 60 second analysis time
- 3-part differential, 21 parameters with 3 histograms and eosinophil flag
- Species specific reference ranges
- Impedance technology
- Stores up to 100,000 results
- Fully automated
- Touch screen
- Low maintenance
- Results transmitted via InSight p-LINK to Practice Management Systems



InSight V5°

Veterinary Haematology Analyser

A fully automated multi-species • Quick and easy to use **Haematology Analyser**

5 PART DIFFERENTIAL WITH LASER TECHNOLOGY.

Parameters Available

WBC	LYM#	HGB	HCT
NEU%	EOS%	MCV	PLT
NEU#	EOS#	MCH	MPV
MON%	BASO%	MCHC	PDW
MON#	BASO#	RDW-CV	PCT
LVN/10/	DDC	DD/W/ CD	

- Small sample volume 15μl whole blood
- 60 second analysis time
- 5-part differential, 23 parameters
- Species specific reference ranges
- Independent channel for Basophil measurement
- Powerful capability to flag abnormal cells
- Fully automated
- Touch screen
- Low maintenance
- Results transmitted via InSight p-LINK to Practice Management Systems



Thisight 5 Diff Retic

THSIRM'S Diff Retio

nSight 5 Diff Retic

5 Part Differential Laser Haematology DUAL LASER

IMPEDANCE TECHNOLOGY. **Parameters Available**

White Blood Cells Red Blood Cells

		rica biood cons		
WBC		RBC	Retics#*	
NE#	NE%	Hb	Retics%*	
LY#	LY%	HCT	RDW	
MO#	MO%	MCV	RSD	
EO#	EO%	MCH	NRBC#*	
BA#	BA%	MCHC	NRBC%*	
BAND#	BAND%	MCVr	RDWr	

Platelets Printed WBC, RBC, PIT PDW **PLT Cyotgrams** MPV PCT *certain species only

Analyser with Reticulocytes • Compact, benchtop size

- Small sample size 20µl whole blood
- True 31 parameter system, 2 minute analysis time
- Reticulocytes on every sample
- Quick and easy operation, low maintenance
- Over 80 species on board
- 5-part differential with additional reticulocytes, band neutrophils & NRBC
- 2 new Retic measurements for maturity index and anisocytosis
- Touch screen
- Results transmitted via InSight p-LINK to Practice Management Systems



To facilitate connectivity and results integration with Practice Management Systems

- InSight p-LINK is a small computer which enables sample results from in-practice analysers to be transferred to a Vet Practice network folder with results in Text, CSV, PDF or VetXML formats
- Now includes cloud based access to lab reports
- Quick and easy to set up via a wired LAN point or Wi-Fi
- Works in the background no user interaction required

RapidVet®- H Blood Typing Kits - Canine & Feline

In-practice tests for rapid determination of Canine DEA 1.1 and Feline A, B or AB blood types.

These fast in-house blood typing card kits for Dogs and Cats are essential for transfusion and breeding decisions. They are fast, easy to use, economical and can be stored at room temperature.


RapidVet-H Canine

- For Canine DEA 1.1
- Fast, easy, economical and VITAL
- Results in 2 minutes
- Auto-agglutination screen
- Requires 3 drops of whole blood
- Available in kits of 5



RapidVet-H Feline

- For Feline A,B & AB
- Fast, easy, economical and VITAL
- Results in 2 minutes
- Auto-agglutination screen
- Requires 3 drops of whole blood
- Available in kits of 5



"Every Dog & Cat in need of a transfusion should be blood typed" Dr. Urs Giger, Dr.med.vet., Dip ACVIM, ECVIM

RapidVet®- H Canine DAL

Canine blood types other than DEA 1 can also be important in transfusion decisions. RapidVet-H Canine DAL is intended for use to classify dogs as DAL positive or negative. This information can be essential in donor selection, especially in the Dalmation, Doberman, Shih Tzu and Lhasa Apso breeds, which reportedly have a higher percentage of DAL-negative blood types. It is also useful in determining the cause of an incompatible crossmatch between DEA 1 compatible dogs.

As many as 98% of dogs, including blood donors, carry the DAL antigen. This antigen elicits strong immunogenicity. Sensitisation through transfusion can occur in dogs that lack it.

DAL-negative blood types are more common in certain breeds and they are more at risk of delayed and acute haemolytic transfusion reactions if transfused with DAL positive blood. Doberman Pinschers are especially at risk since the prevalence of both DAL-negative blood and von Willebrand disease is relatively high.



- Enable more complete, in-practice, pre-transfusion work-up and donor selection
- Useful in determining possible cause of an incompatible crossmatch between DEA 1 compatible dogs
- Clear, easy to follow instructions
- Visible results within 1-2 minutes
- 2 or 5 card kits available

RapidVet®- H Crossmatch Kits

In-practice crossmatch test kits for use in tranfusion medicine on either Canine or Feline patients.

- Will detect all major Canine and Feline blood group incompatibilities
- Performed in-practice in 10 minutes
- · Easy to use and results are clear
- Available in kits of 3 tests
- Minor Crossmatch Kits also available

Essential for good transfusion practice



Use alongside Canine & Feline Blood Typing Kits

InSight® Roller Mixer

- Gentle rocking and rolling motion for optimum mixing
- Fixed speed of 30 RPM with speed accuracy of +/-2 RPM
- Continuous/infinite mode of operation
- Tilt angle of +/- 3°
- Tool free removal/installation of rollers for replacement and cleaning
- Adjustable rollers spacing to hold different diameter tubes
- Suitable for mixing of blood samples, viscous substances and for liquid-solid suspensions

InSight® Veterinary Blood Urea Test Strips

A quick and easy, low cost test for in-practice quantitative determination of urea in whole blood.

- One drop of venous or capillary whole blood
- Test time only 90 seconds
- Colour chart for reading test results is printed on the bottle label



InSight qLabs Veterinary Coagulation Analyser

- Requires 1 drop of non-anticoagulated whole blood (10µl)
- PT/aPTT tested simultaneously on 1 test strip
- Room temperature test strip storage
- Fill test strip directly from syringe
- QC material available
- No need to transfer sample into citrate tube
- Low cost analyser and consumables
- 'Test count up' for immediate assessment
- 2 year warranty

In-practice Point of Care coagulation testing

(InSight qLabs

HemoCue® WBC Analyser

HemoCue WBC provides a total white cell count. The system offers an innovative, quick and reliable way to assess the health of the animal.

- Easy to use
- Results in 2 minutes
- Small sample volume 10μl
- WBC count available using whole blood and synovial fluid



Simple to use & results in 2 minutes

Partner with InSight HCT Meter for **Rapid Haematological Assessment**

InSight® Microscope Eyepiece **Adapter for Smartphones**



- Mounts any smartphone to any microscope in minutes!
- Smartphone camera can then capture photos & videos for real-time viewing on smartphone screen, via video conferencing, emailing or sharing online

Capture & share microscopic images



The InSight HCT Meter is a low cost handheld meter for the measurement of Canine, Feline and Equine haematocrit and total haemoglobin in capillary and venous whole blood.

- Calibrated for Dogs, Cats and Horses
- 30 second test time
- One drop of whole blood (Venous, Capillary or EDTA)
- HCT results reported as %



Calibrated and validated specifically for Dogs, Cats and Horses

Vision V5000 Microscope

- LED illuminator
- Excellent optical resolution
- Wide-field eyepieces
- Dioptric compensation
- DIN achomatic objectives 4x, 10x, 40x, 100x
- Coaxial coarse and fine focusing



Haematology Atlas

Haematology Atlas written exclusively for Woodley Equipment by Kostas Papasouliotis DVM PhD DipRCPath MRCVS (Bristol Vet School).



Titles Available - Canine, Feline & Equine These atlases are an excellent microscope-side reference guide offering normal cell morphology.

- Erythrocytes Neutrophils
- Basophils
- Lymphocytes
- Platelets
- Methodology: Microhaematocrit
- Eosinophils
- Methodology: Blood Smears

InSight Mini Chem°

The InSight Mini Chem is intended for the quantitative determination of Urea (UR), Creatinine (CR) and Uric Acid (UA) in capillary whole blood, venous whole blood, plasma (lithium heparin) and serum.

- Easy to operate
- Portable meter
- Quick and accurate results in under 5 minutes
- Multi function 3 parameters on 1 reagent strip
- Bluetooth built-in to enable transfer of data to your phone, tablet or PC
- Rechargeable built-in Lithium-ion battery
- Small sample requirement 35µl of whole blood, plasma or serum
- Stores up to 500 results



3 Parameters on 1 Reagent Strip





A fully automated clinical chemistry analyser which uses liquid reagents and is approximately 70% cheaper to run than dry chemistry

Chemistry Analyser

- User friendly intuitive veterinary software
- Onboard refrigerated reagent/sample compartment
- Minimal user maintenance
- Fast throughput example batch times: 9 mins 1x GHP, 13 mins 2x GHP, 16 mins 3x GHP
- Continuous loading, no need to wait for batch completion
- Stat function for urgent tests
- Fully automated with auto-reflex dilution of over range results
- Flexible testing, single parameter or up to 20 user defined profiles
- Small footprint 414 (W) x 565 (D) mm
- Online remote assistance directly to your analyser
- Long life ceramic syringes to ensure accuracy and precision
- Environmentally friendly semi-disposable reaction cuvettes with high efficiency wash system

Customer defined profiles and/or tests can be programmed from the following:

Albumin, Alkaline Phosphatase, ALT, AST, Amylase, Bilirubin (Total), Calcium, Cholesterol, CK, Creatinine, Glucose, GGT, Phosphate, Protein (Total) & Urea.

BHB, Bile Acids, Chloride, Copper, Fibrinogen, Fructosamine, GLDH, GSH-Px (as Selenium), Haemoglobin, IgG, LDH, Lipase, Magnesium, NEFA, SAA, T4 (Total Thyroxin), Triglycerides, Uric Acid, Urine Creatinine, Zinc, Phenobarbitone, Iron, Ferritin and TIBC



Skyla VB1 Veterinary Chemistry Analyser



- Compact, benchtop size
- Performs up to 15 biochemistry tests simultaneously
- Performs up to 3 single/dual tests simultaneously
- 200µl whole blood, serum or plasma
- LCD colour touch screen & built-in printer
- Fully automatic system with walk away operation
- Built-in QC and calibration for each test
- Species specific reference ranges
- Stores up to 50,000 results

Parameters Available:

ALB, ALP, ALT, AMY, AST, BUN, Ca, CREA, CPK, GGT, GLU, K+, Na+, PHOS, TBIL, CHOL, TP, B Acids, T4, Mg, UA, TCO2, TRIG, LAC, FRU, LIPA, UCRE, UPRO, GLOB*, A/G*, Na+/K+, BUN/CREA*, UPC*, UREA* * Calculated

For both Rotor Profiles & Single/Dual Tests

Highly accurate test results in under 15 minutes

skyla VB1 & Solution **Panel Test Portfolio:**

- Preanaesthetic Liver Plus
- Renal
- Avian & Reptile
- Diabetes
- Diagnosis Diagnosis Plus

- Electrolytes
- Equine
- Diagnosis II
- Thyroid/T4
- Critical Care
- Large Animal



skyla

Carly la

Veterinary Chemistry + Immunoassay Analyser



- Can run 3 rotors simultaneously, 2 Chemistry + 1 Immunoassay
- Chemistry rotors and parameters identical to skyla VB1⁺ Analyser
- Immunoassay parameters currently available: cCRP, SAA, PHBR, cPROG, TT4, cTSH, cCOR, cA1C
- Easy and quick operation
- Easy to use LCD colour touch screen
- Stores up to 50,000 results

Single/Dual Chemistry Parameters for skyla VB1 & skyla Solution

Single Test

- CREA
- BUN
- TBIL
- LAC
- Ammonia
- **Dual Test** Na⁺ & K⁺
- PHOS & Ca
- CHOL & TRIG
- GGT & CPK
- ALT & ALP
- BA & CHOL
- AST & CPK
- LIPA & AMY
- UPC Kit
- (UPRO & UCRE)



Veterinary Glucose Meter

g-Pet PLUS Veterinary Glucose Meter is calibrated specifically for use in Dogs, Cats and Horses, ensuring an accurate result is reported for each species.

- Biosensor auto code technology
- Accurate results
- Easy to use
- Small sample of whole blood (0.7µl)
- Accurate results displayed within 5 seconds
- Readable range 1.1 33.3 mmol/L
- 500 result automatic memory
- HCT range 20-60%
- Lifetime warranty



Calibrated for Canine, Feline & **Equine use**

Pet

Veterinary Lactate Meter

L-Pet only requires a small drop of blood (3µl). The test strips draw blood in and then measure using an advanced biosensor specific for lactate concentration in blood.

- Advanced Biosensor Test Strips
- Small sample size (3µl)
- Results in less than 45 seconds
- Large LCD display
- Stores up to 300 results
- Easy operation
- Portable
- Specific Lactate Analyser QC Material available
- 5 year warranty



An easy to use handheld meter providing blood lactate results in just 45 seconds

i-SMART 30 VET ELECTROLYTE ANALYSER

Fast & Accurate

- Fast, accurate analysis
- Tests electrolytes Na+, K+, Cl-& iCa++
- 60μl sample volume
- Results in 35 seconds
- Species specific reference ranges

All-in-One Cartridge

- Cartridge contains electrodes, all consumables & waste
- No user maintenance
- Automatic calibration



Measured electrolytes

20 - 250 mmol/L

0.5 - 20.0 mmol/L

Cl-20 - 250 mmol/L

iCa++ 0.25 - 5.00 mmol/L

Fast, accurate & maintenance free In-Practice electrolyte analysis - now with Ionised Calcium

PocketChem™ BA Blood Ammonia

Analyser

PocketChem BA measures blood ammonia levels. The device enables immediate testing and delivers fast, accurate results.

- Whole blood measurement
- Ready for use 24 hours a day
- Small 20μl sample volume
- Battery powered
- Easy to maintain
- Temperature compensation
- Lightweight and compact design
- Detachable printer available
- Supports the diagnosis of hepatic cephalopathy
- Helps detect and evaluate urea cycle defects
- Monitor patients with liver disease



- Accurate, quantitative results in 13 seconds
- Easy to use and portable
- Small 0.6µl sample size
- · Venous or arterial whole blood
- Correct results for HCT
- No calibration required
- Eliminates interfering substances
- Stores 400 results

StatStrip® Creatinine Xpress Meter

- Rapid assessment of kidney function
- Accurate, quantitative results in 30 seconds
- Easy to use and portable
- Small 1.2µl sample size
- Capillary, venous or arterial whole blood
- Correct results for HCT
- No calibration required
- Eliminates interfering substances
- Stores 400 results

NovaVet Blood Ketone Meter

- Measures blood beta-hydroxybutrate (BHBA), the gold standard for detecting subclinical ketosis
- · Accurate, quantitative results in 10 seconds
- Easy to use and portable
- No calibration required
- Simple test procedure in the field
- Small 0.8µl sample size
- Stores 400 results



021 BS

$In Sight^{\circ}$ Vital Signs + ECG

InSight Vital Signs + ECG Handheld Monitor is ideal for use in Veterinary Critical Care, General Health Screening and Anaesthesia Monitoring.

- Small handheld monitor
- Measures ECG + NIBP + SpO₃ + TEMP + Heart Rate + Mean Arterial Pressure (MAP)
- Includes cage clip & charging dock
- Easy to use
- Accurate measurement in 0.05% **Low Perfusion Index**
- Audio and visual alarm
- Records up to 48 hours of measurement data

Simple, Accurate, Low Cost, **Vital Signs Monitoring**

Measures ECG + NIBP + SpO₂ + TEMP + Heart Rate + Mean Arterial Pressure (MAP)



InSight Blood Pressure **Monitor**

Portable veterinary blood pressure monitor suitable for spot checks or continuous monitoring.

- Accurate and reliable uses 2 patented, algorithms specifically developed for animals based on animal vital signs
- Easy to use large buttons or intuitive touch screen
- Silent measures BP quietly and alarms can be silenced to avoid distressing animals
- Fast take multiple measurements quickly
- Memory 1000 results in Monitor Mode, 100 patients and 200 Data sets per patient in Spot Mode
- Portable rechargeable lithium or alkaline battery operation
- Data transfer APP available to transfer data to mobile phone or Tablet/PC via Blutooth

Measures NIBP



InSight® Vet Vital Signs

InSight Vet Vital Signs Handheld Monitor is ideal for use in Veterinary Critical Care, General Health Screening and Anaesthesia Monitoring.

- Small handheld monitor
- Measures NIBP + SpO₂ +TEMP + Heart Rate + Mean Arterial Pressure (MAP)
- Includes cage clip & charging dock
- Easy to use
- Accurate measurement in 0.05% **Low Perfusion Index**
- Audio and visual alarm

Simple, Accurate, Low Cost,

• Records up to 48 hours of measurement data

Vital Signs Monitoring

Measures NIBP + SpO₂ + TEMP + Heart Rate + Mean Arterial Pressure (MAP)



InSight Vital Signs + ECG Handheld Monitor is suitable for small animal vital signs monitoring of blood oxygen saturation (SpO₂).

- Small handheld monitor
- Measures SpO₃
- Includes cage clip & charging dock
- Easy to use
- Accurate measurement in 0.05% Low Perfusion Index
- Audio and visual alarm
- Records up to 48 hours of measurement data

Simple, Accurate, Low Cost, **Vital Signs Monitoring**

Measures SpO₃



Veterinary iPhone ECG/Heart Monitor Universal Adaptor is designed to fix onto any iPhone® and wirelessly communicates with the AliveECG Vet app. Once partnered with the ECG monitor, it transforms any iPhone® into a clinical-quality, single-lead ECG recorder.

- Easy to use No stress or shaving required in most animals
- Powerful Leverage the power, display and communication capabilities of your iPhone/Smartphone
- Accurate Veterinary testing demonstrates accuracy in rate and rhythm compared to traditional devices
- Accessible Securely store ECG notes and patient information on the device and in the cloud
- Convenient Screen a wide variety of animals quicker and easier during wellness and sickness exams, in addition to pre-operative and diagnostic work-ups

An ECG device as mobile as you are

ImmunoComb[®] VacciCheck[®] In-Practice Antibody Test Kits

ImmunoComb VacciCheck Antibody Test Kit monitors serological status, reduces vaccination failure and prevents over vaccination consequences for Dogs and Cats.

VacciCheck is recommended by the WSAVA to comply with regulations to reduce unnecessary vaccinations

VACCICheck Canine

The Canine VacciCheck Antibody Test Kit is designed to determine Dog serum antibody titre to:

- Infectious Canine Hepatitis (ICH)
- Canine Parvovirus (CPV)
- Canine Distemper Virus (CDV)

The main purpose of this kit is to provide a useful tool to assess the immune status of Dogs to the three core Canine infections.



VACCICheck Feline

Quick & easy to use


The Feline VacciCheck Antibody Test Kit is designed to determine Cat serum IgG antibody titre to:

- Feline Panleukopenia (FPLV)
- Feline Herpes Virus (FHV)
- Feline Calici Virus (FCV)

The main purpose of this kit is to provide a useful tool to assess the immune status of Cats to the three core Feline infections.

WSAVA Guidelines state that titre testing for the 3 core Canine and Feline pathogens should be performed instead of automatic annual revaccination





Fully supported

Nationwide by our

Specialist Team of

Field Service Technicians



Take a 360° look at your endocrinology testing with the TOSOH AIA-360



AIA-360 Analyser

Offers the automation powers and performance of a large analyser in a benchtop format. Random access allows any combination of tests, in any order at any time and a throughput of 36 results per hour with the first test result in 20 minutes ensures a fast reliable service.

- Simple to use
- Fast, reliable results
- Simple maintenance
- Touchscreen with intuitive menus

Key Assays Available:

- Free T4
- ACTH
- Total T4
- Insulin
- Cortisol
- cTnl
- Progesterone

Integrates with Practice Management Systems



Immunofluorescence Analyser

Simplify In-Practice Point of Care Testing with Reference Lab Quality - Avoiding Test Referrals

Achieve true analysis with fluorescence technology

- Fully quantitative automated analyser
- Accurate diagnosis for Dogs and Cats in minutes
- Analyser provides print out and result storage



Canine Tests Available:

- cPL
- T4
- TSH
- Cortisol
- CRP
- Progesterone
- D-Dimer
- NT-proBNP
- Troponin I
- Core Antibody: Parvovirus, Distemper & Hepatitis

Feline Tests Available:

- T4
- fPL
- SAA
- NT-proBNP
- Troponin I
- Core Antibody:
 Calicivirus, Herpes
 Parvovirus

Specifications

Testing Time: 5 minutes (Based on CRP, SAA)

Weight: 2.5kg

Screen: Touchscreen 7" Colour TFT

Size: 20 x 24 x 20 cm Printer: Built-in Printer

Sample Required: Serum or Plasma (Heparin)

Low Cost, In-Practice, Reference Lab Quality Immunoassay Testing





Anigen™ Rapid Test Kits







Feline Tests Include:

Anigen Rapid FIV Ab/FeLV Ag Feline Immunodeficiency virus antibody and Feline Leukemia virus antigen in Feline whole blood, plasma or serum Anigen Rapid FeliD-3 Ag

Feline Parvo, Corona and Giardia antigen in Feline faeces

Anigen Rapid FCoV Ag/Ab

Feline Corona virus antigen/antibody in Feline faeces (Ag) and whole blood, plasma or serum (Ab)

High Sensitivity, High Specificity Rapid Tests Kits

Anigen Rapid FIV Ab

Immunodeficiency virus antibody in Feline whole blood, plasma or serum

Anigen Rapid FCoV Ab

Feline Corona virus antibody in whole blood, plasma or serum

Anigen Rapid FPV Ag

Feline Parvo virus antigen in Feline faeces

Anigen Rapid FHW Ag/Ab

Feline heartworm antigen and antibody in Feline whole blood, plasma or serum

Anigen Rapid FeLV Ag

Feline Leukemia virus antigen in Feline whole blood, plasma or serum

Anigen Rapid FCoV Ag

For the detection of Feline Corona virus antigen in Feline faeces

Anigen Rapid Toxoplasma Ab

Toxoplasma gondii antibody in Feline whole blood, plasma or serum

Anigen Rapid FeliCheck-3

Feline immunodeficiency virus antibody, feline leukemia virus antigen and feline heartworm antibody in Feline whole blood, plasma or serum

Available exclusively from Woodley Canine Tests Include:

Anigen Rapid CPV/CCV/Giardia Ag Canine Parvo, Corona and Giardia antigen in Canine faeces

Anigen Rapid CPV Ag

Canine Parvo virus antigen in Canine faeces

Anigen Rapid CCV Ag

Canine Corona virus antigen in Canine faeces

Anigen Rapid CIRD-3 Ag

Canine Infectious Respiratory Disease (CDV, CAV and CIV) in conjunctiva (CDV Ag) and nasal fluid (CAV and CIV Ag)

Anigen Rapid Rabies Ag

Rabies virus antigen

Anigen Rapid CPV Ab

For the titration of Canine Parvo virus in whole blood, plasma or serum

Anigen Rapid CDV Ab

For the titration of Canine Distemper virus antibody in whole blood, plasma or serum

Anigen Rapid CHW Ag

Canine Dirofilaria immitis antigen in Canine serum, plasma or whole blood

Anigen Rapid CDV Ag

Canine Distemper virus antigen in conjunctiva (recommended), urine, serum or plasma

Anigen Rapid CDV Ag/CAV Ag

Canine Distemper virus and Adeno virus antigen in conjunctiva (CDV Ag) and nasal fluid (CAV Ag)

Anigen Rapid C.Brucella Ab

Canine Brucella antibody in whole blood, plasma or serum

Anigen Rapid Leishmania Ab

Leishmania imfantum antibody in whole blood, plasma or serum

Anigen Rapid E.Canis Ab

Ehrlicia Canis antibody in whole blood, plasma or serum

Anigen Rapid Anaplasma Ab

Anaplasma phagocytophilum/ Anaplasma platys antibody in whole blood plasma or serum

Anigen Rapid Lyme Ab

Borrelia burgdorferi in Canine whole blood, plasma or serum

Anigen Rapid CPV Ag/CCV Ag

Canine Parvo virus and Corona virus antigen in Canine faeces

Anigen Rapid Giardia Ag

Giardia antigen in Canine faeces

Anigen Rapid CaniV-4

Dirofilaria immitis antigen, Ehrlichia canis antibody, Leishmania imfantum antibody, Anaplasma, phagocytophilum/ Anaplasma platys antibody in whole blood, plasma or serum

Anigen Rapid E.canis/

Anaplasma Ab

Ehrlichia canis and Anaplasma phagocytophilum/Anaplasma platys in whole blood, plasma or serum

Anigen Rapid Total IgE

Canine Total IgE in serum

Anigen Rapid Leishmania Ab

Leishmania infantum antibodies in Canine whole blood, plasma or serum

This list is not exhaustive - please contact us with your requirements



epoc Critical Care/Blood Gas Analyser

A global leader in handheld Critical Care/Blood Gas!

Test Card Contains All Parameters

Measured Analytes:

pH, pCO₂, pO₂, Na⁺, K⁺, Cl⁻, Ca⁺⁺, Glu, Lac, Hct, BUN, Creatinine, tCO₃

Calculated Values:

HCO₃, BE (ecf), BE (b), sO₃, Hgb, AGapK

Features

- Accurate blood gas, electrolyte and critical care biochemistry results in 3 minutes
- Room temperature card storage
- 92µl sample, fresh whole blood (arterial, venous or capillary)
- Veterinary species specific software
- Portable handheld system with rechargeable battery & Bluetooth printer
- Ready to use anytime, anywhere: in-practice, theatre, patient side or equine yard
- Wi-Fi and Bluetooth enabled easy transfer of data



RAPIDPoint® 500e Blood Gas System

Siemens Healthineers RAPIDPoint® 500e Blood Gas System is designed to meet the challenges of Point of Care testing by providing fast, accurate and comprehensive test results.

Features & Benefits

- Accurate, lab-quality test results in 60 seconds
- Easy to use, reliable, maintenance-free operation with cartridge based system, automatic calibration and QC
- Robust accuracy for sample to sample by leveraging the Siemens Healthineers designed Integri-sense™ technology
- Onboard software defences provide reassuring level of protection against cyber attacks and confidential data theft
- Simple touch screen simply press start after inserting a sample
- Preset custom panels
- Hands-free automatic sampling with universal sample port
- Integrated barcode scanner

RAPIDPoint[®] 500e Tests

Blood Gas: pH, pCO,, pO, Electrolytes: Na+, K+, Ca++, Cl-

Metabolites: Glu, Lac

Co-oximetry: nBili, tHb, SO,, O, Hb, HHb, COHb, MetHb





Ideal for sites with a larger workload

Clinispin Vet Combi Multi Purpose Veterinary Centrifuge



Versatile with Combination Rotor

The Clinispin Vet Combi is a high speed centrifuge equipped with a maintenance free drive, digital display and simple interface for silent and efficient operation.

The combination Rotor provides efficient separation of blood and urine samples and accurate Haematocrit determination.

- 3 pre-defined presets for easy and rapid operation at one press of a button:
 - Serum/Plasma (12,000 RPM for 3 minutes)
 - HCT/PCV (12,000 RPM for 5 minutes)
 - Urine (1,500 RPM for 2 minutes)
 Presets can be varied for each centrifugation if required
- Rotor allows for quick and easy sample loading
- Maintenance free drive
- Digital display
- Easy to use
- Quiet operation
- Displays RPM or RCF (g-force)

Clinispin horizon 6 & 24 Centrifuges

Clinispin horizon 6 & 24 are versatile centrifuges with 3 default preset cycles for chemistry, coagulation and urine.

- Easy To Use simple two-button operation to select a preset cycle or lock to a single cycle
 - 3 preset cycles
 - 6 and 24 tube capacity
 - Erasable label for custom cycles
- Lid Lights flash at the end of each cycle and continue flashing until the lid is opened
- Automatic Lid Lock will not run with lid open and lid cannot be opened while samples are spinning
- Cool Operation unique design of the convection based cooling system protects the samples even under the heaviest workloads
- Clear Shatterproof Lid for easy monitoring
- Brushless DC Motor no routine maintenance
- Reinforced Guard Bowl maximum safety and durability
- 2 Stage Horizontal Rotor patented swing-out design produces ideal horizontally separated samples

Minimise error, reduce training time and deliver the best patient outcome with Clinispin horizon Routine Centrifuges



Clinispin horizon 6
- ideal where space
is at a premium and
consistency is key



Clinispin horizon 24
- ideal for sites with
a larger workload

In-Practice Urinalysis in 24 hours simplified with


InSight Urinalysis

A 2-stage process available 24 hours a day

1. RapidBac™ Vet TEST FOR BACTERIURIA

- * Negative no bacteriuria indicated
- * Positive proceed to Flexicult Vet if antibiotic sensitivity is required





Results in 20 minutes

RapidBac[™] Vet is a lateral flow immunoassay for the detection of Gram-positive and Gram-negative bacteria in Dog and Cat urine.

- Test immediately in-practice
- Results in less than 20 minutes
- Easy to use
- Room temperature storage
- 98.5% accuracy
- No special apparatus required

Proceed with Flexicult® Vet if RapidBac™ is positive and antibiotic sensitivity is required

2. Flexicult® Vet

Chromogenic In-Practice
Urine Culture & Sensitivity Test



Flexicult® Vet is a 24 hour chromogenic culture test for diagnosing UTI's in Dogs and Cats. It aids the identification and quantification of bacteria and helps select the most appropriate antibiotic for effective patient management.

Results within 24 hours

Flexicult Vet contains the following antibiotics:

Ampicillin

Perform in-house

without the need

for specialist knowledge

- Amoxicillin/clavulanate
- Enrofloxacin
- Trimethoprim/ sulfamethoxazole

Agar plate contains a chromogenic substance
Causes bacterial colonies on

the agar to appear with different colours depending on the bacterium species





Examples of Flexicult Vet inoculated with E. coli (left) and E. faecalis (right)

Incubator required - available from Woodley





Urinalysis is an important tool in disease detection, as well as monitoring and screening animal health. Abnormalities can be indicative of diseases of the urinary system as well as other organ systems, including liver function, acid base status and carbohydrate metabolism.

InSight® Urine Test Strips

MS-11 Urine Strips Test for:

- Urobilinogen
- Blood
- Glucose
- pH
- Bilirubin
- Protein
- Ketone
- Creatinine
- Specific Gravity Nitrate
- Leucocytes with urine protein/creatinine ratio

MS-2 Urine Strips Test for:

• Microalbumin • Creatinine with urine microalbumin/creatinine ratio



Multistix® 10 SG Reagent Strips

The Multistix® 10 SG reagent strip is the reliable frontline test for detection of a broad range of conditions, from detecting urinary tract infections (UTI1) to diabetes and kidney disorders. Accurately gain a broader clinical insight into a patient's health.

- Reliable testing for UTI, diabetes and kidney disorders
- 10 tests for important markers on a single strip
- Combined leukocyte/nitrile results for better prediction of UTI with fewer false negatives





InSight® Urinalysis Reader

- Semi-automatic urine test strip reader with touch screen & printer to be used with InSight **Urine Test Strips**
- Simply dip the strip into the urine sample and place it on the strip loading plate
- Results displayed in 60 seconds



Urinalysis Atlas

- Sample collections
- Physical characteristics
- Examination of sediment
- Photographic identification of casts
- Contaminants in urine sediment

Refractometer

- Handheld veterinary refractometer with adjustable calibration
- Provides quick and accurate indication of vital fluid levels.
- The triple scale reads: Urine Specific Gravity, Total Serum Protein and Refractive Index
- Refractometer equipped with ATC (Automatic Temperature Compensation)



EZ-SmartTemp™ Min Max Refrigeration Thermometer

- Easy to use simple and small design with easy to read digital display
- This calibrated thermometer attaches to the exterior of the refrigeration unit and measures the internal temperature via a probe which is sealed in a non toxic glycol solution
- Bottle insulates sensor from rapid temperature changes when the door is opened
- Triple display simultaneously shows the miniumum, maximum and current temperature
- Alarm feature alerts if temperature rises above or falls below a set point programmable in 1°C increments
- Supplied with velcro and magnetic strips for mounting
- Supplied with calibration certificate

EZ-SmartLog[™] One USB Temperature Datalogger

- Measures and stores up to 16,382 temperature readings over a -35°C to +80°C range
- Direct USB connection for easy set-up, functionality and data download
- User-programmable alarm thresholds
- Status indication via red, green and orange LEDs
- Complete with Windows control software which is easy to download and install
- Supplied with replaceable lithium battery typically allows logging for 1 year





Hematek® 3000 System Slide Stainer

The Hematek 3000 System provides semi-automated slide staining for laboratories that need a robust system capable of staining up to 60 slides per hour. The system helps maximise productivity and deliver quality results through user-defined features and precision controls.

- Continuous-feed loading platform eliminates wasted time associated with batch staining methods
- Individual, consistent slide staining optimises staining quality and enables user-adaptable intensity for standardised slide staining
- Minimal system maintenance helps maximise uptime in your lab
- Compact, RoHS-compliant system guarantees incorporation of environmentally friendly, sustainable substances and supports 'green' initiatives
- All-in-one stain pack includes enough stain for up to 900 slides and offers easy load-and-go operation for optimised daily workflow with minimal hands-on time
- Pumps deliver a fresh, precise volume of reagent for every slide, while the platen provides uniform stain from field to field and slide to slide with accurate mixing and uniform stain and buffer delivery
- Simple pump release mechanism streamlines maintenance
- Slides are immediately available for viewing upon completion of staining and the rinse reagent contains the optimal alcohol concentration for proper slide drying
- System warnings alert users of low stain-level and full waste bin to help ensure uninterrupted, quality results

Approved Partner Calmentaly Completed as Healthineers



Practical Automation in Slide Staining

Approved Partner SIEMENS ... Labeling Degreeous Healthingers ...

RAL Diff-Quik™

RAL Diff-Quik is a simplified and cost-effective alternative to Wright-Giemsa staining and is faster than typical staining procedures. The speed of execution makes it an ideal stain for emergency diagnostics of medullary and bloods smears with a **process time of just 15 seconds**.

Easy

- 3-step method is very easy
- Consists of 3 brief immersions in each of the RAL Diff-Quik solutions
- This manual protocol allows the user to customise and vary the staining intensity
- The different degrees of shading allow preferential staining for every lab

Applications:

- Haematology: Blood smears
- Veterinary Parasitology: Piroplasma
- Haematology: Bone marrow
- Cytology: Spermocytograms
- Cyto-bacteriology: Cytology of biological fluids, urines, cytopunctures and CSF

Ready to Use

 A new packaging has been created for direct staining without the need to decant stains



Convenient pack sizes:

- 3 x 100ml
- 3 x 500ml

Process Time of just 15 Seconds

SIEMENS

Healthineers

ADVIA® 2120i Haematology System

Trusted haematology solutions for advanced multispecies testing

The ADVIA 2120i Haematology System is the only haematology system for human and animal samples offering complete multispecies haematology testing capability with practical automation.

- High throughput capacity of 120 samples per hour with sophisticated technologies for analysing White Blood Cell (WBC) Differentials, Red Blood Cells (RBC), Reticulocytes and Platelets
- Accepts up to 150 samples at a time
- 21 pre-programmed species provide ready-to-use cell gating across a wide spectrum of species and an unlimited number of user-definable archetypes allow easy expansion to new species
- Unique playback function allows replay of sample results in alternate species without redrawing or retesting
- Data manager provides full traceability at no additional cost, so you can run customised reports at any time



High Capacity Multispecies Tests

- Autoslide connectivity
- NRBCs on every sample
- Body Fluid Capability
- Unique Unifluidics™ Technology
- Efficient and cost effective workflow
- Fully automated daily maintenance
- "Right first time" results

IMMULITE® 2000 XPi Immunoassay System

The highly reliable IMMULITE® 2000 XPi System is an easy to use, continuous random access immunoassay analyser with one of the largest automated immunoassay menus available and a maximum throughput of 200 tests per hour. Medium to high-volume laboratories can incorporate specialty and allergy testing into routine workloads to improve productivity and efficiency.

Up To 200 Tests/Hour Easy To Use

- Customised automated startup for daily routine maintenance and quality control scheduling
- Reagent onboard stability of 90 days, minimises waste
- Maximum throughput of 200 tests per hour
- No pause sample loading to enhance turnaround time
- Tube-top sample cups for micro-sample handling
- Home screens with at-a-glance status on a 19" touch screen monitor
- Options available for direct water feed and bulk solid waste chute
- Siemens Automation solutions available through VersaCell® Systems



IMMULITE Veterinary Test Menu:

- Total T4
- TLI
- TSH
- Digoxin
- FT4 • ACTH
- Insulin
- Cortisol
- Troponin I

Phenobarbitol

- Progesterone Testosterone
- Vitamin B12
- Folic Acid

This list is not exhaustive - please contact us with your requirements





SRUC



External Quality Assessment (EQA)

External Quality Assessment (EQA) of laboratory equipment is required by the RCVS Practice Standards Scheme and is recommended for any practice with in-practice analysers. In conjunction with the SRUC (Scotland's Rural College), Woodley Equipment are proud to introduce a 'true' External Assessment with...



A certified independent External Quality Assessment scheme enabling you to provide the best possible care

For Biochemistry, Endocrinology, Blood Gases, Haematology & Coagulation



InSight EQA is independent of all analyser manufacturers

Why use EQA?

- · Maintaining and improving analytical quality
- Improving inter-laboratory agreement thus raising standards
- Detecting equipment failures and identifying reagent problems
- Reviewing staff training
- Initiating and evaluating corrective actions
- Comparing different analytical methods
- Compliant with the RCVS Practice Standards Scheme

Why choose EQA?

- InSight EQA is an internet based EQA scheme for biochemistry, endocrinology, blood gases, haematology and coagulation
- InSight EQA provides a 'true' external assessment by comparing results from multiple analysers of the same make and NOT just with a reference laboratory analyser
- InSight EQA is independent of all analyser manufacturers
- Flexible scheme so you choose which analytes to assess

3 Simple Steps

- 1. We send you a sample to analyse each month
- 2. You email back the results
- 3. We send you a report within 7 days comparing your results with all subscribers for both identical analysers and all registered with the scheme



InSight EQA is a 'NOT FOR PROFIT' scheme in conjunction with Scotland's Rural College (SRUC)

For further information on SRUC's activities visit www.sruc.ac.uk





Internal Quality Control Material

Here at Woodley Equipment, we think that accurately verified analyser test results are essential. That is why we are providing Woodley InSight Analyser customers with free monthly Quality Control material.

The veterinary laboratory equipment market has grown considerably over the last 20 years with the majority of UK practices now performing in-house testing. With more surgeries than ever before managing their own laboratories, there has never been a greater need for quality assurance (QA) programs.

Quality controls (QC) form part of a quality assurance (QA) program to guide the work of the practice laboratory and assess its overall quality and performance. This assures the veterinarian that test results generated are both accurate and precise. The Royal College of Veterinary Surgeons Practice Standards states that internal laboratory QC procedures used in conjunction with an external quality assurance scheme must be routinely undertaken to create the uniform high quality that makes laboratory practice results reliable.

In order to assist practice laboratories assess the quality and performance of their analysers, Woodley Equipment will be providing free haematology in-house QC material and training to all of our Woodley InSight analyser customers.









If you would like more information please contact Woodley Equipment Company

Insight Analysis

"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"

































01204 669033 (Option 1)



sales@woodleyequipment.com



www.woodleyequipment.com

Old Station Park Buildings St. John Street Horwich Bolton Lancashire BL6 7NY UK

"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"

