# **WOODLEY EQUIPMENT COMPANY** nelivering etcel **PRODUCT GUIDE**

Issue 2.3









"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"





# WELCOME TO OUR PRODUCT GUIDE

Woodley Equipment Company was **established in 1989** and prides itself on being a market leading supplier of specialist **laboratory equipment**, **diagnostic and critical care** products to the **veterinary industry worldwide**.

A dedicated on-the road Account Manager team provides **in-practice support**, training and on-going advice. This team is supported at our headquarters with a technical support and service department staffed with **Biomedical Scientists** and **Service Engineers** who deliver proven high quality after-sales care.

As well as a **UK nationwide sales force**, Woodley has an international distribution network offering a rapidly growing portfolio of **high quality products** with worldwide distribution rights. Woodley also has a proven track record and outstanding **customer loyalty** in the form of repeat business and client referrals.

Woodley seeks to supply veterinarians and distributors with the equipment and services to best support Practices, thus enabling them to respond to the needs of their clients, while operating a more **efficient and profitable business**.

**Quantum Vet Diagnostics**, a division of Woodley Equipment, is the **UK's number 1** provider of veterinary wet chemistry solutions.



## "PASSIONATE ABOUT LABORATORY DIAGNOSTICS"







01204 669033 (Option 1)

sales@woodleyequipment.com

HCT results reported as %

www.woodleyequipment.com

## **i-SMART 30 VET ELECTROLYTE ANALYSER**

#### Fast & Accurate

- Fast, accurate analysis
- Tests electrolytes Na<sup>+</sup>, K<sup>+</sup> & Cl<sup>-</sup>
- 60µl sample volume
- Results in 35 seconds
- Species specific reference ranges

#### **All-in-One Cartridge**

- No user maintenance
- Autocalibration
- Cartridge contains electrodes, all consumables & waste
- Veterinary specific software

Fast, accurate & maintenance free In-Practice electrolyte analysis

#### PocketChem<sup>™</sup> BA Blood Ammonia Measurement

PocketChem BA measures blood ammonia levels. The device enables immediate testing and delivers fast, accurate results.

- Whole blood measurement
- Ready for use 24 hours a day
- Small sample volume
- Battery powered
- Easy to maintain
- Temperature compensation
- Lightweight and compact design
- Detachable printer available
- Supports the diagnosis of hepatic cephalopathy
- Helps detect and evaluate urea cycle defects
- Monitor patients with liver disease

## InSight<sup>™</sup> Veterinary Blood Urea **Test Strips**

A quick and easy, low cost test for in-practice quantitative determination of urea in whole blood.

- One drop of venous or capillary whole blood
- Test time only 90 seconds
- Colour chart for reading test results is printed on the bottle label



## StatStrip<sup>®</sup> Lactate Xpress Meter

- Results in 13 seconds
- Easy to use and portable
- Small 0.6µl sample size
- Venous or arterial whole blood
- Corrects results for HCT
- No calibration required
- Eliminates interfering substances
- Stores 400 patient results

#### StatSensor<sup>®</sup> Creatinine Xpress Meter

- Rapid assessment of kidney function
- Results in 30 seconds
- Easy to use and portable
- Small 1.2µl sample size
- Capillary, venous or arterial whole blood
- Correct results for HCT
- No calibration required
- Eliminates interfering substances
- Stores 400 patient results

#### **NovaVet Blood Ketone Meter**


- Measures blood beta-hydroxybutrate (BHBA), the gold standard for detecting subclinical ketosis
- Easy to use and portable
- Predict and monitor diabetic ketoacidosis
- No calibration required
- Simple test procedure in the field
- Results in 10 seconds
- Small 0.8µl sample size
- Stores 400 results

## **Poison Test - Anti-Freeze Ethylene Glycol Testing**

Tests for the presence of Ethylene Glycol commonly found in 'anti-freeze', lacquer, detergents & de-icers. These test strips report quantitative results of Ethylene Glycol in plasma.

- Feline & Canine sensitivity at 20 mg/dL
- Only 20µl of Heparin plasma required
- Detects Ethylene Glycol approximately 30 minutes after ingestion
- Results in 10 minutes
- Sample window 12-14 hours post ingestion (varies by size, age and condition of animal)
- Room temperature storage







1021 "BIG"

007



01204 669033 (Option 1)

www.woodleyequipment.com



**Measured electrolytes** 

20 - 250 mmol/L

0.5 - 20.0 mmol/L

20 - 250 mmol/L

Na<sup>+</sup>

K<sup>+</sup>

Cl-

Point of Care Blood Ammonia



## Veterinary Glucose Meter

g-Pet PLUS Veterinary Glucose Meter is calibrated specifically for use in Dogs, Cats and Horses, ensuring an accurate result is reported for each species.

- Biosensor auto code technology
- Accurate results
- Easy to use
- $\bullet$  Small sample of whole blood (7µl)
- Accurate results displayed within 5 seconds
- Readable range 1.1-33.3 mmol/L
- 500 result automatic memory
- HCT range 20-60%
- Lifetime warranty

# **InSight**<sup>™</sup> Vet Vital Signs

InSight Vet Vital Signs Handheld Monitor is ideal for use in Veterinary Critical Care, General Health Screening and Anaesthesia Monitoring.

- Small handheld monitor
- Measures NIBP + SpO<sub>2</sub> +TEMP + Heart Rate
- Includes cage clip & charging dock
- Easy to use
- Accurate measurement in 0.05% Low Perfusion Index
- Audio and visual alarm
- Records up to 48 hours of measurement data



07-15

11: 13\*\*

Calibrated for

Canine, Feline &

Equine use



## Veterinary Lactate Meter

L-Pet only requires a small drop of blood (<3µl). The test strips draw the sample in and then measure using an advanced biosensor specific for lactate concentration in blood.

- Advanced Biosensor Test Strips
- Small sample size (<3µl)
- Results in less than 45 seconds
- Large LCD display
- Stores up to 300 results
- Easy operation
- Portable
- Specific Lactate Analyser QC Material available
- 2 year warranty

# **InSight**<sup>™</sup>Vet Capnograph

InSight Vet Capnograph Monitor is ideal for use in Veterinary Critical Care and Anaesthesia Monitoring.

- Small handheld monitor
- Use to measure SpO<sub>2</sub> + ETCO<sub>2</sub> + Heart Rate
- Includes cage clip & charging dock
- Easy to use
- Accurate measurement in 0.05% Low Perfusion Index
- Audio and visual alarm
- Records up to 48 hours of measurement data



#### AliveCor<sup>®</sup> iPhone<sup>®</sup> Veterinary Heart Monitor An ECG device as mobile as you are

Veterinary iPhone ECG and Heart Monitor is designed to snap onto your iPhone like a case and wirelessly communicates with the AliveECG Vet app. Once partnered with the ECG monitor, it transforms your iPhone into a clinical-quality, single-lead ECG recorder.

• Easy to use - No stress or shaving required in most animals

Measure Blood Pressure,

SpO2, Heart Rate & Temperature

- Powerful Leverages power, display and communication capabilities of your iPhone
- Accurate Veterinary testing demonstrates accuracy in rate and rhythm compared to traditional devices
- Accessible Securely store ECG notes and patient information on the device and in the cloud



#### **Compatible With:**

• iPhone 4/4S/5/5S/6/6S/SE/7



An easy to use handheld meter providing blood lactate results in just 45 seconds

# Vcheck V200 point of care analyser

## Simplify In-Practice Point of Care Testing with Reference Lab Quality - avoiding test referrals

## Canine: cPL, T4, TSH, Cortisol, CRP, Progesterone Feline: T4, fPL, SAA,

## Achieve true analysis with fluorescence technology

- Fully quantitative automated analyser
- Accurate diagnosis for Dogs and Cats in minutes
- Analyser provides print out and result storage



#### **Canine Tests Available:**

BIONOTE

- cPL
- T4
- TSH
- Cortisol
- CRP
- Progesterone

#### Feline Tests Available:

- T4
- fPL
- SAA

Low Cost, In-Practice, Reference Lab Quality Immunoassay Testing

## **Small analyser - low cost per test**



Results available on your smartphone via:

A stand alone Lab Management System to integrate with your Practice Management System





# Rapid Test Kits Available exclusively from Woodley

## **Canine Tests Include:**

#### Anigen Rapid CPV Ag

For the detection of Canine Parvo virus antigen in Canine faeces.

Anigen Rapid CDV Ag

For the detection of Canine Distemper virus antigen in conjunctiva (recommended), urine, serum or plasma.

**Anigen Rapid CCV Ag** For the detection of Canine Corona virus antigen in Canine faeces.

Anigen Rapid CPV Ag/CCV Ag For the detection of Canine Parvo virus and Corona virus antigen in Canine faeces.

Anigen Rapid CDV Ag/CAV Ag For the detection of Canine Distemper virus and Adeno virus antigen in conjunctiva (CDV Ag) and nasal fluid (CAV Ag).

Anigen Rapid CIRD-3 Ag For the detection of Canine Infectious Respiratory Disease (CDV, CAV and CIV) in conjunctiva (CDV Ag) and nasal fluid (CAV and CIV Ag).

Anigen Rapid CPV/CCV/Giardia Ag For the detection of Canine Parvo, Corona and Giardia antigen in Canine faeces.

## Feline Tests Include:

#### Anigen Rapid FHW Ag/Ab

For the detection of Feline heartworm antigen and antibody in whole blood, plasma or serum.

Anigen Rapid FeliCheck-3 For the detection of Feline Immunodeficiency virus antibody, Feline Leukemia Virus antigen and Feline Heartworm antibody in whole blood, plasma or serum.

**Anigen Rapid Toxoplasma Ab** For the detection of Toxoplasma gondii antibody in Feline whole blood, plasma or serum. Anigen Rapid Rabies Ag

For the detection of Rabies virus antigen. Anigen Rapid Giardia Ag For the detection of Giardia antigen

in Canine faeces. Anigen Rapid CPV Ab

For the titration of Canine Parvo virus in whole blood, plasma or serum. Anigen Rapid CDV Ab

For the titration of Canine Distemper virus antibody in whole blood, plasma or serum.

**Anigen Rapid C.Brucella Ab** For the detection of Canine Brucella antibody in whole blood, plasma or serum.

Anigen Rapid Leishmania Ab For the detection of Leishmania imfantum antibody in whole blood, plasma or serum.

Anigen Rapid E.Canis Ab For the detection of Ehrlicia Canis antibody in whole blood, plasma or serum.



#### Anigen Rapid FeLV Ag

For the detection of Feline Leukemia virus antigen in Feline whole blood, plasma or serum.

Anigen Rapid FIV Ab

For the detection of Feline Immunodeficiency virus antibody in Feline whole blood, plasma or serum. Anigen Rapid FCoV Ag/Ab

For the detection of Feline Corona virus antigen/antibody in Feline faeces (Ag) and whole blood, plasma or serum (Ab). **Anigen Rapid FPV Ag** For the detection of Feline Parvo

virus antigen in Feline faeces.



#### Anigen Rapid Anaplasma Ab

For the detection of Anaplasma phagocytophilum/Anaplasma platys antibody in whole blood, plasma or serum.

#### Anigen Rapid CaniV-4

For the detection of Dirofilaria immitis antigen, Ehrlichia canis antibody, Leishmania imfantum antibody, Anaplasma, phagocytophilum/ Anaplasma

platys antibody in whole blood, plasma or serum.

#### Anigen Rapid E.canis/ Anaplasma Ab

For the detection of Ehrlichia canis and Anaplasma phagocytophilum/ Anaplasma platys in whole blood, plasma or serum.

#### Anigen Rapid Total IgE

For the detection of Canine Total IgE in serum.



#### Anigen Rapid FIV Ab/FeLV Ag

For the simultaneous detection of Feline Immunodeficiency virus antibody and Feline Leukemia virus antigen in Feline whole blood, plasma or serum.

Anigen Rapid FCoV Ag

For the detection of Feline Corona virus antigen in Feline faeces.

#### Anigen Rapid FCoV Ab

For the detection of Feline Corona virus antibody in whole blood, plasma or serum. **Anigen Rapid FeliD-3 Ag** For the detection of Feline Parvo, Corona and Giardia antigen in Feline faeces.

## This list is not exhaustive - please contact us with your requirements

01204 669033 (Option 1)

# **NOW YOU HAVE A CHOICE** The 'Woodley Lab Systems'

Haematology, Clinical Chemistry, Immunofluorescence, Coagulation & Urinalysis Analysers



for one low cost payment each month

- Includes Technical Support for 5 years
- Includes all repairs and servicing for 5 years
- Includes scheduled internal QC
- Includes 'Track your delivery'
- Includes InSight Lab Manager

- Includes software updates to improve analyser functionality and test range when appropriate
- Wet Chemistry, Endocrinology and Blood Gas options available

#### **Contact Woodley or your Local Area Account Manager to discuss**

01204 669033 (Option 1) sales@woodleyequipment.com www.woodleyequipment.com

## A stand alone Lab Management System to integrate with your Practice Management System

**InSight**<sup>™</sup>

Lab Manager





- InSight Lab Manager integrates your Lab Equipment with the host Practice Management System
- Alternatively, InSight Lab Manager can be used as a stand alone Lab Management System
- Results are displayed alongside a Reference Range Bar Chart for easy identification of abnormal results

### **Contact Woodley or your Local Area Account Manager to discuss**

Wet Chemistry, Endocrinology & Blood Gas Options





01204 669033 (Option 1)

sales@woodleyequipment.com

www.woodleyequipment.com

# In-Practice Urinalysis in 24 hours simplified with

**Insight**<sup>™</sup>Urinalysis

## A 2-stage process available 24 hours a day

### **1. RapidBac<sup>™</sup> Vet** TEST FOR BACTERIURIA

10

\* Negative - no bacteriuria indicated\* Positive - proceed to Flexicult Vet





#### Results in 20 minutes

RapidBac<sup>™</sup> Vet is a lateral flow immunoassay for the detection of Gram-positive and Gram-negative bacteria in Dog and Cat urine.

- Test immediately in-practice
- Results in less than 20 minutes
- Easy to use
- Room temperature storage
- 98.5% accuracy
- No special apparatus required

#### Proceed with Flexicult<sup>®</sup> Vet if RapidBac<sup>™</sup> is positive or if antibiotic confirmation is needed

## 2. Flexicult<sup>®</sup> Vet

Chromogenic In-Practice Urine Culture & Sensitivity Test



Flexicult<sup>®</sup> Vet is a 24 hour chromogenic culture test for diagnosing UTI's in Dogs and Cats. It aids the identification and quantification of bacteria and helps select the most appropriate antibiotic for effective patient management.

#### Results within 24 hours

Flexicult Vet contains the following antibiotics:

- Ampicillin
- Amoxicillin/clavulanate
- Cephalotin

Perform in-house

without the need

for specialist knowledge

- Enrofloxacin
- Trimethoprim/ sulfamethoxazole

## Agar plate contains a chromogenic substance

Causes bacterial colonies on the agar to appear with different colours depending on the bacterium species



Examples of Flexicult Vet inoculated with E. coli (left) and E. faecalis (right)



Incubator required available from Woodley

01204 669033 (Option 1)



Urinalysis is an important tool in disease detection, as well as monitoring and screening animal health. Abnormalities can be indicative of diseases of the urinary system as well as other organ systems, including liver function, acid base status and carbohydrate metabolism.

#### **InSight™ Urine Test Strips**

Blood

Protein

• pH

#### MS-11 Urine Strips Test for:

- Urobilinogen
- Glucose
- Bilirubin
- Ketone

ne • Creatinine ific Gravity • Nitrate

Specific Gravity
N
Leucocytes

with urine protein/creatinine ratio

#### MS-2 Urine Strips Test for:

• Microalbumin • Creatinine with urine microalbumin/creatinine ratio



WOODLEY

### **InSight™ Urinalysis Reader**

- Semi-automatic urine test strip reader with touch screen & printer to be used with InSight Urine Test Strips
- Simply dip the strip into the urine sample and place it on the strip loading plate
- Results displayed in 60 seconds

#### **Urinalysis Atlas**

- Sample collections
- Physical characteristics
- Examination of sediment
- Photographic identification of casts
- Contaminants in urine sediment

Use in conjunction with RapidBac Vet & Flexicult Vet for complete In-Practice Urinalysis

### InSight<sup>™</sup> Microscope Eyepiece Adapter for Smartphones



- Mounts any smartphone to any microscope in minutes!
- Smartphone camera can then capture photos & videos for real-time viewing on smartphone screen, via video conferencing, emailing or sharing online

Capture & share microscopic images

#### Vision V5000 Microscope

- LED illuminator
- Excellent optical resolution
- Wide-field eyepieces
- Dioptric compensation
- DIN achomatic objectives 4x, 10x, 40x, 100x
- Coaxial coarse and fine focusing

#### InSight<sup>™</sup> 50x Microscope Lens

- Mounts onto any Smartphone
- Enhances the Smartphone camera
- Captures high quality microscopic images & video
- No batteries required





# ImmunoComb<sup>®</sup> VacciCheck<sup>®</sup> In-Practice Antibody Test Kits

ImmunoComb VacciCheck Antibody Test Kit monitors serological status, reduces vaccination failure and prevents over vaccination consequences for Dogs and Cats.

VacciCheck is recommended by the WSAVA to comply with regulations to reduce unnecessary vaccinations

# ACCICheck Canine

The Canine VacciCheck Antibody Test Kit is designed to determine Dog serum antibody titre to:

- Infectious Canine Hepatitis (ICH)
- Canine Parvovirus (CPV)
- Canine Distemper Virus (CDV)

The main purpose of this kit is to provide a useful tool to assess the immune status of Dogs to the three core Canine infections.

# ACCICheck Feline

The Feline VacciCheck Antibody Test Kit is designed to determine Cat serum IgG antibody titre to:


- Feline Panleukopenia (FPLV)
- Feline Herpes Virus (FHV)
- Feline Calici Virus (FCV)

The main purpose of this kit is to provide a useful tool to assess the immune status of Cats to the three core Feline infections.

WSAVA Guidelines state that titre testing for the 3 core Canine and Feline pathogens should be performed instead of automatic annual revaccination







# **Clinispin MPC Combi Veterinary Centrifuge**

Quick & easy to use



Four Tube Holders



The Clinispin MPC Combi is a high speed centrifuge equipped with a maintenance free drive, digital display and simple interface for silent & efficient operation. The combination rotor provides efficient separation of blood and urine samples and accurate Haematocrit determination.

- 3 pre-defined presets for easy and rapid operation at the press of a button:
  - Serum/Plasma (12,000 RPM for 2 minutes)
  - HCT/QBC<sup>®</sup> (12,000 RPM for 5 minutes)
  - Urine (1,000 RPM for 1 minute)
- Presets can be varied for each centrifugation if required
- Rotor allows for quick and easy sample loading
- Maintenance free drive
- Digital display
- Easy to use
- Quiet operation
- Displays RPM or RCF (g-force)

## RapidVet<sup>®</sup>- H Blood Typing Kits - Canine & Feline

In-practice tests for rapid determination of Canine DEA 1.1 and Feline A, B or AB blood types.

These fast in-house blood typing card kits for Dogs and Cats are essential for transfusion and breeding decisions. They are fast, easy to use, economical and can be stored at room temperature.


#### **RapidVet-H Canine**

- For Canine DEA 1.1
- Fast, easy, economical and VITAL
- Results in 2 minutes
- Auto-agglutination screen
- Requires 3 drops of whole blood
- Available in kits of 5



#### RapidVet-H Feline

- For Feline A,B & AB
- Fast, easy, economical and VITAL
- Results in 2 minutes
- Auto-agglutination screen
- Requires 3 drops of whole blood
- Available in kits of 5



"Every Dog & Cat in need of a transfusion should be blood typed" Dr. Urs Giger, Dr.med.vet., Dip ACVIM, ECVIM

## RapidVet<sup>®</sup>- H Crossmatch Kits

In-practice crossmatch test kits for use in tranfusion medicine on either Canine or Feline patients.

- Will detect all major Canine and Feline blood group incompatibilities
- Performed in-practice in 10 minutes
- Easy to use and results are clear
- Available in kits of 3 tests
- Minor Crossmatch Kits also available

Essential for good transfusion



#### Use alongside Canine & Feline Blood Typing Kits

### Veterinary Parasep<sup>®</sup> Faecal Filters

- Determine worm egg counts in Cattle, Sheep and Horses
- Cleaner, safer method
- Separates eggs from faecal filter in safe, enclosed environment
- No possibility of cross contamination
- More accurate and consistent results
- No laborious, time consuming 'clean up' of sinks, shaker bottles and sieves

## **Haematology Atlas**

Haematology Atlas written exclusively for Woodley Equipment by Kostas Papasouliotis DVM PhD DipRCPath MRCVS (Bristol Vet School).



Titles Available - Canine, Feline and Equine

These atlases are an excellent microscope-side reference guide offering normal cell morphology.

• Erythrocytes

• Lymphocytes

• Eosinophils

- Neutrophils
- Basophils
- Platelets
- Methodology: Microhaematocrit
  - Methodology: Blood Smears

# **InSight**<sup>™</sup>DB



A fully automated clinical chemistry analyser which uses liquid reagents and is approximately 70% cheaper to run than dry chemistry

## **Chemistry Analyser**

- User friendly intuitive veterinary software
- Onboard refrigerated reagent/sample compartment
- Minimal user maintenance
- Fast throughput example batch times: 9 mins 1x GHP, 13mins 2x GHP, 16mins 3x GHP
- Continuous loading, no need to wait for batch completion
- Stat function for urgent tests
- Fully automated with auto-reflex dilution of over range results
- Flexible testing, single parameter or up to 20 user defined profiles
- Small footprint W 414mm x D 565mm
- Online remote assistance directly to your analyser
- Long life ceramic syringes to ensure accuracy and precision
- Environmentally friendly semi-disposable reaction cuvettes with high efficiency wash system

## Customer defined profiles and/or tests can be programmed from the following:

Albumin, Alkaline Phosphatase, ALT, AST, Amylase, Bilirubin (Total), Calcium, Cholesterol, CK, Creatinine, Glucose, GGT, Phosphate, Protein (Total) & Urea.

BHB, Bile Acids, Chloride, Copper, Fibrinogen, Fructosamine, GLDH, GSH-Px (as Selenium), Haemoglobin, IgG, LDH, Lipase, Magnesium, NEFA, SAA, T4 (Total Thyroxin), Triglycerides, Uric Acid, Urine Creatinine, Zinc, Phenobarbitone, Iron, Ferritin and TIBC

#### Fully supported Nationwide by our Specialist Team of Field Service Technicians





## A stand alone Lab Management System which can integrate to your Practice Management System

- InSight Lab Manager integrates your Lab Equipment with the host Practice Management System
- Alternatively, InSight Lab Manager can be used as a stand alone Lab Management System
- Results are displayed alongside a Reference Range Bar Chart for easy identification of abnormal results
- Smartphone and tablet friendly access Lab Results anytime, anywhere



Take a 360° look at your endocrinology testing with the TOSOH AIA-360

## AIA-360 Analyser

Offers the automation powers and performance of a large analyser in a benchtop format. Random access allows any combination of tests, in any order at any time and a throughput of 36 results per hour with the first test result in 20 minutes ensures a fast reliable service.

- Simple to use
- Fast, reliable results
- Simple maintenance
- Touchscreen with intuitive menus

• ACTH

• Insulin

#### Key Assays Available:

- Free T4
- Total T4
- Cortisol
   CTnl
- Progesterone



01204 669033 (Option 1)

www.woodleyequipment.com

TOSOH



## **Haematology Analysers**

A fully automated multi-species Haematology Analyser

5 Part Differential with Laser Technology



A true 5 part differential laser Haematology Analyser including Reticulocytes

Uses dual laser and impedance technology

- Quick and easy to use
- Large user-friendly colour touch screen
- Low maintenance
- Fully automated
- 23 parameters
- 60 second analysis time
- Small sample volume 15µl whole blood
- Stores up to 50,000 results
- Dog, Cat, Rabbit, Horse & Cow Settings

#### Both analysers Integrate with Practice Management Systems

- Compact, benchtop size
- Small sample size 20µl whole blood
- 27 parameters
- 2 minute cycle time
- Quick and easy operation
- Low maintenance
- Touch screen
- Over 50 species on board.



(Tesigi



## Skyla VB1 Veterinary Chemistry Analyser

For both Rotor Profiles & Single/Dual Tests

Highly accurate test results in under 15 minutes

Integrates with InSight Lab Manager - a stand alone Lab Management System which can integrate to your Practice Management System





- Fast and reliable
- Compact, benchtop size
- Performs up to 15 biochemistry tests simultaneously
- Performs up to 3 single/dual tests simultaneously
- 200µl whole blood, serum or plasma
- LCD colour touch screen & built-in printer
- Fully automatic system with walk away operation
- Built-in QC and calibration for each test
- Species specific reference ranges
- Stores up to 50,000 results

#### **Parameters Available**

ALB, ALP, ALT, AMY, AST, BUN, Ca, Cl<sup>-</sup>, CREA, CPK, GGT, GLU, K<sup>+</sup>, Na<sup>+</sup>, PHOS, TBIL, CHOL, TP, B Acids, T4, Mg, UA,  $TCO_2$ , TRIG, LAC, FRU, LIPA, UCRE, UPRO, GLOB<sup>\*</sup>, A/G<sup>\*</sup>, Na<sup>+</sup>/K<sup>+</sup>, BUN/CREA<sup>\*</sup>, UPC<sup>\*</sup>, UREA<sup>\*</sup> \* Calculated



Singl	e Test	Dual	Test

- Na⁺ & K⁺
- CREA • Cl<sup>-</sup>
- AMY
- BUN
- TBIL • LAC
  - LAC
- AST & CPKLIPA & AMY

• ALB & TP

PHOS & Ca

GGT & CPK

ALT & ALP

CHOL & TRIG

- UPC Kit
  - (UPRO & UCRE)

#### Integrates with Practice Management Systems







01204 669033 (Option 1)



sales@woodleyequipment.com



www.woodleyequipment.com

 Old Station Park Buildings St. John Street Horwich Bolton Lancashire BL6 7NY UK "PASSIONATE ABOUT LABORATORY DIAGNOSTICS"

