



Vcheck F (V200)

Achieve True Analysis with Fluorescence Technology



Rapid Test Kits



"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"

WOODLEY
EQUIPMENT COMPANY LTD.

BIONOTE Vcheck F (V200)

Immunofluorescence Analyser

Simplify In-Practice Point of Care Testing with Reference Lab Quality - Avoiding Test Referrals

Achieve true analysis with fluorescence technology

- Fully quantitative automated analyser
- Accurate diagnosis for Dogs, Cats and Horses in minutes
- Analyser provides print out and result storage



SPECIFICATIONS

Testing Time:	5 minutes (Based on CRP, SAA)
Weight:	2.5kg
Screen:	Touchscreen 7" Colour TFT
Size:	20 x 24 x 20 cm
Printer:	Built-in Printer
Sample Required:	Serum or Plasma (Heparin)

Canine Tests Available:

- cPL
- T4
- TSH
- Cortisol
- CRP
- Progesterone
- D-Dimer
- NT-proBNP
- Troponin I
- **Antibody Titre Testing For:** Parvovirus, Distemper & Hepatitis

Feline Tests Available:

- T4
- fPL
- SAA
- NT-proBNP
- Troponin I
- **Antibody Titre Testing For:** Calicivirus, Herpes & Parvovirus

Equine Tests Available:

- Progesterone
- SAA
- Foal IgG

**Low Cost, In-Practice,
Reference Lab Quality Immunoassay Testing**



Available exclusively from Woodley

**HIGH SENSITIVITY,
HIGH SPECIFICITY.**

Feline Tests Include:

Anigen Rapid FIV Ab/FelV Ag

Feline Immunodeficiency Virus antibody and Feline Leukemia virus antigen in Feline whole blood, plasma or serum

Anigen Rapid FeliD-3 Ag

Feline Parvo, Corona and Giardia antigen in Feline faeces

Anigen Rapid FCoV Ag/Ab

Feline Coronavirus antigen/antibody in Feline faeces (Ag) and whole blood, plasma or serum (Ab)

Anigen Rapid FHV Ag

Feline Herpes Virus antigen in Feline conjunctival

Canine Tests Include:

Anigen Rapid CPV/CCV/Giardia Ag

Canine Parvo, Corona and Giardia antigen in Canine faeces

Anigen Rapid CPV Ag

Canine Parvovirus antigen in Canine faeces

Anigen Rapid CCV Ag

Canine Coronavirus antigen in Canine faeces

Anigen Rapid CIRDC-3 Ag

Canine Infectious Respiratory Disease (CDV, CAV and CIV) in conjunctiva (CDV Ag) and nasal fluid (CAV and CIV Ag)

Anigen Rapid Rabies Ag

Rabies Virus antigen

Anigen Rapid CPV Ab 3.0

Detection and titration of Canine Parvovirus in whole blood, plasma or serum

Anigen Rapid CDV Ab 3.0

Detection and titration of Canine Distemper Virus antibody in whole blood, plasma or serum

Anigen Rapid CHW Ag

Canine Heartworm antigen in Canine serum, plasma or whole blood

Anigen Rapid Total IgE

Canine Total IgE allergy indicator test in serum

Anigen Rapid FIV Ab

Immunodeficiency Virus antibody in Feline whole blood, plasma or serum

Anigen Rapid FCoV Ab

Feline Coronavirus antibody in whole blood, plasma or serum

Anigen Rapid FPV Ag

Feline Parvovirus antigen in Feline faeces

Anigen Rapid FHW Ag/Ab

Feline Heartworm antigen and antibody in Feline whole blood, plasma or serum

Anigen Rapid FCV Ag

Feline Calicivirus antigen in Feline conjunctival

Anigen Rapid CDV Ag

Canine Distemper Virus antigen in conjunctiva (recommended), urine, serum or plasma

Anigen Rapid CDV Ag/CAV Ag

Canine Distemper Virus and Canine Adenovirus antigen in conjunctiva (CDV Ag) and nasal fluid (CAV Ag)

Anigen Rapid C.Brucella Ab

Canine Brucella antibody in whole blood, plasma or serum

Anigen Rapid Leishmania Ab

Leishmania infantum antibody in whole blood, plasma or serum

Anigen Rapid E.Canis Ab

Ehrlichia Canis antibody in whole blood, plasma or serum

Anigen Rapid Anaplasma Ab

Anaplasma phagocytophilum/Anaplasma platys antibody in whole blood plasma or serum

Anigen Rapid Lyme Ab

Borrelia burgdorferi in Canine whole blood, plasma or serum

Anigen Rapid FeLV Ag

Feline Leukemia Virus antigen in Feline whole blood, plasma or serum

Anigen Rapid FCoV Ag

For the detection of Feline Coronavirus antigen in Feline faeces

Anigen Rapid Toxoplasma Ab

Toxoplasma gondii antibody in Feline whole blood, plasma or serum

Anigen Rapid FeliCheck-3

Feline Immunodeficiency Virus antibody, Feline Leukemia Virus antigen and Feline Heartworm antibody in Feline whole blood, plasma or serum

Anigen Rapid CPV Ag/CCV Ag

Canine Parvovirus and Coronavirus antigen in Canine faeces

Anigen Rapid Giardia Ag

Giardia antigen in Canine faeces

Anigen Rapid CaniV-4

Dirofilaria immitis antigen, Ehrlichia canis antibody, Leishmania infantum antibody, Anaplasma, phagocytophilum/Anaplasma platys antibody in whole blood, plasma or serum

Anigen Rapid E.canis/

Anaplasma Ab

Ehrlichia canis and Anaplasma phagocytophilum/Anaplasma platys in whole blood, plasma or serum

Bovine Tests:

Anigen Rapid BoviD-6 Ag

Cryptosporidium antigen, Rotavirus antigen, Coronavirus antigen, Escherichia coli K99, F17, Giardia lamblia and Clostridium perfringens antigen in Bovine faeces

This list is not exhaustive - please contact us with your requirements

BIONOTE Vcheck F (V200) Immunofluorescence Analyser



LOW COST,
REFERENCE LAB QUALITY,
IMMUNOASSAY TESTING.

Anigen™ Rapid Test Kits



HIGH SENSITIVITY,
HIGH SPECIFICITY.

**Available exclusively
from Woodley**



+44 (0) 1204 669033 (Option 1)



sales@woodleyequipment.com



www.woodleyequipment.com



Old Station Park Buildings
St. John Street
Horwich
Bolton
Lancashire
BL6 7NY
UK

***"PASSIONATE ABOUT
LABORATORY DIAGNOSTICS"***

WOODLEY
EQUIPMENT COMPANY LTD.