

CAT		
PT		
Sample	qLabs	Idexx
1	15	20
2	17	21
3	15	18
4	14	16
5	17	16
6	16	18
7	15	20
8	14	19
9	16	18
10	15	18



CAT Reference Ranges
 PT qLabs 13 - 20 s
 PT Idexx 15 - 22 s

aPTT	
Sample	qLabs
1	111
2	109
3	109
4	111
5	126
6	108
7	107
8	106
9	110
10	115

DOG		
PT		
Sample	qLabs	QuickVet
1	17	17
2	17	19
3	15	17
4	14	17
5	17	19
6	15	16
7	14	17
8	18	16
9	16	15
10	17	18
11	31	30
12	15	16

DOG Reference Ranges
 PT qLabs 14 - 20 s
 PT QuickVet 14 - 19 s

aPTT	
Sample	qLabs
1	128
2	128
3	116
4	106
5	103
6	103
7	109
8	109
9	169
10	111
11	125
12	110

 Signifies result interpretation in agreement
 Signifies result interpretation not in agreement

Idexx
89
78
88
83
125
83
83
97
83
98

CAT Reference Ranges

aPTT qLabs 96 - 122 s

aPTT Idexx 65 – 119 s

QuickVet
111
159
98
105
102
115
99
108
114
101
107
97

DOG Reference Ranges

aPTT qLabs 94 - 123 s

aPTT QuickVet 75 - 105 s

?QV False prolong