



Blood Pressure Monitor *PLUS*

A portable veterinary blood pressure monitor, suitable for spot checks or continuous monitoring.



SIMPLE, ACCURATE
BLOOD PRESSURE MONITORING.

"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"

WOODLEY
EQUIPMENT COMPANY LTD.



Blood Pressure Monitor *PLUS*

A portable veterinary blood pressure monitor, suitable for spot checks or continuous monitoring.

FEATURES

- Accurate and reliable - uses 2 patented algorithms specifically developed for animals based on animal vital signs
- Easy to use - large buttons or intuitive touch screen
- Silent - measures BP quietly and alarms can be silenced to avoid distressing animals
- Fast - take multiple measurements quickly
- Memory - 1000 results in Monitor Mode, 100 patients and 200 Data sets per patient in Spot Mode
- Portable - rechargeable lithium or alkaline battery operation
- Data transfer - APP available to transfer data to mobile phone or Tablet/PC via Bluetooth
- Includes real-time waveform display on screen

SPECIFICATIONS

Method of Measurement:	Schematic method
Blood Pressure Range:	SYS: 40-270 mmHg DIA: 20-230 mmHg MEAN: 10-200 mmHg
Pulse Rate Range:	BPM 40-240
Cuff Accuracy:	0-300 mmHg (± 3 mmHg)
Conditions of Use:	0-50°C, 15%-90% humidity, no condensation
Storage Conditions:	-20-65°C, 15%-90% humidity, no condensation
Power:	Lithium ion batteries and/or AC adapter
Input Voltage Range:	Universal (100-240V AC)
Input Frequency Range:	General purpose (50-60Hz)
Adapter Output:	DC: 9 V-2A
Lithium Batteries:	7.4V, 2.6Ah, 19.24Wh
Dimensions:	15.9 x 12.7 x 13.3 cm
Weight:	0.95 kg (with battery)



Includes Protective Sleeve in choice of colours - blue, green or pink

Includes Real-Time Waveform Display On Screen



+44 (0) 1204 669033 (Option 1)



sales@woodleyequipment.com



www.woodleyequipment.com



Old Station Park Buildings
St. John Street
Horwich
Bolton
Lancashire
BL6 7NY
UK

"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"

**WOODLEY
EQUIPMENT COMPANY LTD.**