

# **InSight<sup>®</sup> &** **InSight<sup>®</sup> BP<sup>PLUS</sup>**

## **Blood Pressure Monitors**

**Portable veterinary blood pressure monitors,  
suitable for spot checks or continuous monitoring.**



**SIMPLE, ACCURATE**  
**BLOOD PRESSURE MONITORING.**

# InSight® & InSight® BP PLUS

## Blood Pressure Monitors

Portable veterinary blood pressure monitors,  
suitable for spot checks or continuous monitoring.

### FEATURES

- Accurate and reliable - uses 2 patented algorithms specifically developed for animals based on animal vital signs
- Easy to use - large buttons or intuitive touch screen
- Silent - measures BP quietly and alarms can be silenced to avoid distressing animals
- Fast - take multiple measurements quickly
- Memory - 1000 results in Monitor Mode, 100 patients and 200 Data sets per patient in Spot Mode
- Portable - rechargeable lithium or alkaline battery operation
- Data transfer - APP available to transfer data to mobile phone or Tablet/PC via Bluetooth
- InSight BP Plus includes real-time waveform display on screen



Includes Protective Sleeve  
in choice of colours -  
blue, green or pink

### SPECIFICATIONS

<b>Method of Measurement:</b>	Schematic method
<b>Blood Pressure Range:</b>	SYS: 40-270 mmHg DIA: 20-230 mmHg MEAN: 10-200 mmHg
<b>Pulse Rate Range:</b>	BPM 40-240
<b>Cuff Accuracy:</b>	0-300 mmHg ( $\pm 3$ mmHg)
<b>Conditions of Use:</b>	0-50°C, 15%-90% humidity, no condensation
<b>Storage Conditions:</b>	-20-65°C, 15%-90% humidity, no condensation
<b>Power:</b>	Lithium ion batteries and/or AC adapter
<b>Input Voltage Range:</b>	Universal (100-240V AC)
<b>Input Frequency Range:</b>	General purpose (50-60Hz)
<b>Adapter Output:</b>	DC: 9 V-2A
<b>Lithium Batteries:</b>	7.4V, 2.6Ah, 19.24Wh
<b>Dimensions:</b>	15.9 x 12.7 x 13.3 cm
<b>Weight:</b>	0.95 kg (with battery)

**InSight BP Plus Includes  
Real-Time Waveform  
Display On Screen**



+44 (0) 1204 669033 (Option 1)



sales@woodleyequipment.com



www.woodleyequipment.com



Old Station Park Buildings  
St. John Street  
Horwich  
Bolton  
Lancashire  
BL6 7NY  
UK

**"PASSIONATE ABOUT  
LABORATORY DIAGNOSTICS"**

**WOODLEY  
EQUIPMENT COMPANY LTD.**