

InSight[®] Equine-IA

Immunoassay Analyser

A portable, handheld, easy to use Equine Immunoassay Analyser providing quick, accurate and reliable results.



**LOW COST, PORTABLE,
EQUINE IMMUNOASSAY TESTING**

eACTH, eSAA, eProg, Foal IgG, Insulin



"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"

WOODLEY
EQUIPMENT COMPANY LTD.

Immunoassay Analyser

The InSight Equine-IA is a portable, handheld, easy to use Equine Immunoassay Analyser providing quick, accurate and reliable results.

FEATURES

- Immunofluorescence technology
- Testing time 3-15 minutes (depending on analyte)
- Portable, with built-in rechargeable battery - ready to use anytime, anywhere
- Easy to use with touch screen technology
- Wi-Fi and Bluetooth connectivity - easy transfer of data
- Accurate and reliable
- Multi-language operation system
- Wi-Fi connectivity
- Reagents stored at room temperature
- Touch screen
- Maintenance free
- Results transmitted via InSight p-LINK to Practice Management Systems



SPECIFICATIONS

Operating System:	Built-in Android test management software
Sample Type:	Serum and Plasma (EDTA)
Data Storage:	ROM Capacity: 8GB RAM Capacity: 1GB
Built-in Battery:	Removable rechargeable battery Standby Time: >72 hours Operating Time: >4 hours
Data Output:	USB data exporting Built-in thermal printer Data Transfer
Wireless Communication:	4G, Wi-Fi, Bluetooth
Image Processing:	Built-in 500W autofocus camera
Working Environment:	Temperature: +10°C to +30°C Humidity: ≤80% RH
Input Voltage:	DC 12V/5A
Dimensions:	240 (W) x 90 (D) x 65 (H) mm
Weight:	0.8kg

TESTS AVAILABLE

- Equine ACTH
- Equine SAA
- Equine Progesterone
- Foal IgG

TESTS IN DEVELOPMENT

- Equine Insulin - launching later this year

Software Languages

English, Chinese, Spanish, French, Russian, Portuguese, Turkish, Italian, Polish, German



+44 (0) 1204 669033 (Option 1)



sales@woodleyequipment.com



www.woodleyequipment.com



Old Station Park Buildings
St. John Street
Horwich
Bolton
Lancashire
BL6 7NY
UK

**"PASSIONATE ABOUT
LABORATORY DIAGNOSTICS"**

WOODLEY
EQUIPMENT COMPANY LTD.