

# **InSight V-IA<sup>®</sup>**

## **Veterinary Immunoassay Analyser**

**An easy to use Immunoassay Analyser providing accurate and reliable results in 3-15 minutes.**



**LOW COST,  
REFERENCE LAB QUALITY,  
IMMUNOASSAY TESTING.**



**"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"**

**WOODLEY  
EQUIPMENT COMPANY LTD.**

An easy to use Immunoassay Analyser providing accurate and reliable results in 3-15 minutes.

## FEATURES

- Immunofluorescence technology
- Testing time 3-15 minutes (depending on analyte)
- Easy to use
- Accurate and reliable
- Multi-language operation system
- Reagents stored at room temperature
- Touch screen
- Maintenance free
- Results transmitted via InSight p-LINK to Practice Management Systems



## SPECIFICATIONS

- Dimensions:** 213 (W) x 303 (D) x 210 (H) mm  
**Weight:** 4kg  
**Power Adapter:** 100~240V, 50~60Hz

### Canine Tests Available:

- cPL
  - T4
  - TSH
  - Cortisol
  - Progesterone
  - NT-proBNP/Troponin I (combined test)
  - Cystatin C (SDMA alternative)
  - cCRP
  - D-Dimer
  - HbA1C
  - NT-proBNP
  - Troponin I
- **Core Antibody:** Distemper, Parvovirus, Distemper/ Parvovirus/Infectious Hepatitis Virus combined test
- **Infectious Disease:** Coronavirus, Distemper, Parvovirus, Parvovirus/ Coronavirus combined test

### Feline Tests Available:

- fPL
  - T4
  - TSH
  - Cortisol
  - Cystatin C (SDMA alternative)
  - D-Dimer
  - fSAA
  - NT-proBNP
  - Troponin I
- **Core Antibody:** Calicivirus, Herpes, Panleukopenia Virus, Herpes/Calicivirus/ Panleukopenia Virus combined test
- **Infectious Disease:** Coronavirus, Calicivirus, Herpes, Panleukopenia Virus

### Equine Tests Available:

- eSAA

### Tests In Development:

- fTT4 (Feline Total Thyroxine) Quantitative Test
- HP (Helicobacter Pylori) Quantitative Test
- CPSE
- TOXO Ab
- RB virus Ab
- AFP (Alpha-fetoprotein)
- Fructosamine
- Equine ACTH
- Equine Progesterone
- Foal IgG



+44 (0) 1204 669033 (Option 1)



sales@woodleyequipment.com



www.woodleyequipment.com



Old Station Park Buildings  
St. John Street  
Horwich  
Bolton  
Lancashire  
BL6 7NY  
UK

**"PASSIONATE ABOUT  
LABORATORY DIAGNOSTICS"**

**WOODLEY  
EQUIPMENT COMPANY LTD.**