

# **InSight V-IA<sup>®</sup> PLUS**

## **Veterinary Immunoassay Analyser**

**An easy to use Immunoassay Analyser providing accurate and reliable results in 3-15 minutes.**



**LOW COST,  
REFERENCE LAB QUALITY,  
IMMUNOASSAY TESTING.**



**"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"**

**WOODLEY  
EQUIPMENT COMPANY LTD.**

# InSight V-IA<sup>PLUS</sup> Veterinary Immunoassay Analyser

An easy to use Immunoassay Analyser providing accurate and reliable results in 3-15 minutes

## FEATURES

- Immunofluorescence technology
- Testing time 3-15 minutes (depending on analyte)
- Easy to use
- Accurate and reliable
- Multi-language operation system
- Wi-Fi connectivity
- Reagents stored at room temperature
- Touch screen
- Maintenance free
- Results transmitted via InSight p-LINK to Practice Management Systems



Scan QR Code for more product information and resources

## Canine Tests Available:

- cPL
- T4
- TSH
- Cortisol
- Progesterone
- Alpha-fetoprotein
- Relaxin
- Distemper
- Parvovirus
- Rabies Virus
- Ehrlichia
- Toxoplasma Gondii
- Parvovirus/Coronavirus combined test
- Cystatin C (SDMA alternative)
- Anaplasma Phagocytophilum
- NT-proBNP/Troponin I combined test
- Distemper/Parvovirus/Infectious Hepatitis Virus combined test
- cCRP
- D-Dimer
- HbA1c
- NT-proBNP
- Troponin I
- CPSE
- HP (Helicobacter Pylori)
- Lyme Disease
- Coronavirus
- Giardia
- Rotavirus

## Feline Tests Available:

- fPL
- fTT4
- TSH
- Cortisol
- Troponin I
- Calicivirus
- Herpes
- Rabies Virus
- Toxoplasma Gondii
- NT-proBNP/Troponin I combined test
- Cystatin C (SDMA alternative)
- Herpes/Calicivirus/Panleukopenia Virus combined test
- fSAA
- D-Dimer
- HbA1c
- NT-proBNP
- HP (Helicobacter Pylori)
- Coronavirus
- Panleukopenia Virus
- Giardia
- Rotavirus

## Equine Tests Available:

- eSAA
- Progesterone
- Foal IgG

## Software Languages

English, Chinese, Spanish, French, Russian, Portuguese, Turkish, Italian, Polish, German

## SPECIFICATIONS

**Dimensions:** 213 (W) x 243 (D) x 195 (H) mm

**Weight:** 4kg

**Power Adapter:** 100~240V, 50~60 Hz

## Tests In Development:

- Fructosamine
- Equine ACTH
- FIV Ab
- FeLV Ag



+44 (0) 1204 669033 (Option 1)



sales@woodleyequipment.com



www.woodleyequipment.com



Old Station Park Buildings  
St. John Street  
Horwich  
Bolton  
Lancashire  
BL6 7NY  
UK

**"PASSIONATE ABOUT  
LABORATORY DIAGNOSTICS"**

**WOODLEY  
EQUIPMENT COMPANY LTD.**