

InSight V-IA[®] Testing Information

Canine Cystatin C				
Species	Sample Volume	Dilution Volume	Test Time	Sample Type
Canine	10µl	75µl	10 mins	Serum or Plasma (EDTA)

Feline Cystatin C				
Species	Sample Volume	Dilution Volume	Test Time	Sample Type
Feline	75µl	75µl	10 mins	Serum or Plasma (EDTA)

Canine NT-proBNP				
Species	Sample Volume	Dilution Volume	Test Time	Sample Type
Canine	75µl	75µl	15 mins	Serum or Plasma (Lithium Heparin)

Feline NT-proBNP				
Species	Sample Volume	Dilution Volume	Test Time	Sample Type
Feline	75µl	75µl	15 mins	Serum or Plasma (EDTA)

cPL				
Species	Sample Volume	Dilution Volume	Test Time	Sample Type
Canine	75µl	75µl	15 mins	Serum or Plasma (Lithium Heparin)

fPL				
Species	Sample Volume	Dilution Volume	Test Time	Sample Type
Feline	75µl	75µl	15 mins	Serum or Plasma (Lithium Heparin)

CRP				
Species	Sample Volume	Dilution Volume	Test Time	Sample Type
Canine	10µl	75µl	3 mins	Serum or Plasma (Lithium Heparin)

Feline Troponin I				
Species	Sample Volume	Dilution Volume	Test Time	Sample Type
Feline	75µl	75µl	15 mins	Serum or Plasma (Lithium Heparin)

Cortisol				
Species	Sample Volume	Dilution Volume	Test Time	Sample Type
Canine or Feline	75µl	75µl	15 mins	Serum (Plain) or Plasma (Lithium Heparin)

fSAA				
Species	Sample Volume	Dilution Volume	Test Time	Sample Type
Feline	10µl	75µl	3 mins	Serum or Plasma (Lithium Heparin)

Progesterone				
Species	Sample Volume	Dilution Volume	Test Time	Sample Type
Canine	75µl	75µl	15 mins	Serum (Plain) or Plasma (Lithium Heparin)

HbA1C				
Species	Sample Volume	Dilution Volume	Test Time	Sample Type
Canine	10µl	75µl	5 mins	Whole Blood (EDTA)

D-Dimer				
Species	Sample Volume	Dilution Volume	Test Time	Sample Type
Canine and Feline	10µl	75µl	5 mins	Plasma (Lithium Heparin)

T4				
Species	Sample Volume	Dilution Volume	Test Time	Sample Type
Canine and Feline	75µl	75µl	15 mins	Serum (Plain) or Plasma (Lithium Heparin)

eSAA				
Species	Sample Volume	Dilution Volume	Test Time	Sample Type
Equine	10µl	75µl	5 mins	Serum (Plain) or Plasma (Lithium Heparin)

TSH				
Species	Sample Volume	Dilution Volume	Test Time	Sample Type
Canine and Feline	75µl	75µl	15 mins	Serum (Plain) or Plasma (Lithium Heparin)

P/S-3232 Rev: 1

Tel: +44 (0) 1204 699033 (Option 1)

"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"

Email: sales@woodleyequipment.com

Web: www.woodleyequipment.com

WOODLEY
EQUIPMENT COMPANY LTD.