

[www.woodleyvetdiagnostics.com](http://www.woodleyvetdiagnostics.com)

**ANNOUNCING**

# THE MYTHIC 18 VET

**:: THE ULTIMATE TECHNOLOGY IN HAEMATOLOGY ::**



**FULLY AUTOMATED 18 PARAMETER  
HAEMATOLOGY ANALYSER DEDICATED  
FOR VET APPLICATIONS**



Swiss Made

- Low reagent consumption
- Low maintenance
- Less servicing required
- Fully automated

- Light weight
- Low power consumption
- Low blood sample volume
- Quick and easy to use

Tel: +44(0)1204 669033

Fax: +44(0)1204 669034

Email: [sales@woodleyvetdiagnostics.com](mailto:sales@woodleyvetdiagnostics.com)


Web: [www.woodleyvetdiagnostics.com](http://www.woodleyvetdiagnostics.com)

**WOODLEY**  
**VETERINARY DIAGNOSTICS**

OLD STATION PARK BUILDINGS, ST. JOHNS STREET, HORWICH, BOLTON, LANCASHIRE, BL6 7NY, UK

# THE MYTHIC 18 VET HAEMATOTOLOGY

**EASY TO USE.....EASY TO SERVICE.....**



## QVGA Touch Screen Portrait

Position identical to print out for easy reading

## Built-In Reagent Tray

Bench space saving and easy to replace

## Only 3 Hydraulic Modules

Cost saving and easy to replace

## Needle and Chambers Dismountable without tools

Simple and easy maintenance

**Parameters** 18 parameters + 3 curves: WBC, LYM% and #, MON% and #, GRA% and #, EOS flag, RBC, HGB, HCT, MCV, MCH, MCHC, RDW, PLT, MPV, PCT, PDW, RBC curve, WBC curve, PLT curve

**Throughput** 60 samples per hour

**Sample Volume** < 10 µl

**Output** RS232, Ethernet (UDP/IP, TCP/IP protocol) 2 USB ports, Centronics port, Barcode scanner port

**Reagents** Cyanide free (up to 165 cycles per pack)

**Printers** External printer Dot matrix printer / thermal printer / monochrome injet printer / monochrome laser printer

**Power (voltage & consumption)** 90 to 250 V, 50 to 60 Hz, Less 50 VA

**Graphic Interface** Menu with icons for direct access to main species

**Technology** Impedance (WBC, RBC, PLT) Spectrophotometry (HGB)

**User Interface** TFT colour touchscreen  
Built-in numerical keyboard

10 Species Available

## Calibration

Flags

## Quality Control

Memory

## Multi-language

Traceability

Dog, Cat, Horse, Cow, Pig, Sheep, Goat, Rabbit, Mouse, Rat.  
10 other species can be programmed by the user

## Automatic or manual

Pathological flags (programmable)  
Lab limits (programmable)  
Reagent alerts  
Instrument alerts

## Liquid QC available Table and Levey Jennings graphs

1,500 patient internal memory  
Mass storage USB up to 60,000 patients

**French, English, Spanish, Portuguese, Italian, German, Chinese, Russian, Polish, Serbian, Hungarian**

Patient identification: name, patient ID, sample ID  
Daily log  
Reagent monitoring: lot number, expiry date, volume  
Service log

Reagents - M-Pack for Mythic 18 Vet  
Options - external bar code scanner

**Tel:** +44(0)1204 669033

**Fax:** +44(0)1204 669034

**Email:** [sales@woodleyvetdiagnostics.com](mailto:sales@woodleyvetdiagnostics.com)

**Web:** [www.woodleyvetdiagnostics.com](http://www.woodleyvetdiagnostics.com)

OLD STATION PARK BUILDINGS, ST. JOHNS STREET, HORWICH, BOLTON, LANCASHIRE, BL6 7NY, UK

**WOODLEY**  
**VETERINARY DIAGNOSTICS**