



# skyla VB1

## Veterinary Chemistry Analyser

For both Rotor Profiles and Single/Dual Tests



## FAST, ACCURATE CHEMISTRY

## PANELS & SINGLE TESTS.

**5 YEAR WARRANTY.**



**"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"**

**WOODLEY**  
EQUIPMENT COMPANY LTD.

# Veterinary Chemistry Analyser

The skyla VB1 Veterinary Chemistry Analyser is a fast, compact and accurate veterinary chemistry analyser ideal for any vet practice or small laboratory.

**Highly accurate test results in under 15 minutes**

## VB1 Panel Test Portfolio

Analytes Profile	Items	ALP	ALT	AST	BUN	CREA	GLU	ALB	TP	TBIL	AMY	Ca	PHOS	Na	K	Cl	BA	GGT	CPK	CHOL	TRIG	tCO2	Others
Preanesthetic Panel	7	ALP	ALT		BUN	CREA	GLU	ALB	TP														
Renal Panel	9				BUN	CREA		ALB				Ca	PHOS	Na	K	Cl						tCO2	
Diagnosis Panel	13	ALP	ALT	AST	BUN	CREA	GLU	ALB	TP	TBIL	AMY	Ca	PHOS							CHOL			
Diagnosis Plus Panel	13	ALP	ALT		BUN	CREA	GLU	ALB	TP	TBIL	AMY	Ca	PHOS	Na	K								
Equine Panel	14	ALP		AST	BUN	CREA	GLU	ALB	TP	TBIL		Ca		Na	K		GGT	CPK				tCO2	
Thyroid/T4 Panel	2																			CHOL			T4
Large Animal Panel	14	ALP		AST	BUN		GLU	ALB	TP			Ca	PHOS	Na	K	Cl	GGT	CPK					Mg
Liver Plus Panel	10	ALP	ALT	AST	BUN			ALB	TP	TBIL							BA	GGT		CHOL			
Avian & Reptile Panel	12			AST			GLU	ALB	TP			Ca	PHOS	Na	K	Cl	BA		CPK				UA
Diabetes Panel	4						GLU	ALB	TP														FRU
Electrolytes Panel	4													Na	K	Cl						tCO2	
Diagnosis-II Panel	14	ALP	ALT		BUN	CREA	GLU	ALB	TP	TBIL	AMY	Ca	PHOS					GGT		CHOL			LIPA
Critical Care Panel	14	ALP	ALT		BUN	CREA	GLU	ALB	TP			Ca		Na	K	Cl						tCO2	LAC, CPK

## FEATURES

- Compact, benchtop size
- Performs up to 15 biochemistry tests simultaneously
- Performs up to 3 single/dual tests simultaneously
- Highly accurate test results in under 15 minutes
- Precision equal to larger laboratory analyser
- Easy and quick operation
- Integrates to Practice Management Systems
- LCD colour touch screen
- Stores up to 50,000 results
- Species specific reference ranges
- No maintenance required
- **5 year warranty**

## SINGLE/DUAL PARAMETERS AVAILABLE

- |                    |                                    |               |
|--------------------|------------------------------------|---------------|
| <b>Single Test</b> | <b>Dual Test</b>                   |               |
| • CREA             | • Na <sup>+</sup> & K <sup>+</sup> | • BA & CHOL   |
| • BUN              | • PHOS & Ca                        | • AST & CPK   |
| • TBIL             | • CHOL + TRIG                      | • LIPA & AMY  |
| • LAC              | • GGT & CPK                        | • UPC Kit     |
| • Ammonia          | • ALT & ALP                        | (UPRO & UCRE) |

## FULLY AUTOMATIC TEST

- Built-in printer
- Automatic system with walk away operation
- Built-in QC and calibration for each test
- Rotors are used straight from the fridge

## NO SAMPLE PRE-TREATMENT

- Single-use reagent discs
- Panels and single tests available
- 200µl whole blood, serum or plasma



1. Add sample



2. Place in analyser



3. Read results



+44 (0) 1204 669033 (Option 1)



sales@woodleyequipment.com



www.woodleyequipment.com



Old Station Park Buildings  
St. John Street  
Horwich  
Bolton  
Lancashire  
BL6 7NY  
UK

**"PASSIONATE ABOUT  
LABORATORY DIAGNOSTICS"**

**WOODLEY  
EQUIPMENT COMPANY LTD.**