


Vcheck Analyzer

Fluorescent immunoassay for quantitative results

EUROPIUM, A POWERFUL FLUORESCENT MATERIAL

Vcheck Analyzer utilizes Europium, which has exceptional sensitivity over other general fluorescence materials.

Europium, which continues to fluoresce even when stopping the light, provides precise results without any interference of the noise signals.



Quantitative Test

Antibody Titer Test

Infectious Test

* Biomarkers under development



Inflammatory
CRP, SAA



Pancreatitis
cPL, fPL



Hormone
Cortisol, T4, TSH
Progesterone



Ab Titers
CPV, CDV, CAV Ab
FPV, FHV, FCV Ab



Coagulation
D-dimer



Cardiac
NT-proBNP (Canine, Feline)
TnI (Canine, Feline*)



Tumor
TK1*, AFP*

In-House Quantitative
Diagnosis Possible

Up to 20+ Biomarkers
Available

User Friendly
3-Step Procedures

Ordering Information

Product No.	Product Name	Product Type	Packing Unit
VC7402EA	Vcheck V200	Equipment	1 EA
VC7403EA	Vcheck V2400	Equipment	1 EA



Manufactured by

BioNote, Inc.

22 Samsung1ro 4-gil, Hwaseong-si, Gyeonggi-do 18449, Republic of Korea
TEL: 82-31-211-0516 | FAX: 82-31-8003-0618 | www.bionote.co.kr