

3301 Pulse Oximeter

Available for hire or purchase for your clinical trial.



Overview

The 3301 pulse oximeter is a lightweight handheld device delivering pulse oximetry wherever spot-checking is needed.

It has an easy to read display and is ergonomically designed to fit in your hand. It delivers fast reliable measurements of SpO2 and pulse rate.

Ideal for use in a variety of settings it accurately delivers measurements from neonatal to adult patients.

Operating Features

- Durable and lightweight
- Ergonomic design
- Stores data for 99 different patients
- Simple on/off control with auto shut off to preserve battery life
- Quantitative pulse strength bar on screen
- Monitor patients from neonatal to adult
- Bright, easy to read display, even in low light settings
- Reliable for use in various locations
- Durable power and conveniently portable

Specification

- Pulse Strength: logarithmically scaled, 8-segment bar graph
- SpO2: range 0 - 99%, accuracy $\pm 2\%$ at 70 - 99%, averaging 8 beat
- Pulse Rate: range 30 - 254 bpm, accuracy $\pm 2\%$ at 30 - 254 bpm, averaging 8 second
- Memory: 1,000 spot checks, up to 99 patients
- Display: LED numeric display
- Power Requirements: 3 x "C" cell alkaline or NiCAD batteries
- Battery Life: approximately 24 hours in continuous use mode - 1,500 spot checks
- Operating Temperature: 0°C to +55°C
- Storage Temperature: -34°C to +70°C
- Relative Humidity: 10 - 95% RH (storage)
- Weight: 462g (with batteries)
- Dimensions: 160h x 82w x 31d mm



Please contact us for a quotation for the hire or purchase of this product

Product description: - - 3301 Pulse Oximeter

Product code - WD7790

Product specifications can be subject to change without prior notice. Version: V2 10/19

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company