

3303 Pulse Oximeter

Available for hire or purchase for your clinical trial.



Overview

The 3303 pulse oximeter is cost-effective, portable and light-weight. It is a continuous or spot check pulse oximetry monitoring system with alarms and an internal rechargeable battery.

The 3303 pulse oximeter provides fast, reliable SpO₂ and pulse rate measurements on patients from neonatal to adult.

Stores data for later sleep screening analysis via a PC or commercially available printer.

Operating Features

- Portable and lightweight
- Spot checking or continuous use
- 3 operating modes - clinic, home or sleep
- Menu free operations using a one touch, direct function key
- Adjustable LED display brightness
- Pulse strength bar on screen
- Adjustable alarm limits for SpO₂ and pulse rate
- Auxiliary serial printer/computer output - data saved every 30 seconds can be downloaded or printed
- Distinguishable audible patient and system alarms
- On-screen Patient History: recall up to 99 patients readings
- 24 hour patient reading memory
- Protective rubber boot included
- Pole mount available

Specification

- Pulse Strength: logarithmically scaled, 8-segment LED bar graph
- SpO₂: range: 0 - 100%, accuracy: $\pm 2\%$ at 70 - 100%, averaging: 4, 8 or 16 pulse beat, alarm range: high 50 - 100% or OFF, low 50 - 99% or OFF
- Pulse Rate: range: 30 - 254 bpm, accuracy: $\pm 2\%$ at 30 - 254, resolution: 1 bpm, averaging: 8 or 16 sec, alarm: high 5 - 250 bpm or OFF, low 5 - 250 bpm or off
- Memory: up to 99 patients and 24 hours of run time
- Display: LED numeric display
- Power Requirements: 207 - 253 VAC, 50/60 Hz or battery: standard NiMH
- Battery Life: internal rechargeable not user replaceable - charges in 6 hours for 24 hours continuous use
- Weight: 539g (with batteries)
- Dimensions: 184h x 84w x 47d mm



Please contact us for a quotation for the hire or purchase of this product

Product description: -3303 Pulse Oximeter

Product code - WD7793

Product specifications can be subject to change without prior notice. Version: V2 10/19

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company