

3402 Pulse Oximeter

Available for hire or purchase for your clinical trial.



Overview

The 3402 pulse oximeter quickly and reliably delivers pulse oximetry wherever spot-checking or continuous monitoring is required - ideal for bedside monitoring as well as during sleep studies.

It includes pulse rate, pulse strength and SpO₂ - with a built in infrared printer output port.

The 3402 pulse oximeter has Serial Autocorrelation (SA) to inform the operator of excess artifacts and improper sensor positioning.

Operating Features

- SpO₂, pulse rate and pulse strength
- Infrared link to printer
- Serial Autocorrelation (SAC) technology informing of excess artifacts and improper sensor positioning
- Portable and lightweight
- LED numeric display
- Spot checking or continuous use
- 3 Operating Modes: clinic, home or sleep
- Easy to use 5 key operation
- Battery powered or battery bypass
- Pulse strength bar on screen
- Adjustable alarm limits
- Protective boot available

Specification

- Pulse Strength: logarithmically scaled, 8-segment LED bar graph
- SpO₂: range: 0 - 99% functional SpO₂, accuracy: $\pm 2\%$ at 70 - 99%, averaging: 4, 8 or 16 pulse beat, alarm range: high - 50 - 99% or OFF, low 50 - 99% or OFF
- Pulse Rate: range: 30 - 254 bpm, accuracy: $\pm 2\%$ at 30 - 254 bpm, alarm range: high - 5 - 250 bpm or OFF, low - 5 - 250 bpm or OFF, averaging: 8 / 16 sec
- Trend Storage: pulse rate, pulse strength, SpO₂ and relative time: 4 - 30 second intervals
- Memory: up to 90 hours collected from up to 99 patients
- Printer Output: data log mode or trend data mode
- Power Requirements: AC power 12VDC supply or 6 standard AA alkaline batteries
- Battery Life: 24 hours without printing
- Weight: 454g (with batteries)
- Dimensions: 167h x 70w x 36d mm



Please contact us for a quotation for the hire or purchase of this product

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

Product description: -3402 Pulse Oximeter

Product code - WD7792

Product specifications can be subject to change without prior notice. Version: V2 10/19

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company