

## ► Eco Temp Basic Digital Thermometer

*available for hire or purchase for your clinical trial.*

### Overview

The Eco Temp Basic Digital Thermometer has a round tip designed for the patient to use safely and conveniently during all measurement methods.

The Eco Temp Basic is water resistant and has an exchangeable battery. The on/off button is easy to use.

For oral, axillary and rectal use.

### Specification

- Measurement Range: 32°C to 42°C
- Accuracy of measurement: +/- 0.1°C
- Power Supply: 1.5V DC, 1 LR41 Alkaline-Magnesium Button Battery
- Power Consumption: 0.1 mW
- Battery Life: Approximately 2 years or more (3 times a day)
- Temperature Display: 3 digit, +°C display in 0.1 degree increments
- Operating Temperature: +10°C to +40°C
- Relative Humidity: 30% to 85% RH
- Weight: approximately 12g (including battery)
- Dimensions: 19.4 x 132.5 x 10 mm

### Operating Features

- Digital fever thermometer for oral, axillary and rectal use
- Normal response in 60 seconds
- Water resistant
- 3 digit display (1 decimal place)
- Reading display in Celsius or Fahrenheit
- Beeper signal when measurement is finished
- Memory for last reading
- Automatic switch off saves on battery-life
- Supplied with a protective case
- Exchangeable battery



*Please contact us for a quotation for the hire or purchase of this product.*

**Product description** - Eco Temp Basic Digital Thermometer  
**Product code** -



#### UK Office

T: +44(0)8456 777001

F: +44(0)8456 777002

E: [sales@woodleyequipment.com](mailto:sales@woodleyequipment.com)

#### USA Office

T: 1 800 471 9200 / 1 508 625 1693

F: 1 508 625 1721

[www.woodleyequipment.com](http://www.woodleyequipment.com)

**WOODLEY**  
EQUIPMENT COMPANY LTD.