

## ► Microman E Pipette 100 – 1000ul

*available for hire or purchase for your clinical trial.*

### Overview

Microman E pipette is the ideal solution to improve accuracy and precision when pipetting hot or cold liquids and is also ideal for pipetting viscous and dense samples. The pipette works like a syringe and capillary pistons eliminate the air cushion between the sample and the disposable piston, providing 100% protection against pipette contamination and producing accurate, reproducible, and precise results every time. The clip fitting system allows for a capillary piston to fit on to Microman E as simply as fitting a standard tip to a standard pipette.

### Specification

- Volume (µL):           Min. 100   500   Max. 1000
- Systematic Error (µL):     ± 3     ± 5     ± 8
- Random Error (µL):       ≤ 1.6   ≤ 2.5   ≤ 4



### Operating Features

- The volume securing button allows for easy volume adjustment and minimizes the risk of accidentally changing the volume during pipetting cycles.
- The Microman E pipette is comfortable and easy to use.
- The shape of the handle, with an added finger hook, rests comfortably in your hand.
- Fully visible volumeter on the front of the pipette so that the volume can be seen while pipetting.
- A large push button allows you to aspirate and dispense with ease.

*Please contact us for a quotation for the hire or purchase of this product.*

**Product description** -Microman E Pipette 100 – 1000ul

**Product Code:** WD7LAB009



#### UK Office

T: +44(0)8456 777001

F: +44(0)8456 777002

E: enquiries@woodleyequipment.com

#### USA Office

T: 1 800 471 9200 /1 718 606 0516

F: 1 718 606 0955

[www.woodleyequipment.com](http://www.woodleyequipment.com)

**WOODLEY**  
 EQUIPMENT COMPANY LTD.