

► Gilson Pipette Classic P5000 1-5ml



available for hire or purchase for your clinical trial.

Overview

The Gilson Pipette Classic P5000 model is robust, accurate and precise pipette. It is ideal for large-volume dispensing of aqueous fluids of moderate viscosity and density. This model is equipped with a polyethylene filter that prolongs piston life and prevents contamination.

Gilson Pipette Classic P5000 model is a fully adjustable, air-displacement pipette with the selected volume shown on the digital indicator (Volumeter).

Please note that this model does not have a tip ejector.

Specification

- Volume (μL): Min. 1000 2000 Max. 5000
- Systematic Error (μL): ± 12 , ± 12 ± 30
- Random Error (μL): ≤ 3 ≤ 5 ≤ 8



Please contact us for a quotation for the hire or purchase of this product.

Product description - Gilson Pipette Classic P5000 1-5ml

Product Code: WD7MISC195 (US)



UK Office

T: +44(0)8456 777001

F: +44(0)8456 777002

E: enquiries@woodleyequipment.com

USA Office

T: 1 800 471 9200 / 1 718 606 0516

F: 1 718 606 0955

www.woodleyequipment.com

WOODLEY
EQUIPMENT COMPANY LTD.