

-100°C Digital Thermometer

Available for hire or purchase for your clinical trial.



Overview

This -100°C digital thermometer has a wide temperature range with Traceable® Platinum Freezer technology. With a temperature range of -99.9°C to +199.9°C this digital thermometer can be used in a wide range of situations.

It monitors temperatures in freezers, water baths, heating blocks, incubators and refrigerators.

Operating Features

- RTD platinum probe
- 3m ultra-thin micro-cable permits freezer doors to open and close on it without affected the seal
- Fast response 100-Ω platinum probe is triple purpose for use with liquids, air/gas and frozen material
- High-impact, chemical resistant ABS case
- Adjustable probe holder, Velcro™, and magnetic strips are provided for mounting of both thermometer and probe
- Supplied with probe, 9V battery, carry case and Traceable® Certificate ISO 17025 A2LA to NIST

Specification

- Range: -99.9°C to +199.9°C
- Resolution: 0.1°C
- Accuracy: ± 2°C
- RTD platinum probe: 3.2 x 127mm
- Weight: 120g
- Dimensions: 108h x 70w x 19d mm



Please contact us for a quotation for the hire or purchase of this product

Product description: -100°C Digital Thermometer

Product Code: -

Product specifications can be subject to change without prior notice. Version: V2 10/19

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company