

-121 °F Chest Freezer, 2.6 cu.ft

Available for hire or purchase for your clinical trial.



Overview

This is one of the most simple ultra low temperature freezers on the market.

With a single compressor the -121°F chest freezer provides low noise, low energy consumption and a compact design.

The -121°F chest freezer demonstrates high performance freezing with an integrated cooling system.

Operating Features

- Single compressor
- Digital display
- Castors and key lock
- Low heat dissipation
- Low energy consumption
- Low noise level
- Visual and acoustic alarms -
 - power failure
 - adjustable high/low temperature
 - probe failure
 - Instrument failure alarm
- Contact for remote alarm
- Prepared for GSM alarm module
- Advanced battery back-up for alarms and temperature display (72 hours)

Specification

- Capacity: 2.6 cu. ft
- Temperature Range: -121°F to -40°F
- Noise Level: 50 dBA
- Power Supply: 120V, 60Hz
- Power Consumption: 7.7 KWH/24H
- Weight: 108lbs
- External Dimensions: 33.5h x 21.7w x 25.5d inches



Please contact us for a quotation for the hire or purchase of this product

Product description: -121F, 2.6 cu. ft Chest Freeze

Product code: - WD7ULT010

Version: V2 10/19

Product specifications can be subject to change without prior notice.

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company