

-121 °F/-85 °C Chest Freezer, 368 Litre

Available for hire or purchase for your clinical trial.



Overview

This is one of the most simple ultra low temperature freezers on the market.

With a single compressor the -121 °F / -85°C, 368 litre chest freezer provides low noise, low energy consumption and a compact design.

The -121 °F / -85°C chest freezer demonstrates high performance freezing with an integrated cooling system.

Operating Features

- Single compressor
- Digital display
- Castors and key lock
- Low heat dissipation
- Low energy consumption
- Low noise level
- Visual and acoustic alarms -
 - power failure
 - adjustable high/low temperature
 - probe failure
 - Instrument failure alarm
- Contact for remote alarm
- Prepared for GSM alarm module
- Advanced battery back-up for alarms and temperature display (72 hours)

Specification

- Capacity: 13.5 cu. ft / 368 litres
- Temperature Range: -121°F to -40°F / -85°C to -40°C
- Noise Level: 50 dBA
- Power Supply: 120V, 60Hz
- Power Consumption: 7.7 KWH/24H
- Weight: 194 lbs / 88 Kgs
- External Dimensions:
33.5 x 61.4 x 27.5 inches H x D x W
854 x 1562 x 698 mm H x D x W



Please contact us for a quotation for the hire or purchase of this product

Product description: -121°F Chest Freezer - 368 Litre

Product code: - WD7ULT044

Version: V2 10/19

Product specifications can be subject to change without prior notice.

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company