

# -130°F, 0.25 cu.ft Upright Ultra Low Temperature Freezer

Available for hire or purchase for your clinical trial.



## Overview

This compact 130°F 0.25 cu.ft capacity upright ultra low temperature freezer is ideal when space is limited.

It is environmentally safe product with no compressor or refrigerant. All features such as alarms and datalogging functions are included as standard. The freezer has been produced with the true and original single compress or technology.

With microprocessor controller and digital display. Direct download of logged data on USB memory stick and direct upload of new software via a USB memory stick.

## Operating Features

- Battery backup for alarm (approx 72 hours)
- Contact for remote alarm & prepared for GSM alarm module
- Auto cycle in case of probe failure
- 3 level password protected
- Ambient temperature display
- Temperature graph
- RS485/RS232 interface
- Integrated memory for 30 years
- Visual and acoustic alarm
- Shows all alarms
- Alarms: power failure, adjustable high/low temperature, open door and probe failure
- Alarm display as text not codes
- Manual defrost
- Porthole for installation of PT100CO2 back up
- Inner material is stainless steel
- Lighting
- Refrigerant: Ep88

### UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

### USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

[www.woodleytrialsolutions.com](http://www.woodleytrialsolutions.com)

## Specification

- Net Capacity: 0.25 cu. ft.
- Temperature Range: -130°F to -40°F
- Max. Ambient Temperature: 77°F
- Power Supply: 110V
- Frequency: 60Hz
- Power Consumption: 4.6 KWH/24H
- Weight: 83 lbs
- Internal Dimensions: 12.2 x 5.9 x 5.6 inches (HxWxD)
- External Dimensions: 26.2 x 15.7 x 24 inches (HxWxD)



**Please contact us for a quotation for the hire or purchase of this product**

Product description: -130°F Medical Upright Freezer

Product code: WD7ULT071

Version: V2 10/19

Product specifications can be subject to change without prior notice.

**woodley**  
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company