

-76 °F / -60 °C Chest Freezer 133 Litre



Available for hire or purchase for your clinical trial.

Overview

This is one of the most simple ultra low temperature freezers on the market.

With a single compressor the -76 °F / -60° C, 133 litre chest freezer provides low noise, low energy consumption and a compact design.

The chest freezer gives reliable low temperature freezing, and good temperature uniformity. Digital display, alarms, contact for remote monitoring, basket(s) and castors are some of the standard features included.

Operating Features

- Single compressor
- Digital display
- Castors, baskets and key lock
- Low heat dissipation
- Low energy consumption
- Low noise level
- Visual and acoustic alarms -
 - power failure
 - adjustable high/low temperature
 - probe failure
 - Instrument failure alarm
- Contact for remote alarm

Specification

- Capacity: 4.9 cu. ft / 133 litres
- Temperature Range: - 76°F to -40°F / -60°C to -40°C
- Noise Level: 50 dBA
- Power Supply: 110V, 60Hz
- Power Consumption: 6.6KWH/24H
- Weight: 130 lbs / 59 Kgs
- External Dimensions: 3
4.8 x 28.1 x 27.51 inches H x W x D
385 x 720 x 695 mm H x W x D



Please contact us for a quotation for the hire or purchase of this product

Product description: - 76°F Chest Freezer 133 Litre

Product code: -WD7ULT033

Version: V2 10/19

Product specifications can be subject to change without prior notice.

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company