

-76°F /-60°C Chest Freezer - 368 Litre

Available for hire or purchase for your clinical trial.



Overview

This counter height medical laboratory freezer has excellent temperature stability and management and is suitable for use in medical, scientific, pharmaceutical and laboratory facilities where precise temperature conditions are required.

The 5.0 cu. ft. 22" wide counter height freezer is designed for free standing use. Two open shelves include cooling coils for optimum performance, with additional door racks for easy storage.

Operating Features

- Single compressor
- Digital display
- Castors, baskets and key lock
- Low heat dissipation
- Low energy consumption
- Low noise level
- Visual and acoustic alarms -
 - power failure
 - adjustable high/low temperature
 - probe failure
 - Instrument failure alarm
- Contact for remote alarm

Specification

- Capacity: 13.5 cu. ft / 368 litres
- Temperature Range: - 76°F to -40°F /-60°C to -40°C
- Noise Level: 50 dBA
- Power Supply: 110V, 60Hz
- Power Consumption: 7.7 KWH/24H
- Weight: 194 lbs / 88 Kgs
- External Dimensions:
34.8 x 61.4 x 27.51 inches / H x W x D
895 x 1560 x 700 mm H x W x D



Please contact us for a quotation for the hire or purchase of this product

Product description: - 76°F, 368 Litre Chest Freezer 368 Litre:

Product code: - WD7ULT043

Version: V2 10/19

Product specifications can be subject to change without prior notice.

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company