

-86°C Upright Medical Freezer 54L

Available for hire or purchase for your clinical trial.



Overview

The 54L -86°C/-40°F freezer is compact in design making it suitable for the storage of samples in either the laboratory or office.

This versatile, low energy consumption freezer has approximately 72 hours battery back up for acoustic alarm and temperature display in case of power cut.

Operating Features

- Single compressor
- Low energy consumption
- Low noise level
- Low heat dissipation
- Castors
- Insulated inner doors
- Lock on door
- Heated door frame and heated vacuum valve
- Porthole for external temperature probe
- RS 232 & RS 485 Interface
- Power failure alarm
- Digital display
- Visual and acoustic alarms
- Adjustable high/low temperature alarm
- Open door alarm
- Contact for remote alarm
- Manual defrost
- USB
- Data logger

Specification

- Temperature Range: -40°C to -86°C
- Net Capacity: 54 litres/1.9 cu.ft
- Insulation: 4.7 Inches
- Maximum Ambient Temperature: 77°F
- Power Requirements: 110V / 60Hz
- AMP: 3.2 (A)
- Weight: 123 lbs
- Shipping Weight: 451 lbs
- 100% HCFC/CFC free
- Internal Dimensions: 12.6h x 14.2w x 18.4d inches
- External Dimensions: 31.9h x 23.6w x 27.6d inches



Please contact us for a quotation for the hire or purchase of this product

Product description: Upright -86°C Medical Freeze

Product Code: WD7ULT072

Version: V3 10/19

Product specifications can be subject to change without prior notice.

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company