

-86°C Chest Freezer - 71 Litres

Available for hire or purchase for your clinical trial.



Overview

This is one of the most simple ultra low temperature freezers on the market.

With a single compressor the -86°C chest freezer provides low noise, low energy consumption and a compact design.

The -86°C chest freezer demonstrates high performance freezing with an integrated cooling system.

Operating Features

- Single compressor
- Porthole for installation of PT100/CO2 back up
- Digital display
- Manual defrost
- Castors and key lock
- Low heat dissipation
- Low energy consumption
- Low noise level
- Visual and acoustic alarms -
 - power failure
 - adjustable high/low temperature
 - probe failure
- Contact for remote alarm
- Prepared for GSM alarm module
- Advanced battery back-up for alarms and temperature display (72 hours)
- 100% HCFC/CFC free
- Chart recorder (optional)

Specification

- Net Capacity: 71 litres
- Temperature Range: -85°C to -40°C
- Max. Ambient Temperature: 25°C
- Insulation: 80 mm
- Noise Level: 50 dBA
- Power Supply: 230V, 50/60Hz (also available in 110V)
- Power Consumption: 5.65 KWH/24H
- Weight: 58 kgs
- Internal Dimensions: 450h x 390w x 390d mm
- External Dimensions: 850h x 552w x 648d mm



Please contact us for a quotation for the hire or purchase of this product

Product description: -86°C Medical Chest Freezer

Product code - WD7ULT009

Version: V3 10/19

Product specifications can be subject to change without prior notice.

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company