

-90°C, 7L Upright Ultra Low Temp Freezer

Available for hire or purchase for your clinical trial.



Overview

This compact -90°C, 7 litre capacity upright ultra low temperature freezer is ideal when space is limited.

It is environmentally safe product with no compressor or refrigerant.

Operating Features

- Low energy consumption
- Low noise level
- Low heat dissipation
- Insulated inner doors
- Door key lock
- Ergonomic loading
- Heated door frame
- Heated vacuum valve
- 100% HCFC/CFC free
- Porthole for external temperature probes
- Visual and acoustic alarms
- Alarms: power failure, adjustable high/low temperature, open door and probe failure
- Alarm display as text not codes
- Contact for remote alarm and prepared for GSM alarm module
- Battery backup for alarms (approx 72 hours)
- Porthole for installation of PT100CO2 back up
- Inner material is stainless steel

Specification

- Net Capacity: 7 litres
- Temperature Range: -90°C to -40°C
- Max. Ambient Temperature: 25°C
- Power Supply: 230V
- Frequency: 50/60Hz
- Power Consumption: 4.6 KWH/24H
- Weight: 55 kgs
- Internal Dimensions: 310h x 150w x 143d mm
- External Dimensions: 665h x 400w x 630d mm



Please contact us for a quotation for the hire or purchase of this product

Product description: -90°C Medical Upright Freeze

Product code: WD7ULT047

Version: V2 10/19

Product specifications can be subject to change without prior notice.

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company