

Accu-Chek® Aviva Glucometer

Available for hire or purchase for your clinical trial.



Overview

The combination of features in the Accu-Chek® Aviva glucometer makes testing fast and easy; the system's technology allows you to perform an accurate test every time.

The Accu-Chek® Aviva system includes:

- Fast, convenient testing
- Reminders & averaging help you manage testing
- 500-result memory capacity
- Rubber grips for easy use

Operating Feature

- User-friendly design: small and slender, with easy-to-hold rubber grips, the Accu-Chek® Aviva meter offers a comfortable, discreet choice for blood glucose testing
- Small sample volume: quick fill, reliable Accu-Chek® Aviva test strips
- Test reminders: the user can set up to 4 automatic reminders for different times throughout the day, to help remember when to test
- Battery backup: an additional small battery in the meter keeps the time and memory for up to 72 hours when the main battery is removed or out of power
- Preset time and date: the meter now comes with the time and date preset, so startup time is faster for the user, as well as the healthcare professional helping with initial setup

Specification

- Display: LCD
- Automatic power:
 - auto power on with strip insertion
 - auto off 90 seconds after test if strip is not removed
 - auto off 5 seconds after strip removed
 - auto off in 2 minutes if meter not used
- Sample size: 0.6 µL
- Power supply: one 3-volt lithium battery (type CR2032)
- Battery life: approx. 1,000 readings (1 year with 3 tests per day)
- Measurement time: approx. 5 seconds
- Operating temperature: 6°C to 44°C
- Storage conditions: -25°C to +70°C
- Operating humidity: 10 - 90% RH
- Memory: 500 results and 20 control tests with date and time
- Test averages: 7, 14, 30 day
- Weight: 60g with battery
- Dimensions: 94h x 53w x 22d mm



Please contact us for a quotation for the hire or purchase of this product

Product description: Accu-Chek® Aviva Glucometer

Product Code:

Product specifications can be subject to change without prior notice. Version: V2 10/19

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company