

AIA-360 Analyser

Take a 360° look at your endocrinology testing with the Tosoh AIA-360



COMPACT BENCHTOP
ENDOCRINOLOGY TESTING.



"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"

WOODLEY
EQUIPMENT COMPANY LTD.


AIA-360 Analyser

The AIA-360 offers the automation powers and performance of a large analyser in a benchtop format. Random access allows any combination of tests, in any order, at any time.

A throughput of 36 results per hour, with the first result in 20 minutes ensures a fast, reliable service.

FEATURES

- Will accept a variety of primary tubes with positive patient identification as well as positive test cup ID to eliminate errors
- Sample level sensing to alert operator to short samples
- Dual clot detection to maintain result integrity
- Long calibration stability
- Simple to use
- Fast, reliable results
- Simple maintenance
- Touchscreen with intuitive menus
- Integrates with Practice Management Systems



KEY ASSAYS AVAILABLE

- Free T4
- Total T4
- Cortisol
- Progesterone
- ACTH
- Insulin
- cTnl
- iPTH

OPERATIONAL SIMPLICITY

- Place the samples/cups on-board and press START
- Touch screen menus
- Automated maintenance

TRACEABILITY

- On-board storage of results with export possibilities
- Real-time printing for fast result reviews

FLEXIBILITY

- Adapted to low volume use
- Long reagent shelf life

SAFETY/RELIABILITY

- Positive identification of samples and reagents
- Quantitative measurement
- Simple maintenance

DIMENSIONS

- 52 (H) x 40 (W) x 40 (D) cm



**A Ready to Use
"All-In-One" Test Cup**



+44 (0) 1204 669033 (Option 1)



sales@woodleyequipment.com



www.woodleyequipment.com



Old Station Park Buildings
St. John Street
Horwich
Bolton
Lancashire
BL6 7NY
UK

**"PASSIONATE ABOUT
LABORATORY DIAGNOSTICS"**

**WOODLEY
EQUIPMENT COMPANY LTD.**