

Alaris® GP Volumetric Pump

Available for hire or purchase for your clinical trial.



Overview

The Alaris® GP volumetric pump provides a versatile and flexible platform for infusion therapy which is suitable across a broad range of applications.

The Alaris® GP volumetric pump offers a range of features suited to general fluid, blood transfusion & drug therapy including: intuitive colour coded set loading, easy to clean with fast release door, a wide range of alarms, prompts & warnings.

Infusion rate: 1 - 1200 ml/hr.

Volume Infused range: 1 - 50,000 ml.

Operating Features

- Alarm beacon with alarm (red) & warning (amber) indicators
- Audible alarm and display messaging with comprehensive protection
- Mechanism: four finger displacement peristaltic
- Occlusion Settings: 3 levels from L2, L5 (default), L8
- Designed to fit into the Alaris® Gateway Workstation
- Easy to clean
- Fast release door
- Large display
- User friendly interface
- Door Removal: for cleaning and maintenance
- Various clamps, valves, filters and lines available

Please contact us for a quotation for the hire or purchase of this product

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

Specification

- Infusion Rate Range: 1 - 1200 ml/h
- VTBI Range Primary: up to 10000 ml
- Volume Infused Range: 1 ml to 50,000 ml
- Volumetric Accuracy: $\pm 5\%$ (nominal)
- Display: mono 240 x 240 pixels
- KVO Rates: 5 ml/h or set infusion rate if lower
- Battery Type: rechargeable NiMH 9.6V 2.4Ah nominal
- Battery Run Time: operates for 4 hours from full charge at 125 ml/h
- Recharge Time: 2.5 hour from discharge to 95% charge
- AC Power Requirements: 100-230 VAC, 50-60Hz, 60VA (maximum)
- Weight: 2.5kgs
- Dimensions: 225h x 140w x 148d mm



Product description: Alaris® GP Volumetric Pump

Product code - WD7INF024

Version: V3 10/19

Product specifications can be subject to change without prior notice

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company