

# Ambulatory Blood Pressure Monitor

Available for hire or purchase for your clinical trial.



## Overview

These Ambulatory Blood Pressure Monitors provide fast, accurate results enabling users to fine-tune their patients treatment. Extensive testing by independent clinicians has validated their accuracy time and time again.

Their reliability has made them one of the most popular ambulatory monitors in the world. They have a compact design, are lightweight and have a quiet operation to maximise patient comfort.

The monitor has advanced reporting software that allows rapid analysis of patient information.

## Operating Features

- Choice of 90207 or UltraLite™ 90217 ABP monitors
- The 90207 monitor has 24 hour monitoring
- The UltraLite™ 90217 monitor has 48 monitoring
- Real-time evaluation of ABPM data
- Individual site configuration
- Compact, lightweight design
- Clinically validated by official organisations in the US, UK, France and Germany
- 5 cuff sizes available (optional)
- ABP report generator/printer (optional)
- ABP report management system (optional):
  - Interfaces with 90207 and 90217 ABP monitors
  - Software offers quick and simple method of retrieving, displaying, editing and printing ABP data
  - Configurable report format and enhanced graphics
  - Advanced search capability
  - Remote communications feature for programming and data retrieval via modem
  - Supports e-mail transfer of reports as PDF files
  - Supports compliance with FDA electronic records requirements and HIPAA

### UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

### USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

[www.woodleytrialsolutions.com](http://www.woodleytrialsolutions.com)

## Specification

- Various cuff sizes available from 12 - 50 cm
- Disposable cuff covers (optional)
- Various Accessories Available:  
Variety of cuffs, report generator/printer, shoulder strap, belt, pouch, cloth and T-tube calibration accessory
- ABP Report Management System (optional):  
Windows compatible, remote communication capability, monitor programming via modem
- 90207 Monitor:
  - Power supply: 4 x AA batteries
  - Weight: 347g including batteries
  - Dimensions: 28h x 86w x 114d mm
- 90217 UltraLite™ Monitor:
  - Power supply: 3 x AA batteries
  - Weight: 255g including batteries
  - Dimensions: 25h x 70w x 100d mm



**Please contact us for a quotation for the hire or purchase of this product**

Product description: Ambulatory BP Monitor

Product Code: WD7BPM029

Version 4 10/19



Product specifications can be subject to change without prior notice

**woodley**  
TRIAL SOLUTIONS

Woodley Trial Solutions is a divisions of Woodley Equipment Company