

# -Basic Counter Top Freezer

## 1.4 Cu.Ft. / 40 L

Available for hire or purchase for your clinical trial.



### Overview

This compact size freezer is ideal for use on a counter top. The factory installed lock provides security when it is needed.

Designed for general purpose use in medical and scientific facilities.

The reversible door offers flexible placement of the freezer.

### Operating Features

- Compact Size
- Factory installed lock
- Manual defrost with a static cooling system
- One piece interior liner
- 1 removable wire shelf
- Adjustable thermostat for easy temperature management
- Flat door liner
- 100% CFC Free
- 2 level legs
- US Electrical Safety - ETL
- Canadian Electrical Safety - ETL-C

### Specification

- Capacity: 1.4 Cu. Ft. / 40 Litres
- Temperature Range: -20 to -15 °C / -5 to +5 °F
- Amps: 0.8
- Input Voltage: 115 V AC/60 Hz
- Weight: 31 lbs
- Dimensions External: 19.63" h x 18.5" w x 17.63" d
- Dimensions Internal: 14.38" h x 13.5" w x 13.25" d



**Please contact us for a quotation for the hire or purchase of this product**

Product description: Basic Counter Top Freezer  
Freezer 1.4 cu. ft.

Product Code: WD7FRZ086

Version: V2 10/19

Product specifications can be subject to change without prior notice.

#### UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

#### USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: [hello@woodleytrialsolutions.com](mailto:hello@woodleytrialsolutions.com)

[www.woodleytrialsolutions.com](http://www.woodleytrialsolutions.com)

**woodley**  
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company