



Vcheck V200

ACHIEVE TRUE ANALYSIS WITH FLUORESCENCE TECHNOLOGY

- In-Practice Point of Care Testing
- Fully quantitative fluorescence analyser



Canine Tests Available:

- cPL
- T4
- TSH
- Cortisol
- CRP
- Progesterone

Feline Tests Available:

- T4
- fPL
- SAA

**LOW COST,
REFERENCE LAB QUALITY,
IMMUNOASSAY TESTING.**



"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"

WOODLEY
EQUIPMENT COMPANY LTD.

BIONOTE

Vcheck V200

Simplify In-Practice Point of Care Testing with Reference Lab Quality - avoiding test referrals

- Fully quantitative automated analyser
- Accurate diagnosis for Dogs and Cats in minutes
- Analyser provides print out and result storage



Canine Tests Available:

- cPL
- T4
- TSH
- Cortisol
- CRP
- Progesterone

Feline Tests Available:

- T4
- fPL
- SAA

BIONOTE Anigen™ Rapid Test Kits

ALSO AVAILABLE
EXCLUSIVELY
FROM WOODLEY



SPECIFICATIONS

Testing Time:	5 minutes (Based on CRP, SAA)
Weight:	2.5kg
Screen:	Touchscreen 7" Colour TFT
Size:	20 x 24 x 20 cm
Printer:	Built-in Printer
Sample Required:	Serum or Plasma (Heparin)

+44 (0) 1204 669033 (Option 1)

sales@woodleyequipment.com

www.woodleyequipment.com

Old Station Park Buildings
St. John Street
Horwich
Bolton
Lancashire
BL6 7NY
UK

**"PASSIONATE ABOUT
LABORATORY DIAGNOSTICS"**

**WOODLEY
EQUIPMENT COMPANY LTD.**