

Clearview® hCG Pregnancy Test



Available for hire or purchase for your clinical trial.

Overview


The Clearview® hCG pregnancy test is the one-step test for high sensitivity detection of hCG in urine. It provides rapid, reliable results in 3 minutes.

It has the ability to detect pregnancy on the day of a missed period and has an inbuilt procedural control for confirmation that the test has been performed correctly.

The test can be stored at room temperature so it is ready to use and leaves valuable refrigerator space free. It also has a long shelf life so small volume users need not worry that the product will expire before use.


Operating Feature

- Results in 3 minutes
- Sensitive to 25ml U/ml
- One step - no additional reagents, training or equipment required
- Proven accuracy of > 99%
- Inbuilt procedural control
- Room temperature storage
- Long shelf life
- 'One shot' pipette for exact specimen volume
- CE marked
- Available in 12 languages



Specification

- Method: lateral flow
- Time To Result: 3 minutes
- Sample Size: 100 µl
- Sample Type: urine
- Sample Procedure: pipette
- Storage Temperature: 2°C to 30°C
- Shelf Life: 24 months from manufacture
- Sensitivity: 25 ml/Uml
- Specificity: 100%
- Accuracy: > 99%
- Kit Size: 20 tests
- Kit Contents:
 - 20 pregnancy test units
 - 20 disposable pipettes
 - 1 package insert



Please contact us for a quotation for the hire or purchase of this product

Product description : Clearview® hCG Pregnancy Test
Product specifications can be subject to change without prior notice. Version: V2 10/19

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company