

# Clinispin Horizon 755VES 110V Centrifuge



Available for hire or purchase for your clinical trial.

## Overview

The Clinispin horizon 755VES features the largest capacity of the horizon range. The versatile standard horizontal rotor can accommodate up to 24 tubes.

The 755VES offers control over almost every aspect of the unit's operation from variable acceleration and deceleration to easy time and speed input.

The unique Performance Plus horizontal rotor accommodates 24 x 10ml tubes, 12 x 15ml tubes or 6 x 50ml tubes.

## Operating Feature

- Electronic timer: selectable from 1 to 99 minutes ( $\pm 1\%$ )
- Brushless DC Motor: for smooth and quiet operation with no routine maintenance
- Variable Acceleration and Deceleration: customise the run from start to finish
- Save settings in one of ten memory locations and recall them at the touch of a button
- End-of-run Audible Indicator: sounds when the rotor has stopped to inform when cycle is complete
- Quiet Operation and Low Vibration: double-encased motor and unique "stability mount" assure quiet performance and low out-of-balance vibration
- Automatic Lid-Lock: centrifuge will not run with the lid open; lid will not open while samples are spinning
- Positive Lid-Lock: firm mechanical lid-lock prevents the escape of material in case of breakage or accident
- Welded Steel Construction: for strength and durability
- Optional Rotors available:  
6 place fixed angle rotor & 20 place fixed angle rotor

## Specification

- Rotor: Performance Plus horizontal rotor
- Variable Speed: 3500 rpm / 2200 x g
- Min. Acceleration time: 15 seconds
- Min. Deceleration time: 15 seconds
- Electronic Timer: 99 minutes
- Weight: 17.7 kg
- Dimensions: 230h x 370w x 430d mm
  
- 4 amp protection breaker
- Max. power requirement: 200 watts
- Voltage: 110V / 230V
- Frequency: 50/60 Hz
- CE, UL and CAN/CSA certified



**Please contact us for a quotation for the hire or purchase of this product**

Product description - Clinispin horizon 755VES Centrifuge  
Performance Plus Horizontal Rotor product code: WD7CEN119  
6 place Fixed Angle Rotor product code: WD7CEN120  
20 place Fixed Angle Rotor product code: WD7CEN121  
Version: V2 10/19

Product specifications can be subject to change without prior notice

### UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

### USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

[www.woodleytrialsolutions.com](http://www.woodleytrialsolutions.com)

**woodley**  
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company